


Dr.Maisch GmbH VisionHT™ Ultra High-Pressure Columns


See Separations with Speed and Clarity

- Ultra-fast separations with superior efficiency, sensitivity and resolution
- Exceptional stability for long column life times
- Comprehensive sub2 µm stationary phase offering
- 12,000psig pressurerating compatible with all ultrahigh-pressure LC systems


Dr. Maisch




VisionHT™ columns offer a new level of performance in HPLC. The powerful combination of high-strength 1.5µm media with ultra-low volume hardware, delivers new clarity to your separations and maintains exceptional column lifetime. Complex samples resolve 95% faster with 4x greater sensitivity when compared to traditional 2.1 x 150mm, 5µm columns. And with a wide variety of phases available, the possibilities are endless.



Efficiency and Retention Times Remain Constant after Exposure to 12,000psig and 2000 Injections



High-pressure competitor sub 2µm column lost efficiency at routine 12,000psig pressures. VisionHT™ columns remain stable.



Before and after chromatograms show a constant level of performance after 2000 injections.



VisionHT™ Ultra High-Pressure Columns

Comprehensive High Purity 1.5µm Selectivity Options


Sub 2µm particles deliver efficiency and speed, but critical to success is having the right stationary phase selectivity. Six VisionHT™ high purity phases are available, each with unique separation benefits. C18-HL, with maximum bonded phase coverage, is ideal for complex hydrophobic samples. Use C18-B for basic compounds at neutral pH, often a requirement for mass spec work in the pharmaceutical industry. Reserve the C18 and C18-P for fastest analysis times. Both offer increased polar interactions to make neutral, non-polar compounds elute faster and retain polar compounds longer. The HILIC and Silica packings are normal phases that typically use near exclusive organic mobile phases; an advantage when seeking highest mass spec sensitivity.

VisionHT™ phases have exceptionally rigid silica structure that withstands routine use of 12,000psig pressure often required in this new technology. Reproducible silica synthesis, bonding and column packing guarantee low bleed and excellent column-to-column consistency.


VisionHT™ Phase Specifications									
Packing	Base Material	Particle Size	Carbon Load	Pore Size	Surface Area	Endcapped	pH Range*	Feature	Benefits
C18-HL	Spherical Silica	1.5µm	10%	120Å	220m ² /g	Yes	1–10	Fully bonded silica, ultra high purity silica	High capacity for hydrophobic compounds. Good peak symmetry.
C18-B	Spherical Silica	1.5µm	5%	120Å	220m ² /g	Proprietary	1–10	High purity silica and unique endcapping.	Improved performance of basic compounds at neutral pH. Better sensitivity and peak shape by mass spec for basic compounds, without the need for acidified mobile phases.
C18	Spherical Silica	1.5µm	6%	100Å	200m ² /g	Yes	1–10	Moderate silica exposure	Classic reversed-phase selectivity. Reduced bonding is optimized for speed and sensitivity.
C18-P	Spherical Silica	1.5µm	5%	100Å	200m ² /g	No	1–10	High silica exposure, low carbon load	Unique polar selectivity. Low carbon load gives fastest reversed-phase elution times while retaining polar compounds longer.
HILIC	Spherical Silica	1.5µm	NA	120Å	220m ² /g	No	2–8	Polar phase with shorter equilibration times. Shipped in ACN/Water.	Peak reversal compared to reversed-phase. Ideal for very polar compounds with high organic mobile phases for improved sensitivity by MS.
Silica	Spherical Silica	1.5µm	NA	120Å	220m ² /g	No	2–8	Traditional normal phase for use in 100% organic mobile phases	For isomeric separation of non-aqueous compatible compounds by absorption chromatography.

*Choice of buffer is critical at pH > 8.

VisionHT Selectivity



Separate Highly Basic Components at neutral pH with Vision HT™ C18-B



Column: VisionHT™ C18-B, 1.5µm, 2.0x50mm
 Mobile Phase: 50mM Ammonium Formate pH7: Methanol(20:80)
 Column Temp: 40°C
 Flow Rate: 0.2mL/min
 Detector: UV at 210nm




VisionHT™ Ultra High-Pressure Columns

Unique Column Hardware Advantages


You can erode the benefits of sub 2µm media if packed in an inefficient column design. VisionHT™ hardware uses thin screens, instead of traditional frits, to retain media and minimize dead volume. A unique insert combination seals end fittings leak-free to 18,000psig. Various dimensions offer even more flexibility over analysis speed and resolution.

VisionHT™ Hardware Specifications

Dead volume:	<15nL
Wetted materials:	316 Stainless Steel, PTFE
Port Geometry:	10–32, for Industry Standard & Waters® connection
Pressure rating, hardware alone:	18,000psig (1250 bar)
Pressure rating, packed column:	12,000psig (830 bar)
i.d.:	1.0, 2.0mm
Length:	20, 30, 50, 100mm



Sulfa Drugs



Column: VisionHT®C18-P, 1.5µm, 50x2.0mm (Part No. 5139606)

Mobile Phase: A: 0.1%Formic Acid B: Methanol

Gradient:

Time:	0	1.4
%B:	20	100

Flow Rate: 0.4mL/min

Detector: UV at 280nm