




# Application - 1

|   | <b>Page</b> |
|---|-------------|
| • Cold medicine .....   | 22          |
| • Phenol mixture EPA 604/625                                      |             |
| • Tricyclic Antidepressants I                                     |             |
| • Fat soluble Vitamins II   |             |
| • Engelhardt test for RP phases .....                             | 23          |
| • Antirheumatic bath ingredient                                   |             |
| • High Speed Separation of Catecholamines                         |             |
| • Organic Acids I   |             |
| • Organic Acids II .....  | 24          |
| • Formic Acid / Acetic Acid                                       |             |
| • Water soluble Vitamins  |             |
| • Peptides I  |             |
| • Proteins I .....  | 25          |
| • Fat soluble Vitamins I  |             |
| • Insulin Derivatives I   |             |
| • Insulin Derivatives II  |             |
| • Isomers of Tocopherol (Vitamin E) .....                         | 26          |
| • Tricyclic Antidepressants II                                    |             |
| • Steroids - UV vs. Evaporative Light Scattering Detection (ELSD) |             |
| • Carotenoids I   |             |
| • N-Methylated Xanthines .....                                    | 27          |
| • Parabens  |             |
| • Carbohydrates I   |             |
| • Carbohydrates II Evaporative Light Scattering Detection         |             |
| • Amino Acids I .....   | 28          |
| • Weak Anions - Ion Pair Chromatography                           |             |
| • Iodate in Galvanic Bath   |             |
| • Proteins II   |             |
| • Peroxide .....  | 29          |
| • 5-Aminolevulinic Acid   |             |
| • Alcohols I  |             |
| • Dihydroxybenzoic Acids  |             |
| • Bitter ingredients of hop I .....                               | 30          |
| • Bitter ingredients of hop II                                    |             |
| • Melatonin   |             |
| • Aldehydes and Ketones I: DNPH Derivatives in Galvanic Bath      |             |
| • Fast Separation of Proteins .....                               | 31          |
| • Fast Separation of Peptides                                     |             |




## Cold medicine ProntoSIL Eurobond C18 5 µm

|             |   |             |  |
|-------------|---|-------------|--|
| Part number | 1204F181PS050   | Detection   | UV 230nm, 0-240 s;<br>UV 261nm, 240-720 s                                |
| Dimension   | 125 x 4.0 mm  | Temperature | 20 °C  |
| Eluent      | A: 5mM Li <sub>2</sub> SO <sub>4</sub> in H <sub>2</sub> O4; pH 2.1<br>B: ACN/50mM H <sub>3</sub> PO <sub>4</sub> | Injection   | 3 µl   |
| Gradient    | 0%B, 0-120 s; 0_12%B, 120-184 s;<br>12-17% B, 184-500 s;<br>17-39% B, 500-720 s                                   | Sample      | 1 : Ascorbic acid<br>2 : Paracetamol<br>3 : Coffein<br>4 : Chlorphenamin |
| Flow        | 1 ml/min  |             |  |




## Phenol mixture EPA 604/625 ProntoSIL EnviroPHE

|             |  |        |  |
|-------------|--|--------|--|
| Part number | 1246F440PS30                                 | Sample | Supelco Calibration Standard<br>(P/N: 4-8859)  |
| Dimension   | 125 x 4.6 mm                                 |        | 1: Phenol  |
| Eluent      | A: H <sub>2</sub> O/1% HAC<br>B: MeOH/1% HAC |        | 2: 4-Nitrophenol<br>3: 2,4-Dinitrophenol<br>4: 2-Chlorophenol<br>5: 2-Nitrophenol<br>6: 2,4-Dimethylphenol<br>7: 4,6-Dinitro-o-Cresol<br>8: 2,4-Dichlorophenol<br>9: 4-Chloro-m-Cresol<br>10: 2,4,6-Trichlorophenol<br>11: Pentachlorophenol |
| Gradient    | 5-100% B, 0-30min                            |        |  |
| Flow        | 1 ml/min                                     |        |  |
| Detection   | UV 280 nm                                    |        |  |
| Temperature | 25 °C  |        |  |
| Injection   | 2 µl   |        |  |




## Tricyclic Antidepressants I ProntoSIL 120-5-CN

|             |   |               |   |
|-------------|---|---------------|---|
| Part number | 2546F200PS050   | Injection     | 5 µl  |
| Dimension   | 250 x 4.6 mm  | Concentration | 50 ppm each   |
| Eluent      | A: 25mM K <sub>2</sub> HPO <sub>4</sub> (pH 7.1)<br>B: MeOH/ACN (15/65) | Sample        | 1 : Trimipramine<br>2 : Doxepin<br>3 : Amitriptyline<br>4 : Imipramine<br>5 : Nortriptyline<br>6 : Desipramine<br>7 : Protriptyline |
| Flow        | 1 ml/min  |               |   |
| Detection   | UV 254nm  |               |   |
| Temperature | 40 °C   |               |   |



## Fat soluble Vitamins II ProntoSIL 120-3-C18 SH


|             |   |               |   |
|-------------|---|---------------|---|
| Part number | 2503F180PS030   | Temperature   | 20 °C   |
| Dimension   | 250 x 3.0 mm  | Injection     | 5 µl  |
| Eluent      | MeOH  | Concentration | 50-300 ppm each   |
| Flow        | 1 ml/min  | Sample        | 1 : Vitamin A<br>2 : Vitamin A acetate<br>3 : Vitamin D <sub>2</sub><br>4 : Vitamin D <sub>3</sub><br>5 : Vitamin E<br>6 : Vitamin K <sub>1</sub> |
| Detection   | Evap. Light Scattering<br>Detector (DDL 31)<br>PMT: 600, T: 33 °C |               |   |






## Engelhardt test for RP phases ProntoSIL Eurobond C18 5 µm

| Part number 1204F181PS050 |                                   | Sample |           |       |       |       |                      |  |
|---------------------------|-----------------------------------|--------|-----------|-------|-------|-------|----------------------|--|
| Dimension                 | 125 x 4 mm                        | Nr.    | Ret.(min) | TP    | TP/m  | Asym. | Name                 |  |
| Eluent                    | H <sub>2</sub> O/MeOH 51/49 (w/w) | 1      | 0.96      | -     | -     | -     | Uracil               |  |
| Flow                      | 1 ml/min                          | 2      | 1.62      | 6756  | 54048 | 1.49  | Aniline              |  |
| Detection                 | UV 254nm                          | 3      | 1.84      | 7117  | 56937 | 1.33  | Phenol               |  |
| Temperature               | 40 °C                             | 4      | 3.37      | 7433  | 59461 | 1.34  | p-Ethylaniline       |  |
| Injection                 | 5 µl                              | 5      | 6.75      | 9450  | 75599 | 1.12  | N,N-Dimethyl-aniline |  |
|                           |                                   | 6      | 9.33      | 10380 | 83037 | 1.13  | Toluene              |  |
|                           |                                   | 7      | 16.24     | 10357 | 82856 | 1.05  | Ethylbenzene         |  |



## Antirheumatic bath ingredient ProntoSIL Eurobond C18 5 µm


| Part number 2504F181PS050 |  |
|---------------------------|--|
| Dimension                 | 250 x 4 mm   |
| Eluent                    | H <sub>2</sub> O/MeOH 30/70 (w/w),<br>0.1% acetic acid   |
| Flow                      | 0.9 ml/min   |
| Detection                 | UV 254nm   |
| Temperature               | 20 °C  |
| Injection                 | 20 µl  |
| Sample                    | 1 : Camphor<br>2 : Nicotinic acid methyl easter<br>3 : Nicotinic acid benzyl easter<br>4 : Ethylenegly colmonosalicylate |



## High Speed Separation of Catecholamines ProntoSIL 120-3-C18 AQ


### ProntoSIL 120-3-C18 AQ

| Part number 0404F184PS030 |   |
|---------------------------|---|
| Dimension                 | 33 x 4 mm   |
| Eluent                    | CAT-A-Phase   |
| Flow                      | 1.0 ml/min  |
| Detection                 | Coulometric   |
| Temperature               | 30 °C   |
| Injection                 | 5 µl  |
| Sample                    | 1 : Norepinephrine<br>2 : Epinephrine<br>3 : Dihydroxybenzylamine<br>4 : Dopamine |




### ProntoSIL 120-3-C18 AQ

| Part number 0203F184PS030 |   |
|---------------------------|---|
| Dimension                 | 25 x 3 mm   |
| Eluent                    | ADAM-B-Phase  |
| Flow                      | 0.95 ml/min   |
| Detection                 | Coulometric   |
| Temperature               | 28 °C   |
| Injection                 | 1 µl  |
| Sample                    | 1 : Norepinephrine<br>2 : Epinephrine<br>3 : Dihydroxybenzylamine<br>4 : Dopamine |



## Organic Acids I ProntoSIL 120-3-C18 AQ


| Part number 2503F184PS030 |                                      | Sample |               |
|---------------------------|--------------------------------------|--------|---------------|
| Dimension                 | 250 x 3.0 mm                         | 1 :    | Glutamic acid |
| Eluent                    | 50 mM H <sub>3</sub> PH <sub>4</sub> | 2 :    | Oxalic acid   |
| Flow                      | 0.7 ml/min                           | 3 :    | Tartaric acid |
| Detection                 | UV 205nm                             | 4 :    | Malic acid    |
| Temperature               | 22 °C                                | 5 :    | Ascorbic acid |
| Injection                 | 5 µl                                 | 6 :    | Acetic acid   |
|                           |                                      | 7 :    | Maleic acid   |
|                           |                                      | 8 :    | Citric acid   |
|                           |                                      | 9 :    | Fumaric acid  |






## Organic Acids II ProntoSIL 120-3-C18 AQ

|             |  |        |  |
|-------------|--|--------|--|
| Part number | 2003F184PS030  | Sample | 1: Oxalic acid                                   |
| Dimension   | 300 x 3.0 mm   |        | 2: Tartaric acid                                 |
| Eluent      | A: H <sub>2</sub> O/ 50 mM H <sub>3</sub> PO <sub>4</sub><br>B: ACN/ 50mM H <sub>3</sub> PO <sub>4</sub> |        | 3: Pyridine-4- carboxylic acid                   |
| Gradient    | 0-100% B, 0-20 min   |        | 4: Pyridine-3- carboxylic acid                   |
| Flow        | 0.6 ml/min   |        | 5: Lactic Acid                                   |
| Detection   | UV 215 nm  |        | 6: Citric Acid                                   |
| Temperature | 25 °C  |        | 7: Pyridine-2,6- dicarboxylic acid               |
| Injection   | 5 µl   |        | 8: 4- hydroxybenzoic acid                        |
|             |  |        | 9: 3- Hydroxybenzoic acid                        |
|             |  |        | 10: Benzoic acid                                 |
|             |  |        | 11: 2- Hydroxybenzoic acid<br>(= Salicylic acid) |




## Formic Acid / Acetic Acid ProntoSIL 120-3-C18 AQ

|             |                                       |
|-------------|---------------------------------------|
| Part number | 2503F184PS030                         |
| Dimension   | 250 x 3.0 mm                          |
| Eluent      | 100 mM H <sub>3</sub> PO <sub>4</sub> |
| Flow        | 0.7 ml/min                            |
| Detection   | UV 202 nm                             |
| Temperature | 22 °C                                 |
| Injection   | 5 µl                                  |
| Sample      | 1: Formic acid<br>2: Acetic acid      |




## Water soluble Vitamins ProntoSIL 120-3-C18 AQ

|             |   |        |   |
|-------------|---|--------|---|
| Part number | 2003F184PS030   | Sample | 1: Pyridoxamine   |
| Dimension   | 200 x 3.0 mm  |        | 2: Thiamine (Vit. B1)+impurity  |
| Eluent      | A: H <sub>2</sub> O/ 50 mM H <sub>3</sub> PO <sub>4</sub><br>B: ACN |        | 3: Ascorbic acid (Vit. C)   |
| Gradient    | 0% B, 0-180 s; 0-30% B,<br>181-400 s; 30% B, 401-790 s              |        | 4: Nicotinic acid (Niacin)  |
| Flow        | 0.7 ml/min  |        | 5: Nicotinamide (Vit. B <sub>3</sub> )                                |
| Detection   | UV 230 nm   |        | 6: Pyridoxal  |
| Temperature | 22 °C   |        | 7: Pyridoxine (Vit. B <sub>6</sub> )+impurity of Vit. B <sub>12</sub> |
| Injection   | 5 µl  |        | 8: Folic acid   |
|             |   |        | 9: Cyanocobalamin (Vit. B <sub>12</sub> )                             |
|             |   |        | 10: Riboflavin (Vit. B <sub>2</sub> )                                 |



## Peptides I ProntoSIL 300-5-C18 H


|             |  |             |                   |
|-------------|--|-------------|-------------------|
| Part number | 2546K185PS050  | Temperature | 30 °C             |
| Dimension   | 250 x 4.6 mm   | Injection   | 5 µl              |
| Eluent      | A: H <sub>2</sub> O/ 0.1% TFA<br>B: ACN/ H <sub>2</sub> O, 70/30 (v/v)<br>0.1% TFA | Sample      | 1: Oxytocin       |
| Gradient    | 20-100% B, 0-30 min  |             | 2: Bradykinin     |
| Flow        | 1.5ml/min  |             | 3: Angiotensin II |
| Detection   | UV 220 nm  |             | 4: Eledoisin      |
|             |  |             | 5: Neurotensin    |
|             |  |             | 6: Angiotensin I  |






## Proteins I ProntoSIL 300-5-C18 H

|             |   |           |                     |
|-------------|---|-----------|---------------------|
| Part number | 2546K185PS050                                     | Injection | 5 µl                |
| Dimension   | 250 x 4.6 mm                                      | Sample    | 1 : Ribonuclease A  |
| Eluent      | A: H <sub>2</sub> O / 0.1% TFA                    |           | 2 : Insulin, bovine |
|             | B: ACN / H <sub>2</sub> O, 95/5 (v/v)<br>0.1% TFA |           | 3 : Lysozyme        |
| Gradient    | 25-80% B, 0-17min                                 |           | 4 : BSA             |
| Flow        | 1.5 ml/min  |           | 5 : Myoglobin       |
| Detection   | UV 220 nm   |           | 6 : Ovalbumin       |
| Temperature | 30 °C   |           |                     |




## Fat soluble Vitamins I ProntoSIL 120-3-C18 H

|             |               |        |                  |
|-------------|---------------|--------|------------------|
| Part number | 2003F185PS030 | Sample | 1 : Vitamin A    |
| Dimension   | 200 x 3.0 mm  |        | 2 : Vitamin D2   |
| Eluent      | Acetonitrile  |        | 3 : Vitamin D3   |
| Flow        | 1 ml/min      |        | 4 : α-Tocopherol |
| Detection   | UV 280 nm     |        |                  |
| Temperature | 30 °C         |        |                  |
| Injection   | 5 µl          |        |                  |




## Insulin Derivatives I ProntoSIL 300-5-C18 H

|             |   |             |                       |
|-------------|---|-------------|-----------------------|
| Part number | 2546K185PS050                                     | Temperature | 30 °C                 |
| Dimension   | 250 x 4.6 mm                                      | Injection   | 5 µl                  |
| Eluent      | A: H <sub>2</sub> O / 0.1% TFA                    | Sample      | 1 : Protamine Sulfate |
|             | B: ACN / H <sub>2</sub> O, 95/5 (v/v)<br>0.1% TFA |             | 2 : Phenol            |
| Gradient    | 25-80% B, 0-17min                                 |             | 3 : Cresol            |
| Flow        | 1.5 ml/min  |             | 4 : Insulin H         |
| Detection   | UV 220 nm   |             | 5 : Protamine insulin |



## Insulin Derivatives II ProntoSIL 300-5-C18 H


|             |   |             |                               |
|-------------|---|-------------|-------------------------------|
| Part number | 2546K185PS050                                     | Temperature | 30 °C                         |
| Dimension   | 250 x 4.6 mm                                      | Injection   | 5 µl                          |
| Eluent      | A: H <sub>2</sub> O / 0.1% TFA                    | Sample      | 1 : Cresol                    |
|             | B: ACN / H <sub>2</sub> O, 95/5 (v/v)<br>0.1% TFA |             | 2 : Insulin H                 |
| Gradient    | 25-80% B, 0-17min                                 |             | 3 : Insulin Chain A, oxidized |
| Flow        | 1.5 ml/min  |             | 4 : Insulin, bovine           |
| Detection   | UV 220 nm   |             | 5 : Insulin Chain B, oxidized |






## Isomers of Tocopherol (Vitamin E) ProntoSIL 200-3-C30

|             |                                   |        |                         |
|-------------|-----------------------------------|--------|-------------------------|
| Part number | 2546H300P5030                     | Sample | 1: $\delta$ -Tocopherol |
| Dimension   | 250 x 4.6 mm                      |        | 2: $\gamma$ -Tocopherol |
| Eluent      | MeOH/H <sub>2</sub> O, 96/4 (v/v) |        | 3: $\beta$ -Tocopherol  |
| Flow        | 0.9 ml/min                        |        | 4: $\alpha$ -Tocopherol |
| Detection   | Coulometric                       |        |                         |
| Temperature | 25 °C                             |        |                         |
| Injection   | 5 $\mu$ l                         |        |                         |



## Tricyclic Antidepressants II ProntoSIL 120-3-C18 AQ


|             |   |             |                  |
|-------------|---|-------------|------------------|
| Part number | 0604F184P5030   | Temperature | 20 °C            |
| Dimension   | 53 x 4.0 mm   | Injection   | 5 $\mu$ l        |
| Eluent      | ACN/MeOH/NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub> (10 mM), 62/13/23 (v/v), titration with TEA to pH 5 | Sample      | 1: Doxepin       |
| Flow        | 0.69 ml/min   |             | 2: Imipramine    |
| Detection   | UV 254 nm   |             | 3: Nortriptyline |
|             |   |             | 4: Amitriptyline |
|             |   |             | 5: Trimipramine  |



## Steroids - UV vs. Evaporative Light Scattering Detection (ELSD)


### ProntoSIL 120-3-C18 H

|             |                                  |               |                    |
|-------------|----------------------------------|---------------|--------------------|
| Part number | 2003F185P5030                    | Concentration | 70 ppm each        |
| Dimension   | 200 x 3 mm                       | Sample        | 1: Hydrocortisone  |
| Eluent      | MeOH                             |               | 2: Progesterone    |
| Flow        | 0.65 ml/min                      |               | 3: Cholecalciferol |
| Detection   | a) UV 200 nm                     |               | 4: Ergocalciferol  |
|             | b) ELSD: PMT gain: 600. T: 40 °C |               | 5: Ergosterol      |
| Temperature | 30 °C                            |               | 6: Cholesterol     |
| Injection   | 5 $\mu$ l                        |               |                    |




## Carotenoids I ProntoSIL 200-3-C30

|             |                         |        |  |
|-------------|-------------------------|--------|--|
| Part number | 2546H300P5030           | Sample | $\beta$ -Carotene, technical mixture, isomerized |
| Dimension   | 250 x 4.6 mm            |        | 1: all-trans                                     |
| Eluent      | MeOH/ TBME, 80/20 (v/v) |        | 2: 9-cis   |
| Flow        | 1.4 ml/min              |        |  |
| Detection   | VIS 450 nm              |        |  |
| Temperature | 20 °C                   |        |  |
| Injection   | 5 $\mu$ l               |        |  |







## N-Methylated Xanthines ProntoSIL 120-5-C18 AQ




Part number 1504F184PS050  
 Dimension 150 x 4.0 mm  
 Eluent 25mM NaH<sub>2</sub>PO<sub>4</sub>, pH 3/  
 MeOH 80/20 (v/v)  
 Gradient 25-80% B, 0-17min  
 Flow 1.0 ml/min  
 Detection UV 254 nm  
 Temperature 35 °C  
 Injection 10 µl  
 Sample 1: Theobromine  
 2: Theophylline  
 3: Caffeine



## Parabens ProntoSIL 120-3-PHENYL

Part number 1204F050PS030  
 Dimension 125 x 4.0 mm  
 Eluent 20mM K<sub>2</sub>HPO<sub>4</sub> / ACN,  
 50/50 (v/v)  
 Flow 1.0 ml/min  
 Detection UV 254 nm  
 Temperature 40 °C


Injection 10 µl  
 Sample Parabens  
 1: Methylparaben  
 2: Ethylparaben  
 3: Propylparaben  
 4: Butylparaben



## Carbohydrates I ProntoSIL 120-3-NH<sub>2</sub>

Part number 1246F190PS030  
 Dimension 125 x 4.0 mm  
 Eluent H<sub>2</sub>O/ACN, 20/80 (v/v)  
 Flow 1.0 ml/min  
 Detection RI  
 Temperature Ambient  
 Injection 6 µl


Sample 1: Fructose  
 2: Mannose  
 3: Glucose  
 4: Galactose  
 5: Saccharose  
 6: Maltose



## Carbohydrates II Evaporative Light Scattering Detection ProntoSIL 120-5-NH<sub>2</sub>

Part number 1246F190PS030  
 Dimension 125 x 4.0 mm  
 Eluent H<sub>2</sub>O/ACN, 20/80 (v/v)  
 Flow 1.0 ml/min  
 Detection RI  
 Temperature Ambient  
 Injection 6 µl


Sample 1: Fructose  
 2: Mannose  
 3: Glucose  
 4: Galactose  
 5: Saccharose  
 6: Maltose






## Amino Acids I ProntoSIL 120-3-C18 H

|             |   |             |                        |
|-------------|---|-------------|------------------------|
| Part number | 2003F185PS030   | Flow        | 0.6 ml/min             |
| Dimension   | 200 x 3.0 mm  | Detection   | UV 340 nm              |
| Eluent      | A: 20 mM CH <sub>3</sub> COONa in H <sub>2</sub> O/ACN, 93/3 (v/v)  | Temperature | 30 °C                  |
|             | B: 20 mM CH <sub>3</sub> COONa in H <sub>2</sub> O/ACN, 50/50 (v/v) | Injection   | 1 µl                   |
| Gradient    | 5-28% B, 0-144 s  | Sample      | SIGMA standard AA-s-   |
|             | 28-45% B, 145-560 s   |             | 18: OPA/ Mercapto-     |
|             | 45-82% B, 561-896 s   |             | propionic acid labeled |
|             | 82-90% B, 897-1215 s  |             |                        |




## Weak Anions - Ion Pair Chromatography ProntoSIL 120-3-C8 SH

|             |   |           |                 |
|-------------|---|-----------|-----------------|
| Part number | 2003F080PS030   | Injection | 5 µl            |
| Dimension   | 200 x 3.0 mm  | Sample    | 1 : acetic acid |
| Eluent      | 10mM KH <sub>2</sub> PO <sub>4</sub> /10 mM TBA-H <sub>3</sub> PO <sub>4</sub> , pH 2.4 |           | 2 : lactic acid |
|             |   |           | 3 : malic acid  |
| Flow        | 0.6 ml/min  |           | 4 : iodate      |
| Detection   | UV 205 nm   |           | 5 : citric acid |
| Temperature | 25 °C   |           |                 |




## Iodate in Galvanic Bath ProntoSIL 120-3-C18 SH

|             |   |           |                          |
|-------------|---|-----------|--------------------------|
| Part number | 2546F180PS050                                     | Injection | 20 µl                    |
| Dimension   | 200 x 4.6 mm                                      | Sample    | Galvanic bath, technical |
| Eluent      | 10mM KH <sub>2</sub> PO <sub>4</sub> /            |           | 1 : iodate               |
|             | 2 mM TBA-H <sub>3</sub> PO <sub>4</sub> , pH 1.95 |           |                          |
| Flow        | 1 ml/min  |           |                          |
| Detection   | UV 196 nm   |           |                          |
| Temperature | 25 °C   |           |                          |



## Proteins II ProntoSIL 300-3-C4


|             |                                     |             |                       |
|-------------|-------------------------------------|-------------|-----------------------|
| Part number | 2003K040PS030                       | Detection   | UV 220 nm             |
| Dimension   | 200 x 3.0 mm                        | Temperature | 30 °C                 |
| Eluent      | A: H <sub>2</sub> O/0.1% TFA        | Injection   | 5 µl                  |
|             | B: ACN/H <sub>2</sub> O, 95/5 (v/v) |             |                       |
| Gradient    | 10-40% B, 0-218 s;                  | Sample      | 1 : Cytochrome c      |
|             | 40-47% B, 219-512 s;                |             | 2 : Fibrinogen, human |
|             | 47-100% B, 513-612 s                |             | 3 : Chymotrypsinogen  |
|             |                                     |             | 4 : β-Lactoglobulin   |
| Flow        | 0.6 ml/min                          |             | 5 : Thyroglobulin     |






## Peroxide ProntoSIL 120-3-C18 H

|             |   |               |                                     |
|-------------|---|---------------|-------------------------------------|
| Part number | 0603F185PS030   | Detection     | Coulometric                         |
| Dimension   | 53 x 3.0 mm   | Temperature   | 25 °C                               |
| Eluent      | A: 35mM citric acid/<br>85mM sodium acetate/<br>0.5mM heptane<br>sulphonate; pH 4.3 | Injection     | 5 µl                                |
|             | B: MeCH<br>A/B: 99.7/0.3 (w/w)  | Concentration | 500 ppb                             |
| Flow        | 0.5 ml/min  | Sample        | 1: Perbenzoic acid t-butyl<br>ester |




## 5-Aminolevulinic Acid ProntoSIL 120-3-C18 H

|             |  |               |                          |
|-------------|--|---------------|--------------------------|
| Part number | 1204F185PS030  | Temperature   | 26 °C                    |
| Dimension   | 125 x 4.6 mm   | Injection     | 10 µl                    |
| Eluent      | H <sub>2</sub> O/MeOH/THF, 90/6/4 (v/v)                | Concentration | 100 ppb                  |
|             | 0.5 mM octane sulphonic<br>acid (sodium salt); pH 3.26 | Sample        | 1: 5-Aminolevulinic acid |
| Flow        | 1 ml/min   |               |                          |
| Detection   | Conductivity; FS: 5 µS                                 |               |                          |




## Alcohols I ProntoSIL 120-5-C18 AQ

|             |                                      |               |                |
|-------------|--------------------------------------|---------------|----------------|
| Part number | 1546F184PS050                        | Injection     | 10 µl          |
| Dimension   | 150 x 4.6 mm                         | Concentration | 10 ppb each    |
| Eluent      | H <sub>2</sub> O/MeOH, 96/4 (v/v)    | Sample        | 1: Ethanol     |
|             | 0.8 mM Caps; pH 2.5                  |               | 2: Isopropanol |
| Flow        | 1.1 ml/min                           |               | 3: n-Propanol  |
| Detection   | Indirect Conductivity;<br>FS: 500 µS |               |                |
| Temperature | 26 °C                                |               |                |



## Dihydroxybenzoic Acids ProntoSIL 120-3-CN


|             |   |               |  |
|-------------|---|---------------|--|
| Part number | 2003F200PS030                                       | Temperature   | 20 °C  |
| Dimension   | 200 x 3 mm  | Injection     | 5 µl   |
| Eluent      | 30mM Sodium acetate/<br>30mM Sodium Citrate; pH 4.6 | Concentration | 1.8 mMol each  |
|             | 0.5 ml/min  | Sample        | 1: 2,3-Dihydroxybenzoic Acid<br>2: 2,5-Dihydroxybenzoic Acid<br>3: Salicylic Acid (1,2-DHBA) |
| Flow        | 0.5 ml/min  |               |  |
| Detection   | UV 220 nm   |               |  |






## Bitter ingredients of hop I ProntoSIL 120-5-C18 SH

|             |   |           |  |
|-------------|---|-----------|--|
| Part number | 0604F180PS050   | Injection | 10 µl  |
| Dimension   | 53 x 4 mm   | Sample    | Hop CO <sub>2</sub> -extract   |
| Eluent      | H <sub>2</sub> O/MeOH 25/75 (v/v)<br>+ 50 mM H <sub>3</sub> PO <sub>4</sub> |           | 1 : not identified<br>2 : Cohumulon<br>3 : Deoxyhumulon<br>4 : n-plus Adhumulon<br>5 : Colupulon<br>6 : n-plus Colupulon |
| Flow        | 1 ml/min  |           |  |
| Detection   | UV 314 nm   |           |  |
| Temperature | 21 °C   |           |  |




## Bitter ingredients of hop II ProntoSIL 120-5-C18 SH

|             |   |           |  |
|-------------|---|-----------|--|
| Part number | 2504F080PS050   | Injection | 10 µl  |
| Dimension   | 250 x 4 mm  | Sample    | Hop CO <sub>2</sub> -extract   |
| Eluent      | H <sub>2</sub> O/MeOH 20/80 (v/v)<br>+ 50 mM H <sub>3</sub> PO <sub>4</sub> |           | 1 : not identified<br>2 : Cohumulon<br>3 : Deoxyhumulon<br>4 : n-plus Adhumulon<br>5 : Colupulon<br>6 : n-plus Colupulon |
| Flow        | 1 ml/min  |           |  |
| Detection   | UV 314 nm   |           |  |
| Temperature | 21 °C   |           |  |




## Melatonin ProntoSIL 120-3-C18 H

|             |                 |               |               |
|-------------|-----------------|---------------|---------------|
| Part number | 1003F185PS030   | Injection     | 10 µl         |
| Dimension   | 100 x 3 mm      | concentration | 10 ppb        |
| Eluent      | ESA CAT-A Phase | Sample        | 1 : Melatonin |
| Flow        | 0.5 ml/min      |               |               |
| Detection   | Coulometric     |               |               |
| Temperature | 25 °C           |               |               |



## Aldehydes and Ketones I: DNPH Derivatives in Galvanic Bath ProntoSIL 120-5-C18 H


|             |  |               |   |
|-------------|--|---------------|---|
| Part number | 2546F185PS050  | Injection     | 50 µl   |
| Dimension   | 250 x 4.6 mm   | concentration | 10 ppb  |
| Eluent      | A : H <sub>2</sub> O; B: ACN                                 | Sample        | 1 : DNPH<br>2 : Formaldehyde<br>3 : Acetaldehyde<br>4 : Acetone<br>5 : Propanal<br>6 : Pentanal |
| Gradient    | 40% B, 0-60 s;<br>40-60% B, 61-390 s;<br>60-80% B, 391-480 s |               |   |
| Flow        | 1 ml/min   |               |   |
| Detection   | UV 360 nm  |               |   |
| Temperature | 30 °C  |               |   |






## Fast Separation of Proteins ProntoSIL 300-3-C18 H

|             |  |             |  |
|-------------|--|-------------|--|
| Part number | 0746K18PS030   | Detection   | UV 220 nm  |
| Dimension   | 75 x 4.6 mm  | Temperature | 30 °C  |
| Eluent      | A: H <sub>2</sub> O/0.1% TFA<br>B: ACN/H <sub>2</sub> O 95/5 (v/v)/0.1% TFA  | Injection   | 5 µl   |
| Gradient    | a) 29-48% B, 0-25 s;<br>48-100% B, 25-65 s<br>b) 29-48% B, 0-49 s;<br>48-100% B, 49-130 s<br>c) 29-100% B, 0-360 s | Sample      | 1: Ribouclease A<br>2: Insulin, bovine<br>3: Lysozyme<br>4: BSA<br>5: Myoglobin<br>6: Ovalbumin<br>7: not identified |
| Flow        | a) 4ml/min; b) 2ml/min;<br>c) 1ml/min  |             |  |



## Fast Separation of Peptides ProntoSIL 300-3-C18 H

|             |  |             |   |
|-------------|--|-------------|---|
| Part number | 0746K185PS030  | Temperature | 30 °C   |
| Dimension   | 75 x 4.6 mm  | Injection   | 5 µl  |
| Eluent      | A: H <sub>2</sub> O/0.1% TFA<br>B: ACN/H <sub>2</sub> O 70/30 (v/v),<br>0.1% TFA | Sample      | 1: Oxytocin<br>2: Bradykinin<br>3: Angiotensin I<br>4: Eledoisin<br>5: Neurotensin<br>6: Angiotensin II |
| Gradient    | 35-80% B, 0-480 s  |             |   |
| Flow        | 1 ml/min   |             |   |
| Detection   | UV 220 nm  |             |   |






# ProntoSIL HPLC Column Application - 2

|   | <b>Page</b> |
|---|-------------|
| • Ketones I .....   | 33          |
| • Organic Acids III   |             |
| • Coconut Fatty Acid - Ethoxylates 10 EO                      |             |
| • Castor oil Ethoxylates 20 EO                                |             |
| • Soybean oil - Ethoxylates 20 EO .....                       | 34          |
| • Calcium-DL-2-Hydroxy-4-methylthiobutyrate in Animal Food    |             |
| • Ubiquinols and Ubiquinones                                  |             |
| • Aldehydes and Ketones as DNPH Derivatives                   |             |
| • Tryptic digest of Casein Peptone                            |             |
| • Aflatoxines with KBr Derivatization Peptone .....           | 35          |
| • Amino Acids with FMOC                                       |             |
| • Chinin and Saccharin in Beverages                           |             |
| • Flavonoids in Red Wine                                      |             |
| • Analysis of Nitro-PAH with Fluorescence Detection (I) ..... | 36          |
| • Analysis of Nitro-PAH with Fluorescence Detection (II)      |             |




## Ketones I ProntoSIL 120-3-C18 AQ

|             |                            |        |                      |
|-------------|----------------------------|--------|----------------------|
| Part number | 0203F184PS030              | Sample | 1: Tetramethylketone |
| Dimension   | 25 x 3 mm                  |        | 2: 2-Pentanone       |
| Eluent      | A: H <sub>2</sub> O B: ACN |        | 3: 2-Hexanone        |
| Gradient    | 10-100% B, 0-18 min        |        | 4: 2-Hexanone        |
| Flow        | 0.35 ml/min                |        | 5: 2-Octanone        |
| Detection   | UV 270 nm                  |        | 6: 2-Nonanone        |
| Temperature | 30 °C                      |        | 7: 2-Decanone        |
| Injection   | 2 µl                       |        | 8: 2-Undecanone      |
|             |                            |        | 9: 2-Dodecanone      |
|             |                            |        | 10: 2-Tridecanone    |
|             |                            |        | 11: 2-Hexadecanone   |




## Organic Acids III ProntoSIL 120-3-C18 AQ

|               |  |        |                          |
|---------------|--|--------|--------------------------|
| Part number   | 3003F184PS030  | Sample | 1: oxalic acid 150 ppm   |
| Dimension     | 300 x 3 mm   |        | 2: not identified        |
| Eluent        | 5 mM Li <sub>2</sub> SO <sub>4</sub> / H <sub>2</sub> SO <sub>4</sub><br>pH 2.81 |        | 3: tartaric acid 150 ppm |
| Flow          | 0.56 ml/min  |        | 4: quinic acid           |
| Detection     | RI   |        | 5: malic acid            |
| Temperature   | 20 °C  |        | 6: malonic acid          |
| Injection     | 10 µl  |        | 7: shikimic acid         |
| Concentration | 300 ppm each, except 1,3   |        | 8: lactic acid           |
|               |  |        | 9: ascorbic acid         |
|               |  |        | 10: acetic acid          |
|               |  |        | 11: citric acid          |
|               |  |        | 12: fumaric acid         |
|               |  |        | 13: succinic acid        |




## Coconat Fatty Acid Ethoxylates 10 EO ProtoSIL 120-5-NH<sub>2</sub>

|             |   |             |   |
|-------------|---|-------------|---|
| Part number | 2546F190PS050   | Flow        | 1 ml/min  |
| Dimension   | 250 x 4.6 mm  | Detection   | Evap. Light Scattering<br>Detector (ELSD), model<br>DDI 31; PMT: 300,<br>T: 42 °C, p(air): 1bar |
| Eluent      | A: n-hexane<br>B: acetone   | Temperature | 30 °C   |
| Gradient    | 3-20% B, 0-770 s;<br>20-38% B, 771-2600 s;<br>38-70% B, 2601-3600 s<br>70% B, 3600-4900 s | Injection   | 50 µl   |



## Castor oil Ethoxylates 20 EO ProtoSIL 120-5-NH<sub>2</sub>


|             |   |             |   |
|-------------|---|-------------|---|
| Part number | 2546F190PS050   | Detection   | Evap. Light Scattering<br>Detector (ELSD), model<br>DDI 31; PMT: 500,<br>T: 42 °C, p(air): 1bar |
| Dimension   | 250 x 4.6 mm  | Temperature | 30 °C   |
| Eluent      | A: n-hexane<br>B: acetone                                       | Injection   | 50 µl   |
| Gradient    | 0-5% B, 0-1200 s;<br>5-40% B, 2601-3600 s<br>40% B, 2601-3600 s |             |   |
| Flow        | 1 ml/min  |             |   |






**Soybean oil  
Ethoxylates 20 EO  
ProtoSIL 120-5-NH<sub>2</sub>**

|             |   |                        |                        |
|-------------|---|------------------------|------------------------|
| Part number | 2546F190PS050   | Flow                   | 1 ml/min               |
| Dimension   | 250 x 4.6 mm  | Detection              | Evap. Light Scattering |
| Eluent      | A: n-hexane<br>B: acetone   | Detector (ELSD), model | DDI 31; PMT: 300,      |
| Gradient    | 3-20% B, 0-770 s;<br>20-38% B, 771-2600 s;<br>38-70% B, 2601-3600 s<br>70% B, 2601-3600 s | T: 42 °C, p(air): 1bar |                        |
|             |   | Temperature            | 30 °C                  |
|             |   | Injection              | 50 µl                  |




**Calcium-DL-2-Hydroxy-4-  
methylthiobutyrate in  
Animal Food  
ProtoSIL 120-5-NH<sub>2</sub>**

|             |   |             |             |
|-------------|---|-------------|-------------|
| Part number | 1003F185PS030   | Flow        | 0.55 ml/min |
| Dimension   | 100 x 3 mm  | Detection   | Coulometric |
| Eluent      | A: 75 mM NaH <sub>2</sub> PO <sub>4</sub> /<br>0.2 mM heptanesulfonic<br>acid (sodium salt), pH 2.9 | Temperature | 25 °C       |
| Gradient    | B: ACN; A/B 98/2 (v/v)  | Injection   | 10 µl       |



**Aldehydes and Ketones  
as DNPH Derivatives  
ProntoSIL 120-5-C18 H**

|  |               |        |  |
|--|---------------|--------|--|
| Part number  | 2504F185PS050 | Sample | 1: Nonanal (internal standard)<br>2: 2,4-Dinitrophenylhydrazine<br>3: Formaldehyde<br>4: Acetaldehyde<br>5: Acetone<br>6: Acrolein<br>7: Propanal<br>8: Crotonic Aldehyde<br>9: n-Butanal<br>10: Benzaldehyde<br>11: Pentanal<br>12: 2-Methylbenzaldehyde<br>13: 4-Methylbenzaldehyde<br>14: 3,3-Dimethylbutanone<br>15: Hexanal<br>16: 2,4-Dimethylbenzaldehyde<br>17: 2,4-dimethylpentanone<br>18: Heptanal<br>19: Octanal |
| Dimension  | 250 x 4.0 mm  |        |  |
| Eluent/Gradient (final compositions of linear steps are given (v/v/v ratios)): |               |        |  |
| 0-5 min isocratic:   |               |        |  |
| THE/ACN/H <sub>2</sub> O (16/25/59)  |               |        |  |
| 10 min. The/ACN/H <sub>2</sub> O (10/40/50)                                    |               |        |  |
| 20 min. The/ACN/H <sub>2</sub> O (0/80/40)                                     |               |        |  |
| 30 min. The/ACN/H <sub>2</sub> O (0/80/20)                                     |               |        |  |
| 40 min. The/ACN/H <sub>2</sub> O (0/100/0)                                     |               |        |  |
| 45 min. The/ACN/H <sub>2</sub> O (16/25/59)                                    |               |        |  |
| Flow   | 1.3 ml/min    |        |  |
| Detection  | UV 360 nm     |        |  |
| Temperature  | 30 °C         |        |  |
| Injection  | 10 µl         |        |  |




by courtesy of Dr. Turok, Landesuniversität Brandenburg, Germany



**Tryptic digest of Casein  
Peptone  
ProntoSIL 120-3-C18 AQ**


|             |   |             |                                     |
|-------------|---|-------------|-------------------------------------|
| Part number | 0746K180PS030   | Flow        | 1 ml/min                            |
| Dimension   | 75 x 4.6 mm   | Detection   | UV 220 nm                           |
| Eluent      | A: H <sub>2</sub> O/0.1% TFA<br>B: ACN/ H <sub>2</sub> O 70/30 (v/v),<br>0.1% TFA | Temperature | 30 °C                               |
| Gradient    | 0-40% B, 0-1400 s<br>40-0% B, 1401-1481 s   | Injection   | 5 µl                                |
|             |   | Sample      | Tryptic digest of<br>casein peptone |






## Aflatoxines with KBr Derivatization ProntoSIL 120-5-C18 AQ

|             |   |             |  |
|-------------|---|-------------|--|
| Part number | 2546F184PS050   | Temperature | 30 °C  |
| Dimension   | 250 x 4.6 mm  | Injection   | 5 µl   |
| Eluent      | ACN/MeOH/H <sub>2</sub> O<br>20/20/60 (v/v/v)           | Sample      | Aflatoxin-Standard<br>G2: 60 fg/ml<br>G1: 100 fg/ml<br>B2: 50 fg/ml<br>B1: 90 fg/ml<br>Injection: 100 µl |
| Flow        | 1.0 ml/min  |             |  |
| Detection   | Fluorescence Detector Model<br>8470 Ex: 365nm Em: 429nm |             |  |




## Amino Acids with FMOc ProntoSIL-AA-FMOc 5.0µm

|             |   |             |  |
|-------------|---|-------------|--|
| Part number | 2504F470PS050   | Detection   | Fluorescence<br>Ex: 263 nm Em: 313 nm                          |
| Dimension   | 250 x 4.0 mm  | Temperature | 30 °C  |
| Eluent      | A: 50mM CH <sub>3</sub> COONa pH 4.2, 5% THF<br>B: ACN  | Injection   | 10 µl  |
| Flow        | 1.0 ml/min  | Sample      | SIGMA standard AA-S-18 + Asn, Gln, IS (β-tAla)<br>FMOc labeled |
| Gradient    | 20-36% B, 0-20 min<br>36-38% B, 20-28 min<br>38% B, 28-34 min<br>38-59% B, 34-51 min<br>59-67% B, 51-56 min |             |  |




## Chinin and Saccharin in Beverages ProntoSIL-120-5-C18 H 5.0 µm

|             |   |             |  |
|-------------|---|-------------|--|
| Part number | 1204F185PS050   | Temperature | 25 °C  |
| Dimension   | 125 x 4.0 mm  | Injection   | 20 µl  |
| Eluent      | MeOH/H <sub>2</sub> O 70/30 (v/v)<br>+ 5ml/l H <sub>3</sub> PO <sub>4</sub> | Sample      | Bitter Lemon Light<br>1: Quinine<br>2: Saccharin |
| Flow        | 1.0 ml/min  |             |  |
| Detection   | UV 210 nm   |             |  |



## Flavonoids in Red Wine ProntoSIL 120-5-C8 ace-EPS 5.0µm


|             |  |             |   |
|-------------|--|-------------|---|
| Part number | 1204F08APS050  | Detection   | UV 280 nm (Stilbenes)<br>370 nm (Flavonoides)   |
| Dimension   | 125 x 4.0 mm   | Temperature | 40 °C   |
| Eluent      | A: 10 mM NaH <sub>2</sub> PO <sub>4</sub> pH 3.0<br>B: ACN   | Injection   | 25 µl   |
| Flow        | 2.5 ml/min   | Sample      | Red Wine<br>1: Myricetin<br>2: cis- Resveratrol<br>3: trans- Resveratrol<br>4: Quercetin<br>5: Caempferol |
| Gradient    | 30% B, 0-2 min<br>30-70% B, 2-4 min<br>70% B, 4-4.4 min<br>70-30% B, 4.4-4.5min<br>30% B, 4.5-5.5min |             |   |






## Analysis of Nitro-PAH with Fluorescence Detection (I) ProntoSIL 120-3-Phenyl


|             |   |             |                |
|-------------|---|-------------|----------------|
| Part number | 1204F050PS030   | Temperature | 25 °C          |
| Dimension   | 125 x 4.0 mm  | Injection   | 20 µl          |
| Eluent      | A: MeOH/H <sub>2</sub> O 60/40 (v/v)<br>B: MeOH 0 % B hold for<br>13 min, 0 – 100 % B in 29 min | Sample      | Nitro-PAH Mix1 |
| Flow        | 0.8 ml/min  |             |                |
| Detection   | Fluorescence Detector Model 8470<br>wavelength program  |             |                |



## Analysis of Nitro-PAH with Fluorescence Detection (II) ProntoSIL 120-3-Phenyl

|             |   |             |                |
|-------------|---|-------------|----------------|
| Part number | 1204F050PS030   | Temperature | 25 °C          |
| Dimension   | 125 x 4.0 mm  | Injection   | 20 µl          |
| Eluent      | A: MeOH/H <sub>2</sub> O 60/40 (v/v)<br>B: MeOH 0 % B hold for<br>13 min, 0 – 100 % B in 29 min | Sample      | Nitro-PAH Mix2 |
| Flow        | 0.8 ml/min  |             |                |
| Detection   | Fluorescence Detector Model 8470<br>wavelength program  |             |                |






# ProntoSIL HPLC Column Application - 3

|  | <b>Page</b> |
|--|-------------|
| • Carbohydrates, Organic Acids and Alcohols in Wine .....                          | 38          |
| • Dimethylsiloxane with Evaporative Light Scattering Detection                     |             |
| • Fast Analysis of a pharmaceutical Teststandard<br>with Multiwavelength Detection |             |
| • Fast Analysis of a Cold Medicine II<br>with Multiwavelength Detection            |             |
| • Standardtestmixture with Multiwavelength Detection .....                         | 39          |
| • Fast Analysis of Parabens II with Multiwavelength Detection                      |             |
| • Fat soluble Vitamins II with Multiwavelength Detection                           |             |
| • Aromatic Hydrocarbons according to prEN 12916                                    |             |
| • Diesel .....   | 40          |
| • Acetaldehyde in Mineralwater   |             |
| • Ergosterol in Grass  |             |
| • Sulfa Drugs  |             |




## Carbohydrates, Organic Acids and Alcohols in Wine Carbohydrate H<sup>+</sup>

Part number 00253776  
 Dimension 300 x 7.8 mm  
 Eluent 1.25 mM H<sub>2</sub>SO<sub>4</sub>  
 Flow 0.6 ml/min  
 Detection RI  
 Temperature 45 °C  
 Injection 50 µl  
 Sample Wine Complet Standard




## Dimethylsiloxane with Evaporative Light Scattering Detection ProntoSIL 120-5-C18 SH

Part number 2546F180PS050  
 Dimension 250 x 4.6 mm  
 Eluent A: ACN/Acetone 70/30 (v/v)  
 B: Ethyl acetate  
 5-70% B in 50 min, 70% B hold for 15 min  
 Flow 0.9 ml/min  
 Detection Evap. Light Scattering detector (DDL31)  
 PMT: 500, Temp: 40 °C  
 Pressure (air): 1.5 bar  
 Temperature 25 °C  
 Sample Silicone-Standard (DMS-T11)  
 Injection: 10 µl




## Fast Analysis of a pharmaceutical Teststandard with Multiwavelength Detection ProntoSIL-120-5-C18 AQ

Part number 0446F184PS050  
 Dimension 33 x 4.6 mm  
 Eluent A: H<sub>2</sub>O + 0.1% TFA; B: ACN  
 Gradient 20-30% B in 60 sec  
 Flow 1.2 ml/min  
 Detection Multiwavelength Detector (DAD-3L)  
 Wavelength 230, 261 and 274 nm; Ref. at 350nm  
 Temperature 25 °C  
 Injection 1 µl  
 Sample 1 Ascorbic acid  
 2 Paracetamol  
 3 Caffeine  
 4 Acetylsalicylic acid



## Fast Analysis of a Cold Medicine II with Multiwavelength Detection ProntoSIL-120-5-C18 AQ


Part number 0446F184PS050  
 Dimension 33 x 4.6 mm  
 Eluent A: H<sub>2</sub>O + 0.1% TFA; B: ACN  
 Gradient 20-30% B in 30 sec  
 Flow 1.2 ml/min  
 Detection Multiwavelength Detector (DAD-3L)  
 Wavelength 230, 254 and 274 nm; Ref. at 350nm  
 Temperature 25 °C  
 Injection 1 µl  
 Sample 1 Paracetamol  
 3 Caffeine  
 4 Propylphenazone






### Standardtestmixture with Multiwavelength Detection ProntoSIL 120-3-C18 SH

|             |                                      |             |            |
|-------------|--------------------------------------|-------------|------------|
| Part number | 0604F180PS030                        | Temperature | 25 °C      |
| Dimension   | 53 x 4.0 mm                          | Injection   | 1 µl       |
| Eluent      | H2O/ ACN; 40/60                      | Sample      | 1 Acetone  |
| Flow        | 1.0 ml/min                           |             | 2 Benzene  |
| Detection   | Multiwavelength Detector<br>(DAD-3L) |             | 3 Toluene  |
|             |                                      |             | 4 m-Xylene |
| Wavelengths | 254, 262 and 273 nm<br>Ref. at 350nm |             |            |




### Fast Analysis of Parabens II with Multiwavelength Detection ProntoSIL 120-5-C18 AQ

|             |                                      |             |                                       |
|-------------|--------------------------------------|-------------|---------------------------------------|
| Part number | 0446F184PS050                        | Wavelengths | 210, 230 and 254 nm<br>Ref. at 350 nm |
| Dimension   | 33 x 4.6 mm                          | Temperature | 25 °C                                 |
| Eluent      | A: H2O B: ACN                        | Injection   | 2 µl                                  |
| Gradient    | 20-30% B in 30 sec                   | Sample      | 1 Uracil                              |
| Flow        | 1.0 ml/min                           |             | 2 Methylparabene                      |
| Detection   | Multiwavelength Detector<br>(DAD-3L) |             | 3 Ethylparabene                       |
|             |                                      |             | 4 Propylparabene                      |




### Fat soluble Vitamins II with Multiwavelength Detection ProntoSIL 120-3-C18 SH

|             |                                      |           |                           |
|-------------|--------------------------------------|-----------|---------------------------|
| Part number | 2503F180PS030                        | Injection | 5 µl                      |
| Dimension   | 250 x 3.0 mm                         | Sample    | 1: Vitamin A              |
| Eluent      | MeOH                                 |           | 2: Vitamin A acetate      |
| Flow        | 1 ml/min                             |           | 3: Vitamin D <sub>2</sub> |
| Detection   | Multiwavelength Detector<br>(DAD-3L) |           | 4: Vitamin D <sub>3</sub> |
|             |                                      |           | 5: Vitamin E              |
|             |                                      |           | 6: Vitamin K <sub>1</sub> |
| Wavelengths | 225, 264, 325 nm                     |           |                           |
| Temperature | 25 °C                                |           |                           |



### Aromatic Hydrocarbons according to prEN 12916 ProntoSIL ENVIRO PTL


|              |               |        |   |
|--------------|---------------|--------|---|
| Part number  | 2546C450PS050 | Sample | Calibration Standard B -<br>EUOPEAN STANDARD<br>prEN 12916 (June 1997)  |
| Dimension    | 250 x 4.6 mm  |        | 1: Cyclohexane (2.0g/100ml)   |
| Eluent       | Heptane       |        | 2: o-Xylene (1.0g/100ml)  |
| Flow         | 1.0 ml/min    |        | 3: 1-Methylnaphthalene<br>(1.0g/100ml)  |
| Backflush on | 16 min        |        | 4: Phenanthrene (0.2g/100ml)  |
| Detection    | RI            |        | 7. 22 between Cyclohexane and<br>o- Xylene<br>(the resolution between Cyclohexane and o-Xylene should not be less than 5.0) |
| Temperature  | 25 °C         |        |   |
| Injection    | 10 µl         |        |   |






## Diesel ProntoSIL ENVIRO PTL

|              |               |        |                        |
|--------------|---------------|--------|------------------------|
| Part number  | 2546C450PS050 | Sample | 1: Cyclohexane         |
| Dimension    | 250 x 4.6 mm  |        | 2: Monocycl. Aromatics |
| Eluent       | Heptane       |        | 3: Dicycl. Aromatics   |
| Flow         | 1.0 ml/min    |        | 4: Dicycl. Aromatics   |
| backflush on | 20 min        |        | 5: Polycycl. Aromatics |
| Detection    | RI            |        |                        |
| Temperature  | 25 °C         |        |                        |
| Injection    | 10 µl         |        |                        |




## Acetaldehyde in Mineralwater ProntoSIL 120-3-C18 ACE- EPS

|             |                                   |        |   |
|-------------|-----------------------------------|--------|---|
| Part number | 1504F18APS030                     | Sample | Mineralwater after Derivatisation with DNPH |
| Dimension   | 150 x 4.0 mm                      |        | 1 DNPH                                      |
| Eluent      | H <sub>2</sub> O/ ACN 50/50 (v/v) |        | 5 Acetaldehyde. 6ppb                        |
| Flow        | 1.0 ml/min                        |        |   |
| Detection   | UV 360 nm                         |        |   |
| Temperature | 25 °C                             |        |   |
| Injection   | 500 µl                            |        |   |




## Ergosterol in Grass ProntoSIL 120-5-C18 H

|             |                                  |        |                                      |
|-------------|----------------------------------|--------|--------------------------------------|
| Part number | 1204F18SPS050                    | Sample | 1 Ergosterol                         |
| Dimension   | 125 x 4.0 mm                     |        | Extracted according to VDLUFA (1993) |
| Eluent      | H <sub>2</sub> O/MeOH 2/98 (v/v) |        |                                      |
| Flow        | 0.75 ml/min                      |        |                                      |
| Detection   | UV 282 nm                        |        |                                      |
| Temperature | 25 °C                            |        |                                      |
| Injection   | 20 µl                            |        |                                      |



## Sulfa Drugs ProntoSIL 120-5-C18 ACE - EPS

|             |   |           |                  |
|-------------|---|-----------|------------------|
| Part number | 1546F18APS050                                     | Injection | 3 µl             |
| Dimension   | 150 x 4.6 mm                                      | Sample    | 1 Sulfanilamide  |
| Eluent      | MeOH/H <sub>2</sub> O/Acetic Acid 20/79/1 (v/v/v) |           | 2 Sulfadiazine   |
| Flow        | 1 ml/min  |           | 3 Sulfathiazole  |
| Detection   | UV 254 nm   |           | 4 Sulfamerazine  |
| Temperature | 25 °C   |           | 5 Sulfamethazine |







# ProntoSIL HPLC Column Application - 4

Page



- Sugaralcohols, Mono-and Disaccharides with BioQuant Carbohydrate Inositol
  - Alcohol, Sugaralcohols, Mono-and Disaccharides with BioQuant Carbohydrate
- .....40





# ProntoSIL HPLC Column Application - 5

- High Speed Separation of Catecholamines





# ProntoSIL HPLC Column Application - 6


Page

- Sugaralcohols, Mono-and Disaccharides with BioQuant Carbohydrate . . 44
- Sugars in Applejuice with BioQuant Carbohydrate
- Sugars in Applejuice
- Sugars in Orangejuice (A) and Diet Orangejuice (B)
- Sugars in Shampoo ..... 45
- Maltooligosaccacharides from Cornsyrup
- Monosaccharides of Glycoproteins
- Plasma Catecholamines
- Catecholamine - 6 - Standard ..... 46
- Catecholamine - 8 - Standard




## Sugaralcohols, Mono- and Disaccharides with BioQuant Carbohydrate

|             |                       |        |                       |
|-------------|-----------------------|--------|-----------------------|
| Column      | PAD Column            | Sample | 1 Inositol            |
| Dimension   | 250 x 4.6 mm PEEK     |        | 2 Arabitol            |
| Eluent      | 100mM NaOH            |        | 3 Sorbitol            |
| Flow        | 0.9 ml/min            |        | 4 Fucose              |
| Detection   | BioQuant Carbohydrate |        | 5 Arabinose           |
| Temperature | 28 °C                 |        | 6 Glucose             |
| Injection   | 20 µl                 |        | 7 Xylose              |
|             |                       |        | 8 Fructose            |
|             |                       |        | 9 Lactose             |
|             |                       |        | 10 Sucrose (1-50mg/l) |




## Sugars in Applejuice with BioQuant Carbohydrate

|             |                         |        |            |
|-------------|-------------------------|--------|------------|
| Column      | PAD Column              | Sample | 1 Sorbitol |
| Dimension   | 250 x 4.6 mm PEEK       |        | 2 Glucose  |
| Eluent      | 100 mM NaOH             |        | 3 Fructose |
| Flow        | 0.9 ml/min              |        | 4 Sucrose  |
| Detection   | BioQuant carbohydrate   |        |            |
| Temperature | 28 °C                   |        |            |
| Injection   | 20 µl (diluted 1: 1000) |        |            |




## Sugars in Applejuice

|             |                       |        |            |
|-------------|-----------------------|--------|------------|
| Column      | hamilton RCX-30       | Sample | 1 Sorbitol |
| Dimension   | 150 x 4.6 mm PEEK     |        | 2 Glucose  |
| Eluent      | 60 mM NaOH            |        | 3 Fructose |
| Flow        | 2.00 ml/min           |        | 4 Sucrose  |
| Detection   | BioQuant Carbohydrate |        |            |
| Temperature | 32 °C                 |        |            |
| Injection   | 1 µl (diluted 1:1000) |        |            |



## Sugars in Oragnejuice (A) and Diet Oragnejuice (B)


|             |                       |           |                        |
|-------------|-----------------------|-----------|------------------------|
| Column      | Hamilton RCX-30       | Injection | 5 µl (diluted 1: 1000) |
| Dimension   | 150 x 4.6 mm PEEK     | Sample    | 1 Glucose              |
| Eluent      | 60 mM NaOH            |           | 2 Fructose             |
| Flow        | 2.00 ml/min           |           | 3 Sucrose              |
| Detection   | BioQuant carbohydrate |           |                        |
| Temperature | 32 °C                 |           |                        |






## Sugars In Shampoo

|             |                       |           |            |
|-------------|-----------------------|-----------|------------|
| Column      | Hamilton RCX-30       | Injection | 1 µl       |
| Dimension   | 150 x 4.6 mm PEEK     | Sample    | 3 Sorbitol |
| Eluent      | 60 mM NaOH            |           | 5 Glucose  |
| Flow        | 2.00 ml/min           |           | 6 Fructose |
| Detection   | BioQuant Carbohydrate |           | 7 Sucrose  |
| Temperature | 32 °C                 |           |            |




## Maltooligosaccharides from Cornsyrup

|           |   |             |                       |
|-----------|---|-------------|-----------------------|
| Column    | Hamilton RCX  | Flow        | 2.00 ml/min           |
| Dimension | 150 x 4.6 mm PEEK   | Detection   | BioQuant carbohydrate |
| Eluent    | Eluent A: 60 mM NaOH<br>Eluent B: 60 mM NaOH:<br>500 mM Na-acetat | Temperature | 32 °C                 |
| Gradient  | 10-60% B in 0-20 min  | Injection   | 5 µl                  |
|           |   | Sample      | G4- G10               |




## Monosaccharides of Glycoproteins

|             |                       |           |                 |
|-------------|-----------------------|-----------|-----------------|
| Column      | Dionex CarboPAC PA1   | Injection | 5 µl            |
| Dimension   | 250 x 4.0 mm PEEK     | Sample    | 1 Galactosamine |
| Eluent      | 20 mM NaOH            |           | 2 Glucosamine   |
| Flow        | 1.00 ml/min           |           |                 |
| Detection   | BioQuant Carbohydrate |           |                 |
| Temperature | 28 °C                 |           |                 |



## Plasma Catacholamines


|             |                        |           |                |
|-------------|------------------------|-----------|----------------|
| Column      | ProntoSIL-120-3-C18 AQ | Injection | 5 µl           |
| Dimension   | 33 x 3.0 mm            | Sample    | 1 Noradrenalin |
| Eluent      | Pronto- Cat Phase      |           | 2 Adrenalin    |
| Flow        | 0.5 ml/min             |           | 3 DHBA         |
| Detection   | BioQuant carbohydrate  |           | 4 Dopamin      |
| Temperature | 30 °C                  |           |                |





## Catecholamine -6 -Standard

|             |                        |        |             |
|-------------|------------------------|--------|-------------|
| Column      | ProntoSIL-120-3-C18 AQ | Sample | 1 Adrenalin |
| Dimension   | 33 x 3.0 mm            |        | 2 DHBA      |
| Eluent      | Pronto- Cat Phase      |        | 3 DOPAC     |
| Flow        | 0.5 ml/min             |        | 4 5-HIAA    |
| Detection   | BioQuant Carbohydrate  |        | 5 HVA       |
| Temperature | 30 °C                  |        | 6 Serotonin |
| Injection   | 5 µl                   |        |             |



## Catecholamine -8 -Standard

|             |                             |           |                |
|-------------|-----------------------------|-----------|----------------|
| Column      | ProntoSIL-120-3-C18 AQ      | Injection | 5 µl           |
| Dimension   | 33 x 3.0 mm                 | Sample    | 1 Noradrenalin |
| Eluent      | Pronto- Cat Phase           |           | 2 Adrenalin    |
| Flow        | 0.5 ml/min                  |           | 3 DHBA         |
| Detection   | BioQuant CAT Amine, +500 mV |           | 4 DOPAC        |
| Temperature | 30 °C                       |           | 5 Dopamin      |
|             |                             |           | 6 5-HIAA       |
|             |                             |           | 7 HVA          |
|             |                             |           | 8 Serotonin    |

