


Pharmacopeia

Para-Hydroxybenzoate (JP)



Conditions


Column : Inertsil ODS-4(5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH
 B) 50 mM KH₂PO₄
 A/B = 13/7, v/v

Flow Rate : 1.3 mL/min
 Col. Temp. : 35 °C
 Detection : UV 272 nm
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LB097, LB098

Sample :

1. *p*-Hydroxybenzoic acid
 2. *p*-Hydroxybenzoic acid methyl ester
 3. *p*-Hydroxybenzoic acid ethyl ester
 4. *p*-Hydroxybenzoic acid *n*-propyl ester

Crospovidone (JP)



Conditions


Column : InertSustain C18 (5 μ m, 250 \times 4.0 mm I.D.)
 Guard Column : InertSustain C18 (5 μ m, 25 \times 4.0 mm I.D.)
 Eluent : A) CH₃CN
 B) H₂O
 A/B = 1/9, v/v

Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 235 nm
 Injection Vol. : 50 μ L
 Data Source : LC InertSearch No. LB167

Sample :

1. 1-Vinyl-2-pyrrolidone
 2. Vinyl acetate

Losartan potassium and Hydrochlorothiazide (USP)



Conditions


Column : Inertsil C8-3 (10 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) 10 mM KH₂PO₄ in H₂O (pH 2.5, H₃PO₄)
 A/B = 40/60, v/v

Flow Rate : 2.3 mL/min
 Col. Temp. : 35 °C
 Detection : UV 230 nm
 Injection Vol. : 20 μ L
 Data Source : LC InertSearch No. LB200

Sample :

1. Hydrochlorothiazide (14 mg/L)
 2. Losartan potassium (55 mg/L)

Carumonam Sodium (JP)



Conditions


Column : Inertsil ODS-2 (5 μ m, 150 \times 4.0 mm I.D.)
 Eluent : A) 0.1 g/L (NH₄)₂SO₄
 B) CH₃OH
 C) CH₃COOH
 A/B = 97/2/1, v/v/v

Flow Rate : 1.4 mL/min
 Col. Temp. : 25 °C
 Detection : UV 254 nm
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LA994

Sample :

1. Resorcinol (1.8 mg/L)
 2. Carumonam sodium (0.4 mg/L)

Glycyrrhetic acid in Glycyrrhiza (JP)



Conditions


Column : InertSustain C18(5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) 2.1 % CH₃COOH in H₂O
 A/B = 40/60, v/v

Flow Rate : 0.45 mL/min
 Col. Temp. : 20 °C
 Detection : UV 254 nm
 Injection Vol. : 20 μ L
 Data Source : LC InertSearch No. LB182

Sample :

1. Glycyrrhizic acid

Terbinafine Hydrochloride Cream (JP)



Conditions

Column : InertSustainSwift C18 (5 μ m, 125 \times 4.0 mm I.D.)
 Eluent : A) CH₃CN
 B) 0.45 % Tetramethylammonium hydroxide in H₂O (pH 8.0, 4 % H₃PO₄ in H₂O)
 C) THF
 A/B/C = 2/2/1, v/v/v

Flow Rate : 1.46 mL/min
 Col. Temp. : 25 °C
 Detection : UV 282 nm
 Injection Vol. : 10 μ L


Sample :

1. *p*-Terphenyl (17.5 mg/L)
 2. Terbinafine hydrochloride (200 mg/L)

Applications

Pharmaceuticals

17 Anti-Depressant Drugs




Conditions

Column : InertSustain C18 HP (3 μ m, 150 \times 2.1 mm I.D.)
 Eluent : A) 0.1 % HCOOH in CH₃CN
 B) 0.1 % HCOOH in H₂O
 A/B = 2/98 - 0.5 min - 40/60
 - 5.5 min - 40/60, v/v
 Flow Rate : 0.4 mL/min
 Col. Temp. : 40 °C
 Detection : LC/MS/MS (4000 QTRAP : ESI, Positive, MRM)
 Injection Vol. : 5 μ L
 Data Source : LC InertSearch No. LA908

Sample :
 1. Sulpiride 10. Imipramine
 2. Milnacipran 11. Nortriptyline
 3. Trazodone 12. Maprotiline
 4. Mianserin 13. Amitriptyline
 5. Amoxapine 14. Trimipramine
 6. Doxepin 15. Fluoxetine
 7. Paroxetine 16. Sertraline
 8. Desipramine 17. Clomipramine
 9. Fluvoxamine (100 ng/mL each)

18 Drugs




Conditions

Column : Inertsil ODS-3 (2 μ m, 50 \times 2.1 mm I.D.)
 Eluent : A) CH₃CN
 B) 0.05 % HCOOH in H₂O
 A/B = 5/95 - 5 min - 95/5, v/v
 Flow Rate : 500 μ L/min
 Col. Temp. : 40 °C
 Detection : LC/MS/MS (4000 QTRAP : ESI, Positive, MRM)
 Injection Vol. : 10 μ L

Sample :
 1. Ranitidine 10. Diphenhydramine
 2. Scopolamine 11. Doxepin
 3. Naltrexone 12. Amitriptyline
 4. Acetaminophen 13. Reserpine
 5. Theophylline 14. Isopropylantipyrene
 6. Metoprolol 15. Ketoprofen
 7. Caffeine 16. Warfarin
 8. Chlorpheniramine 17. Capsaicin
 9. Propranolol 18. Dihydrocapsaicin

Anti-inflammatory Drug




Conditions

Column : InertSustain Phenyl (2 μ m, 50 \times 2.1 mm I.D.)
 Eluent : A) CH₃OH/HCOOH = 100/0.05, v/v
 B) H₂O/HCOOH = 100/0.05, v/v
 A/B = 40/60, v/v
 Flow Rate : 0.6 mL/min
 Col. Temp. : 40 °C
 Detection : LC/MS/MS (4000 QTRAP:ESI, Positive, MRM)
 Injection Vol. : 5 μ L
 Data Source : LC InertSearch No. LB198

Sample :
 1. Hydrocortisone
 2. Prednisolone (0.1 mg/L each)

Histamine Antagonist




Conditions

Column : Inertsil ODS-4 (3 μ m, 150 \times 2.1 mm I.D.)
 Eluent : A) CH₃OH
 B) 2 mM CH₃COONH₄
 A/B = 40/60 - 10 min - 95/5
 - 2 min - 95/5, v/v
 Flow Rate : 0.2 mL/min
 Col. Temp. : 40 °C
 Detection : LC/MS/MS (4000 QTRAP : ESI, Positive, MRM)
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LA678

Sample :
 1. Chlorpheniramine 6. Diphenylpyraline
 2. Cinnarizin 7. Hydroxyzine
 3. Clemastine 8. Promethazine
 4. Difenedol 9. Triprolidine
 5. Diphenhydramine (0.1 mg/L each)

11 Drugs




Conditions

Column : Inertsil ODS-3 (3 μ m, 100 \times 2.1 mm I.D.)
 Eluent : A) CH₃CN
 B) 0.1 % HCOOH in H₂O
 A/B = 0/100 - 6 min - 10/90
 - 12 min - 100/0, v/v
 Flow Rate : 0.2 mL/min
 Col. Temp. : 40 °C
 Detection : LC/MS/MS (4000 QTRAP : ESI, Positive, MRM)
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LA610

Sample :
 1. Sulbutamol (7.5 μ g/L)
 2. Ranitidine (15.0 μ g/L)
 3. Acetaminophen (15.0 μ g/L)
 4. Theophylline (20.0 μ g/L)
 5. Metoprolol (22.5 μ g/L)
 6. Caffeine (15.0 μ g/L)
 7. Propranolol (12.5 μ g/L)
 8. Doxepin (5.0 μ g/L)
 9. Carbamazepine (5.0 μ g/L)
 10. Ketoprofen (15.0 μ g/L)
 11. Indomethacin (20.0 μ g/L)

Cold Medication




Conditions

Column : Inertsil ODS-4 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) 0.1 % H₃PO₄ in H₂O
 A/B = 3/97 - 20 min
 - 75/25, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 210 nm
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LA616

Sample :
 1. Maleic acid 6. Bromovalerylurea
 2. Acetaminophen 7. Apronalide
 3. Caffeine 8. Isopropylantipyrene
 4. Chlorpheniramine 9. Ibuprofen
 5. Ethenamide (50 mg/L each)

Biochemicals

Trypsin Digest of the Protein




Conditions

Column : MonoCap C18 High Resolution 2000 (2000 mm × 0.1 mm I.D.)
 Trap Column: MonoCap C18 Trap Column (50 mm × 0.075 mm I.D.)
 Eluent : A) 0.1 % HCOOH in CH₃CN , B) 0.1 % HCOOH in H₂O
 A/B=10/90 - 600 min - 45/55, v/v
 Flow Rate : 0.5 μL/min

Col. Temp. : ambient
 Detection : LC/MS/MS (JMS-T100-LC, JEOL, TIC, m/z 500-1500)
 Injection Vol. : 5 μL
 Sample : tryptic digest of proteins

Nucleoside and Nucleic acid base




Conditions

Column : InertSustain C18 (5 μm, 150 × 4.6 mm I.D.)
 Eluent : 0.1 M NH₂PO₄, 0.2 M NaClO₄ (pH 2.03, H₃PO₄)
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 260 nm
 Injection Vol. : 1 μL
 Data Source : LC InertSearch No. LA895

Sample :
 1. Cytosine
 2. Uracil
 3. Guanine
 4. Adenine
 5. Cytidine
 6. Uridine
 7. Thymine
 8. Adenosine
 9. Guanosine
 10. Thymidine

Nucleic acid base




Conditions

Column : Inertsil Amide (5 μm, 250 × 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) 10 mM HCOONH₄
 A/B = 90/10, v/v
 Flow Rate : 0.2 mL/min
 Col. Temp. : 35 °C
 Detection : UV 254 nm
 Injection Vol. : 1 μL
 Data Source : LC InertSearch No. LB001

Sample :
 1. Toluene
 2. Thymine
 3. Uracil
 4. Uridine
 5. Adenosine
 6. Cytosine
 7. Cytidine
 8. Guanosine
 (10 mg/L each)

Protein




Conditions

Column : Inertsil WP300 C18 (5 μm, 150 × 3.0 mm I.D.)
 Eluent : A) 0.1 % TFA in CH₃CN
 B) 0.1 % TFA in H₂O
 A/B = 20/80 - 20 min - 70/30, v/v
 Flow Rate : 0.4 mL/min
 Col. Temp. : 40 °C
 Detection : UV 280 nm
 Injection Vol. : 10 μL
 Data Source : LC Technical Note No. 116

Sample :
 1. Ribonuclease B
 2. Cytochrome c
 3. Lysozyme
 4. BSA
 5. Ovalbumin

Catecholamines in Urine



Conditions

Column : Inertsil ODS-4
 (5 μm, 250 × 3.0 mm I.D.)
 Eluent : A) Acetate-citrate buffer *
 B) CH₃CN
 A/B = 100/16, v/v
 Flow Rate : 0.5 mL/min
 Col. Temp. : 35 °C
 Detection : ECD 800 mV vs. Ag/AgCl
 Injection Vol. : 20 μL
 Data Source : LC Technical Note No. 93


Sample :
 1. Norepinephrine (NE)
 2. Epinephrine (E)
 3. 3,4-dihydroxybenzylamine (DHBA, I.S.)
 4. Dopamine (DA)
 (100 ng/mL in 0.1 % Acetic acid solution each)

* : Acetate-citrate buffer :
 Dissolve 0.82 g of anhydrous sodium acetate, 2.10 g of citric acid monohydrate and 0.50 g of sodium 1-octanesulfonate in 500mL of H₂O.

Applications

Foods

Organic acid




Conditions

Column : Inertsil Ph-3 (5 μ m, 250 \times 4.6 mm I.D.)
 + Inertsil CX (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : 3 mM HClO₄ in H₂O
 Reaction Reagent: 0.1 mM BTB + 30 mM Na₂HPO₄ in H₂O
 Flow Rate : 0.5 mL/min
 Col. Temp. : 35 °C
 Detection : VIS 440 nm
 Injection Vol. : 10 μ L
 Data Source : LC Technical Note No. 24

Sample :
 1. Phosphoric acid
 2. Tartaric acid
 3. Malic acid
 4. Formic acid
 5. Citric acid
 6. Lactic acid
 7. Acetic acid
 8. Succinic acid
 9. Pyroglutamic acid
 (1 mg/mL each)

Linear Fatty acids




Conditions

Column : InertSustain C18
 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) 0.1 % H₃PO₄ in CH₃CN
 B) 0.1 % H₃PO₄ in H₂O
 A/B = 10/90 -15 min
 - 90/10 -10 min - 90/10, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 210 nm
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LA901

Sample :
 1. Acetic acid (C2)
 2. Propionic acid (C3)
 3. Butyric acid (C4)
 4. Valeric acid (C5)
 5. Caproic acid (C6)
 6. Enanthic acid (C7)
 7. Caprylic acid (C8)
 8. Pelargonic acid (C9)
 9. Capric acid (C10)
 10. Undecanoic acid (C11)
 11. Lauric acid (C12)
 12. Tridecanoic acid (C13)
 13. Myristic acid (C14)
 14. Pentadecanoic acid (C15)
 (1 mg/mL each)

Preservatives and Sweetener




Conditions

Column : InertSustain Phenyl (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) 0.2 % HCOOH in H₂O
 A/B = 15/85, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 25 °C
 Detection : UV 230 nm
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LB031

Sample :
 1. Acesulfame potassium
 2. Saccharin
 3. Sorbic acid
 4. Benzoic acid
 5. Dehydroacetic acid
 (10 mg/L each)

Phenolic Antioxidants




Conditions

Column : Inertsil Ph-3 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH B) CH₃CN
 C) 5 % Acetic acid
 A/B/C = 20/20/60 - 15 min
 - 50/50/0, v/v/v
 (Mixed by a gradient mixer)
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 280 nm
 Injection : 10 μ L
 Data Source : LC Technical Note No. 65

Sample :
 1. Propyl gallate (PG)
 2. 2,4,5-Trihydroxybutyrophenone (THBP)
 3. tert-Butylhydroquinone (TBHQ)
 4. Nordihydroguaiaretic acid (NDGA)
 5. Butylated Hydroxyanisole (BHA)
 6. Octyl gallate (OG)
 7. 4-Hydroxymethyl-2,6-di-tert-butylphenol (HMBP)
 8. Dodecyl gallate (DG)
 9. Butylated hydroxytoluene (BHT)
 (10 mg/L each)

Cyanoguanidine and Melamine in Food




Conditions

Column : Inertsil HILIC (5 μ m, 150 \times 3.0 mm I.D.)
 Eluent : A) CH₃CN
 B) 10 mM Ammonium acetate
 A/B = 90/10 - 0.5 min - 90/10
 - 5.5 min - 50/50, v/v
 Flow Rate : 0.5 mL/min
 Col. Temp. : 40 °C
 Detection : LC/MS/MS (4000 QTRAP : ESI, Positive, MRM)
 Injection Vol. : 5.0 μ L
 Data Source : LC Technical Note No.132

Sample :
 1. Cyanoguanidine (20 μ g/L)
 2. Melamine (10 μ g/L)

Oxalic acid in Food




Conditions

Column : Inertsil Amide (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) 30 mM Na₂HPO₄ in H₂O (pH 6.8, H₃PO₄)
 A/B = 65/35, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 50 °C
 Detection : UV 220 nm
 Injection : 5 μ L
 Data Source : LC Technical Note No.109

Sample :
 1. Oxalic acid
 (100 mg/L)

Foods

Canthaxanthin in Tea




Conditions

Column : Inertsil ODS-3(4 μm, 150 × 4.6 mm I.D.)
 Eluent : A) CH₃OH
 B) 10 mM NaH₂PO₄ in H₂O
 A/B = 10/90 - 30 min - 50/50, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 280 nm
 Injection Vol. : 10 μL
 Data Source : LC InertSearch No. LA684

Sample :
 1. Galliccatechin (GC)
 2. Epigallocatechin (EGC)
 3. Catechin (C)
 4. Caffeine
 5. Epigallocatechin gallate (EGCg)
 6. Epicatechin (EC)
 7. Galliccatechin gallate (GCg)
 8. Epicatechin gallate (ECg)
 9. Catechin gallate (Cg)

Isoflavones




Conditions

Column : Inertsil ODS-SP (5 μm, 250 × 4.6 mm I.D.)
 Eluent : A) 0.1 % CH₃COOH in CH₃CN
 B) 0.1 % CH₃COOH in H₂O
 A/B = 15/85 - 8 min
 - 15/85 - 42 min - 35/65, v/v
 Flow Rate : 1.5 mL/min
 Col. Temp. : 35 °C
 Detection : PDA 254 nm
 Injection : 10 μL
 Data Source : LC Technical Note No. 66

Sample :
 1. Daidzin (D)
 2. Glycitin (GI)
 3. Genistin (G)
 4. 6"-O-Malonyldaidzin (MD)
 5. 6"-O- Malonylglycitin (MGI)
 6. 6"-O- Acetyldaidzin (AD)
 7. 6"-O- Acetylglycitin (AGI)
 8. 6"-O- Malonylgenistin (MG)
 9. Daizein (De)
 10. Glycitein (Gle)
 11. 6"-O- Acetylgenistin (AG)
 12. Genistein (Ge)
 (10 mg/L each)

Carotene in Food




Conditions

Column : Inertsil ODS-P
 (5 μm, 250 × 4.6 mm I.D.)
 Eluent : A) CH₃OH
 B) Ethanol
 A/B = 90/10, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C

Detection : VIS 455 nm
 Injection Vol. : 20 μL
 Data Source : LC Technical Note No.28

Sample :
 1. α-Carotene
 2. β-Carotene

Tetracycline




Conditions

Column : InertSustain C18
 (5 μm, 150 × 4.6 mm I.D.)
 Eluent : A) Imidazole buffer*
 B) CH₃OH
 A/B = 80/20, v/v (Premix)
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : FL Ex 380 nm Em 520 nm
 Injection Vol. : 20 μL

Sample :
 1. Oxytetracycline
 2. Tetracycline
 3. Chlorotetracycline
 (1 mg/L each)

* Imidazole buffer :
 Dissolve 68.08 g of imidazole, 0.37 g of disodium ethylenediaminetetraacetate and 10.72 g of magnesium acetate in 800 mL of H₂O. Adjust to pH 7.2 with acetic acid and dilute this solution to 1,000 mL with H₂O.

Putrefactive Non-Volatile Amines in Food by Pre-column HPLC



Conditions

Column : Inertsil ODS-SP (5 μm, 250 × 4.6 mm I.D.)
 Guard Column : Inertsil ODS-SP (5 μm, 10 × 4.0 mm I.D.)
 Eluent : A) CH₃CN B) H₂O A/B= 65/35, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : FL Ex 325 nm Em 525 nm
 Injection : 10 μL
 Data Source : LC Technical Note No. 48


Sample :
 1. Putrescine (5 mg/L)
 2. Cadaverine (5 mg/L)
 3. Histamine (100 mg/L)
 4. 1,8-Diaminooctane (10 mg/L)
 5. Tyramine (25 mg/L)
 6. Spermidine (5 mg/L)

Food Dyes

Conditions

Column : Inertsil ODS-3
 (5 μm, 150 × 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) 10 mM Na₂HPO₄ in H₂O
 A/B = 10/90 - 50 min - 35/65, v/v

Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 270 nm
 Injection Vol. : 10 μL
 Data Source : LC InertSearch No. LA509




Sample :

- 1. Tartrazine (Food Yellow No. 4, 7.6 mg/L)
- 2. Amaranth (Food Red No. 2, 3.8 mg/L)
- 3. Ingigocarmine (Food Blue No. 2, 7.6 mg/L)
- 4. New cocine (Food Red No. 102, 3.8 mg/L)
- 5. Sunset Yellow FCF (Food Yellow No. 5, 5.3 mg/L)
- 6. Naphthol Yellow S (7.6 mg/L)
- 7. Uranine (3.8 mg/L)
- 8. Allura red AC (5.3 mg/L)
- 9. Ponceau R (7.6 mg/L)
- 10. Ponceau SX (5.3 mg/L)
- 11. Orange I (5.3 mg/L)
- 12. Fast green FCF (Food Green No. 3, 3.0 mg/L)
- 13. Brilliant blue FCF (Food Blue No. 1, 3.0 mg/L)
- 14. Ponceau 3R (7.6 mg/L)
- 15. Erythrosine (Food Red No. 3, 5.3 mg/L)
- 16. Azure Blue VX (Sulfan blue, 3.0 mg/L)
- 17. Orange II (7.6 mg/L)
- 18. Acid red (Food Red No. 106, 3.0 mg/L)

Applications

Environmental

Non-ionic Surfactant




Conditions

Column : InertSustain C18 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH B) 10 mM Na₂B₄O₇ in H₂O A/B = 38/62, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 510 nm
 Injection Vol. : 20 μ L
 Data Source : LC InertSearch No. LA974

Sample :
 1. Heptaoxyethylene dodecyl ether
 [Deriv.](0.002 mg/L)

Anion Surfactant (Toluene addition sample)




Conditions

Column : Inertsil ODS-3 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : 0.1 M NaClO₄ in CH₃CN/H₂O = 65/35, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : FL Ex 221 nm Em 284 nm
 Injection Vol. : 10 μ L
 Data Source : LC Technical Note No. 102

Sample :
 1. Sodium Decylbenzenesulfonate(C10)
 2. Sodium Undecylbenzenesulfonate(C11)
 3. Sodium Dodecylbenzenesulfonate(C12)
 4. Toluene
 5. Sodium Tridecylbenzenesulfonate(C13)
 6. Sodium Tetradecylbenzenesulfonate(C14)
 (1 mg/L each)

Haloacetic acids




Conditions

Column : InertSustain C18 (3 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH B) 0.2 % HCOOH in H₂O A/B = 5/95 - 38 min - 100/0 -12 min - 100/0, v/v
 Flow Rate : 0.2 mL/min
 Col. Temp. : 30 °C
 Detection : LC/MS/MS (4000 QTRAP : ESI, Negative, MRM)
 Injection Vol. : 100 μ L
 Data Source : LC Technical Note No. 125

Sample :
 1. Monochloroacetic acid (MCAA)(2 μ g/L)
 2. Dichloroacetic acid (DCAA)(4 μ g/L)
 3. Trichloroacetic acid (TCAA)(20 μ g/L)

Chlorophenol




Conditions

Column : InertSustain C18 HP (3 μ m, 100 \times 2.1 mm I.D.)
 Eluent : A) CH₃OH B) CH₃OH/H₂O = 10/90, v/v A/B = 40/60 - 8 min - 90/10 - 0.5 min - 90/10 - 0.1 min - 40/60 - 5 min - 40/60, v/v
 Flow Rate : 0.3 mL/min
 Col. Temp. : 40 °C
 Detection : LC/MS (4000 QTRAP : APCI, Neg, SIM)
 Injection Vol. : 25 μ L

Sample :
 1. Phenol
 2. 2-Chlorophenol
 3. 4-Chlorophenol
 4. 2,6-Chlorophenol
 5. 2,4-Chlorophenol
 6. 2,4,6-Chlorophenol
 (0.83 μ g/L in H₂O each)

13 Kinds of Aldehydes




Conditions

Column : InertSustain C18 HP (3 μ m, 150 \times 3.0 mm I.D.)
 Eluent : A) CH₃CN B) H₂O C) THF A/B/C = 35/55/10, v/v/v
 Flow Rate : 1.5 mL/min
 Col. Temp. : 40 °C
 Detection : UV 360 nm
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LA962

Sample :
 1. DNPH-Formaldehyde
 2. DNPH-Acetaldehyde
 3. DNPH-Acetone
 4. DNPH-Acrolein
 5. DNPH-Propionaldehyde
 6. DNPH-Crotonaldehyde
 7. DNPH-Methylethylketone
 8. DNPH-Methacrolein
 9. DNPH-n-Butyraldehyde
 10. DNPH-Benzaldehyde
 11. DNPH-n-Valeraldehyde
 12. DNPH-m-Tolualdehyde
 13. DNPH-Hexanal
 (150 μ g/L each)

Aromatic Hydrocarbons



Conditions


Column : Inertsil ODS-P (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH B) CH₃OH/H₂O = 70/30, v/v A/B = 0/100 - 5 min - 0/100 - 20 min - 100/0, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 254 nm
 Data Source : InertSearch for LC No. LA336

Sample :
 1. Naphthalene
 2. Acenaphthylene
 3. Acenaphthene
 4. Fluorene
 5. Phenanthrene
 6. Anthracene
 7. Fluorene
 8. Pyrene
 9. Benzo-[a]-anthracene
 10. Chrysene
 11. Benzo-[b]-fluoranthene
 12. Benzo-[k]-fluoranthene
 13. Benzo-[a]-pyrene
 14. Dibenzo-[a]-pyrene
 15. Benzo-[ghi]-pyrene
 16. Indeno-[1,2,3-cd]-pyrene
 (10 ng/mL in CH₃OH each)

Pesticides

160 pesticides

Provided by AB SCIEX




Conditions

Column : InertSustain C18
(2 μ m, 100 \times 2.1 mm I.D.)
Eluent : A) 2 mM CH₃COONH₄
B) CH₃OH
A/B = 5/95 - 0.5 min - 30/70 - 9.5 min
- 95/5 - 5 min - 95/5, v/v
Flow Rate : 0.3 mL/min
Col.Temp. : 40 $^{\circ}$ C
Detection : LC/MS/MS
(4000 QTRAP : ESI, MRM)
Injection Vol. : 10 μ L
Data Source : LC Technical Note No. 129

118 Pesticides

Provided by Doctor KADOKAMI




Conditions


Column : Inertsil ODS-4 HP (3 μ m, 150 \times 2.1 mm I.D.)
Eluent : A) 5 mM CH₃COONH₄ in CH₃OH
B) 5 mM CH₃COONH₄ in H₂O
A/B = 5/95 - 30 min - 95/5 - 20 min - 95/5, v/v
Flow Rate : 0.3 mL/min
Col. Temp. : 40 $^{\circ}$ C
Detection : LC/MS/MS (4000 QTRAP : ESI, Positive, MRM)
Injection Vol. : 2.5 μ L
Data Source : LC InertSearch No. LA 843

Tap Water Pesticides

Positive




Negative



Conditions

Column : InertSustain C18 (3 μ m, 75 \times 2.1 mm I.D.)
Eluent : A) 5 mM CH₃COONH₄ in H₂O B) 5 mM CH₃COONH₄ in CH₃OH
A/B = 95/5 - 4 min - 60/40 - 35 min - 25/75 - 5 min - 0/100 - 6 min - 0/100, v/v
Flow Rate : 0.15 mL/min
Col.Temp. : 40 $^{\circ}$ C
Sample.Temp. : 5 $^{\circ}$ C
Detection : LC/MS/MS (4000 QTRAP : ESI, MRM)
Injection Vol. : 100 μ L
Data Source : LC Technical Note No. 135

4 Pesticides




Conditions

Column : InertSustain Phenyl (3 μ m, 150 \times 2.1 mm I.D.)
Eluent : A) 0.1 % HCOOH in CH₃CN
B) 0.1 % HCOOH in H₂O
A/B = 40/60 - 10 min - 70/30 - 0.01 min - 40/60
- 5 min - 40/60, v/v
Flow Rate : 0.3 mL/min
Col. Temp. : 40 $^{\circ}$ C
Detection : LC/MS/MS (4000 QTRAP : ESI, Positive, MRM)
Injection Vol. : 5 μ L
Data Source : LC InertSearch No. LB077

Sample :
1. Paclobutrazole
2. Diniconazole
3. Propiconazole
4. Difenconazole
(1 μ g/L each)

Carbamate Insecticides



Conditions


Column : InertSustain C18 (5 μ m, 250 \times 4.6 mm I.D.)
Eluent : A) CH₃OH B) H₂O
A/B = 35/65 - 2 min - 35/65 - 0.1 min - 53/47 - 18.4 min - 53/47 - 0.1 min
- 70/30 - 9.4 min - 70/30 - 0.1 min - 35/65 - 9.9 min - 35/65, v/v
Reaction Reagent : OPA reagent
Flow Rate : 1.0 mL/min
Col. Temp. : 40 $^{\circ}$ C
Detection : FL Ex 339 nm Em 455 nm(0 - 18.5 min), Ex 312 nm Em 382 nm(18.6 - 20.1 min),
Ex 339 nm Em 455 nm(20.2 - 30 min)
Injection Vol. : 10 μ L
Data Source : LC InertSearch No. LA916

Sample :
1. Oxamyl 5. Carbaryl 9. Fenobucarb
2. Methomyl 6. Ethiofencarb 10. Methiocarb
3. Aldicarb 7. Pirimicarb (1.0 mg/L each)
4. Bendiocarb 8. Isoprocarb

Applications

Vitamins

Vitamin A in Food




Conditions

Column : Inertsil ODS-3 (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH
 B) H₂O
 A/B = 95/5, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : UV 325 nm
 Injection Vol. : 20 μ L
 Data Source : LC Technical Note No. 32

Sample :
 1. Retinol (50 μ g/L)

Vitamin B1




Conditions

Column : Inertsil ODS-3 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH
 B) 0.01 M NaH₂PO₄,
 0.15 M NaClO₄ in H₂O (pH2.2)
 A/B = 1/9, v/v
 Reaction Reagent : 0.05 w/v % K₃[Fe(CN)₆]
 +15 w/v % NaOH, 0.4 mL/min
 Flow Rate : 0.8 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : FL Ex 375 nm Em 440 nm
 Injection Vol. : 20 μ L
 Data Source : LC Technical Note No. 10

Sample :
 1. Thiamine
 2. Hydroxyethyl thiamine (HET)
 (10 μ g/L each)

Vitamin C in Food




Conditions

Column : Inertsil SIL-100A (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) CH₃COOC₂H₅
 B) n-Hexane
 C) CH₃COOH
 A/B/C = 50/40/10, v/v/v
 Flow Rate : 1.5 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : VIS 495 nm
 Injection Vol. : 20 μ L
 Data Source : LC Technical Note No. 9

Sample :
 1. Ascorbic acid

Vitamin D2, D3




Conditions

Column : Inertsil ODS-P (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : CH₃CN
 Flow Rate : 1.5 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : UV 265 nm
 Injection Vol. : 200 μ L
 Data Source : LC Technical Note No. 33

Sample :
 1. Vitamin D2 (Ergocalciferol)
 2. Vitamin D3 (Cholecalciferol)
 (0.1 mg/L each)

Vitamin E




Conditions

Column : Inertsil NH2 (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) 2-Propanol
 B) n-Hexane
 A/B = 2/98, v/v
 Flow Rate : 1.2 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : FL Ex 298 nm Em 330 nm
 Injection Vol. : 20 μ L
 Data Source : LC InertSearch No. LB199

Sample :
 1. α -Tocopherol
 2. β -Tocopherol
 3. γ -Tocopherol
 4. δ -Tocopherol
 (1.0 mg/L each)

Vitamin K




Conditions

Column : InertSustain C8 (5 μ m, 150 \times 3.0 mm I.D.)
 Eluent : CH₃CN
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : UV 270 nm
 Injection Vol. : 5 μ L
 Data Source : LC InertSearch No. LB030

Sample :
 1. Vitamin K2 (MK-4)
 2. Vitamin K1
 3. Vitamin K2 (MK-7)
 (5 mg/L each)

Others

Sugar Analysis




Conditions

Column : InertSustain NH2 (5 μ m, 250 \times 4.6 mm I.D.)
 Eluent : A) CH₃CN
 B) H₂O
 A/B = 85/15, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : RI
 Injection Vol. : 10 μ L
 Data Source : LC InertSearch No. LB180

Sample :
 1. Rhamnose
 2. Fucose
 3. Fructose
 4. Glucose
 5. Sucrose
 6. Maltose
 7. Lactose
 (10 mg/mL each)

Sugar Analysis Using HPLC-ECD




Conditions

Column : InertSphere Sugar-1 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : 100 mM NaOH
 Flow Rate : 0.5 mL/min
 Col. Temp. : 25 °C
 Detection : ECD Pulse Mode
 Injection Vol. : 10 μ L
 Data Source : LC Technical Note No. 101

Sample :
 1. Fucose
 2. Glucose
 3. Fructose
 4. Lactose
 5. Sucrose
 (10 mg/L each)

Analysis of Pre-column Derivatized Amino Acids




Conditions

Column : Inertsil ODS-4(3 μ m, 150 \times 3.0 mm I.D.)
 Eluent : A) CH₃CN/CH₃OH/H₂O = 45/40/15, v/v/v
 B) 20 mM KH₂PO₄ (pH 6.9, H₃PO₄)
 A/B = 11/89 - 3 min - 11/89 - 9 min
 - 22/78 - 2 min - 28/72 - 9 min - 30/70
 - 4 min - 65/35 - 7 min - 75/25 - 1 min
 - 100/0, v/v
 Flow Rate : 0.7 mL/min
 Col. Temp. : 35 °C
 Detection : FL Ex 350 nm Em 450 nm (0-29 min)
 Ex 266 nm Em 305 nm (29-35 min)
 Injection Vol. : 1 μ L
 Sample : Derivatized Amino Acids
 Data Source : LC InertSearch No. LB088

Sample :
 1. OPA-Aspartic Acid
 2. OPA-Glutamic Acid
 3. OPA-Asparagine
 4. OPA-Serine
 5. OPA-Glutamine
 6. OPA-Histidine
 7. OPA-Glycine
 8. OPA-Threonine
 9. OPA-Citru line
 10. OPA-Arginine
 11. OPA-Alanine
 12. OPA-GABA
 (4-aminobutanoic acid)
 13. OPA-Tyrosine
 14. OPA-Cys-Cys
 15. OPA-Valine
 16. OPA-Methionine
 17. OPA-Tryptophan
 18. OPA-Phenylalanine
 19. OPA-Isoleucine
 20. OPA-Leucine
 21. OPA-Lysine
 22. Fmoc-Proline
 (10 μ g/mL each)

PFOS and PFOA




Conditions

Column : Inertsil ODS-4
 (3 μ m, 100 \times 2.1 mm I.D.)
 Eluent : A) 5 mM CH₃COONH₄ in H₂O
 B) 5 mM CH₃COONH₄ in CH₃OH
 A/B = 40/60 - 8 min - 25/75 - 0.1 min
 - 10/90 - 1.9 min - 10/90 - 0.1 min
 - 40/60 - 4.9 min - 10/90, v/v
 Flow Rate : 0.6 mL/min
 Col. Temp. : 40 °C
 Detection : LC/MS/MS (4000 QTRAP : ESI, Negative, MRM)
 Injection Vol. : 2 μ L
 Data Source : LC InertSearch No. LA864

Sample :
 1. PFOA (Perfluorooctanoic acid)
 2. PFOS (Perfluorooctanesulfonic acid)
 (1 mg/L each)


Analysis of Pre-column Derivatized Aflatoxin



Conditions

Column : InertSustain C18 (5 μ m, 150 \times 4.6 mm I.D.)
 Eluent : A) CH₃OH B) CH₃CN C) H₂O
 A/B/C = 30/10/60, v/v/v (Premix)
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : FL Ex 365 nm Em 450 nm
 Injection Vol. : 20 μ L
 Data Source : LC InertSearch No. LB107

Sample :
 1. Aflatoxin G1
 2. Aflatoxin B1
 3. Aflatoxin G2
 4. Aflatoxin B2
 (5 ng/mL each)

Hydrogen Peroxide (H₂O₂)

Conditions

Column : Inertsil CX (5 μ m, 250 \times 4.6 mm I.D.)
 Flow Rate : 0.8 mL/min
 Detection : ECD
 Injection Vol. : 100 μ L
 Data Source : LC Technical Note No. 49

Sample :
 1. Hydrogen peroxide(H₂O₂)
 (10 μ g/L)