

InertSustain® C18

Physical Properties

- Silica : Newly Developed ES Silica Gel
- Particle Size : 2 µm, 3 µm, 5 µm
- Surface Area : 350 m²/g
- Pore Size : 100 Å (10 nm)
- Pore Volume : 0.85 mL/g
- Bonded Phase : Octadecyl Groups
- End-capping : Complete
- Carbon Loading : 14 %
- USP Code : L1
- pH Range : 1~10



Generally, silica based columns are mechanically stable and provide high efficiencies, however, they cannot be used under alkaline conditions and their residual silanol groups tend to adsorb organic bases. InertSustain C18 employs a radically new type of silica, in which the surface of the silica is uniquely modified, enabling precise control of the silica properties.

InertSustain C18 inherits the advantages of all the current Inertsil HPLC columns (e.g., extremely low operating back pressure, superior inertness to typically any analytes, high efficiency and compatibility with a wide range of solvents), but now can be used for wide pH analysis with consistent performance from column to column and lot to lot.

Figure 1 : Basic Compounds




Figure 2 : Acidic Compounds






Figure 3 : Chelating Compound




Comparison of Performance to Core-Shell Columns

As shown below, core-shell columns show peak tailing due to the presence of trace metals or silanol groups in their silica gel. Quantitative and qualitative analysis will be a source of concern since the adsorption of compounds can be expected.

InertSustain C18 3 μm



Kinetex C18 (1.7 μm)




Conditions

Column : ODS Column(100 × 2.1 mm I.D.)
 Eluent : A) 2 mM CH₃COONH₄ in 95 % CH₃CN
 B) 2 mM CH₃COONH₄ in 10 % CH₃CN
 A / B = 20 / 80 - 2 min - 100 / 0 - 2.5 min
 - 100 / 0 - 0.01 min - 20 / 80, v / v
 (Mixed by a gradient mixer)
 Flow Rate : 0.3 mL / min
 Col. Temp. : 40 °C
 Detection : LC / MS / MS
 (4000 QTRAP® : ESI, Positive, MRM)
 Injection Vol. : 10 μL
 Sample : 1. Nitrofurazone (100 μg / L)
 2. Lasalocid A (100 μg / L)

Wide pH compatibility with Long Column Lifetime

As shown in the experiment below, due to the introduction of Evolved Surface Silica, InertSustain C18 maintained high efficiency and peak shape for 300 hours while other "wide pH" column brands failed.



Purging Conditions

Column Size : 5 μm, 150 × 4.6 mm I.D.
 Eluent : A) CH₃OH
 B) 50 mM Triethylamine (pH 10.0)
 A/B = 30/70, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 50 °C

Analytical Conditions


Eluent : A) CH₃CN
 B) H₂O
 A/B = 65/35, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 40 °C
 Detection : UV 254 nm
 Sample : Naphthalene

Experience the InertSustain! (Inertness and Sustainability)

Highly end-capped ODS column such as InertSustain C18 offers an opportunity to flush out contaminants from the column surface easily using an organic solvent. Coffee melanoidins are brown heterogeneous polymers contained in coffee. Its components are not clarified yet, but it is considered to contain several ionic compounds, which a poorly end-capped column will adsorb those ionic compounds leading to short column lifetime.

As for ODS column, which is commonly used for HPLC and LC/MS/MS, its inertness has an influence not only on peak shape but also detection sensitivity and durability. It is highly recommended to use highly end-capped column which provides good peak shape for both basic and acidic compounds such as InertSustain C18.

The packing material was visually confirmed by removing the column



	Brand A	InertSustain C18
Before Experiment		
Injection of Coffee		
After loading Coffee		
After washing the column with CH ₃ CN 100 % 10 minute washing		


Ionic contaminants were hard to be washed out from the Column

InertSustain® C18

Reliable Reproducibility, Performance and Quality


Rigorous quality control of physical properties and strict chromatographic tests for inertness and selectivity, contribute to the production of InertSustain C18 with an outstanding reproducibility and long column lifetime.

Strong Basic Compound Test




Conditions
 Column Size : 5 µm, 250 × 4.6 mm I.D.
 Eluent : A) CH₃CN B) 25 mM phosphate buffer (pH 7.0)
 A / B = 40 / 60, v / v
 Flow Rate : 1.0 mL / min
 Col. Temp. : 40 °C
 Detection : UV 220 nm
 Sample : 1. Dextromethorphan

Same Retention Behavior between Various Particle Sizes




Conditions
 Column Size : 50 × 2.1 mm I.D.
 Eluent : A) CH₃OH B) 25 mM phosphate buffer (pH 7.0)
 A / B = 30 / 70, v / v
 Flow Rate : 0.2 mL / min
 Col. Temp. : 40 °C
 Detection : UV 230 nm
 Sample : 1. Uracil
 2. Pyridine
 3. Phenol
 4. Berberine chloride
 5. Dextromethorphan

Strong Acidic Compound Test



Conditions
 Column Size : 5 µm, 150 × 4.6 mm I.D.
 Eluent : A) CH₃CN B) 0.1 % H₃PO₄
 A / B = 25 / 75, v / v
 Flow Rate : 1.0 mL / min
 Col. Temp. : 40 °C
 Detection : UV 254 nm
 Sample : 1. Brilliant Blue FCF
 2. Phenol
 3. Salicylic acid

Strong Chelating Compound Test



Conditions
 Column Size : 5 µm, 150 × 4.6 mm I.D.
 Eluent : A) CH₃CN B) 0.1 % H₃PO₄
 A / B = 40 / 60, v / v
 Flow Rate : 1.0 mL / min
 Col. Temp. : 40 °C
 Detection : UV 254 nm
 Sample : 1. Hinokitiol

Analytical Columns

Particle Size: 2 µm	Length \ I.D. (mm)	2.1	3.0
	30	5020-14351	5020-14361
	50	5020-14352	5020-14362
	75	5020-14353	5020-14363
	100	5020-14354	5020-14364
	150	5020-14355	5020-14365

HP Series Particle Size: 3 µm 50 MPa (500 Bar)	Length \ I.D. (mm)	2.1	3.0	4.6
	30	5020-14411	5020-14421	5020-14441
	50	5020-14412	5020-14422	5020-14442
	75	5020-14413	5020-14423	5020-14443
	100	5020-14414	5020-14424	5020-14444
	150	5020-14415	5020-14425	5020-14445
250	5020-14416	5020-14426	5020-14446	

* End-fittings are 1/16" Waters-compatible.
 * UHPLC compatible end-fittings are also available upon request for UHPLC systems (Ex: UPLC) to avoid dead volume.
 * Indicate "UP Type end-fittings" when ordering. (Please note that UP type is not available for a 4.6 mm I.D. column)
 * For maximum operating pressure information, please refer to page 46.

Particle Size: 3 µm	Length \ I.D. (mm)	1.0	1.5		
	30	5020-14301	5020-14311		
	50	5020-14302	5020-14312		
	75	5020-14303	5020-14313		
	100	5020-14304	5020-14314		
	150	5020-14305	5020-14315		
Particle Size: 3 µm	Length \ I.D. (mm)	2.1	3.0	4.0	4.6
	30	5020-07411	5020-07421	5020-07431	5020-07441
	50	5020-07412	5020-07422	5020-07432	5020-07442
	75	5020-07413	5020-07423	5020-07433	5020-07443
	100	5020-07414	5020-07424	5020-07434	5020-07444
	125	5020-07417	5020-07427	5020-07437	5020-07447
	150	5020-07415	5020-07425	5020-07435	5020-07445
	250	5020-07416	5020-07426	5020-07436	5020-07446

Particle Size: 5 µm	Length \ I.D. (mm)	1.0	1.5		
	30	5020-14201	5020-14211		
	50	5020-14202	5020-14212		
	75	5020-14203	5020-14213		
	100	5020-14204	5020-14214		
	150	5020-14205	5020-14215		
Particle Size: 5 µm	Length \ I.D. (mm)	2.1	3.0	4.0	4.6
	30	5020-07311	5020-07321	5020-07331	5020-07341
	50	5020-07312	5020-07322	5020-07332	5020-07342
	75	5020-07313	5020-07323	5020-07333	5020-07343
	100	5020-07314	5020-07324	5020-07334	5020-07344
	125	5020-07317	5020-07327	5020-07337	5020-07348
	150	5020-07315	5020-07325	5020-07335	5020-07345
	250	5020-07316	5020-07326	5020-07336	5020-07346

* End-fittings are 1/16" Waters-compatible.
 * For maximum operating pressure information, please refer to page 46.

Cartridge Guard Column E

I.D. of the Analytical Column Applicable (mm)	Length (mm)	I.D. (mm)	Replacement Cartridge E Guard Column (2 EA.)		Cartridge E Holder / Cartridge Set (2 Cartridge E Guard Columns & 1 Holder)	
			Particle Size		Particle Size	
			3 µm	5 µm	3 µm	5 µm
1.0	10	1.0	5020-19250	5020-19249	5020-19300	5020-19299
1.5, 2.1		1.5	5020-19350	5020-19349	5020-19400	5020-19399
2.1, 3.0		3.0	5020-19150	5020-19149	5020-19200	5020-19199
4.0, 4.6		4.0	5020-19050	5020-19049	5020-19100	5020-19099
2.1, 3.0	20	3.0	5020-19550	5020-19549	5020-19600	5020-19599
4.0, 4.6		4.0	5020-19450	5020-19449	5020-19500	5020-19499
Holder for Cartridge Guard Column E				For 10 mm Length		5020-08500
				For 20 mm Length		5020-08550

* End-fittings are 1/16" Waters-compatible.
 * For maximum operating pressure information, please refer to page 46.

Reversed Phase Columns

HILIC Columns

Normal Phase Columns

SEC Columns

Ion Exchange Columns

Application Specific Columns

Guard Columns

Preparative Columns

Capillary Columns

Applications

Cat. No. Index

InertSustainSwift™ C18



Physical Properties

- Silica : Newly Developed ES Silica Gel
- Particle Size : 1.9 µm, 3 µm, 5 µm
- Surface Area : 200 m²/g
- Pore Size : 200 Å (20 nm)
- Pore Volume : 1.00 mL/g
- Bonded Phase : Octadecyl Groups
- End-capping : Complete
- Carbon Loading : 9.0 %
- USP Code : L1
- pH Range : 1.0~10.0

As shown in figure 1, InertSustainSwift C18 maintains the same extreme inertness, wide pH range and provide rapid separations with symmetric peaks. The optimization of surface area, pore size and chemical bonding delivers superior peak shapes (Figure 2). Figure 3 proves InertSustainSwift C18 is also ideal for LC/MS/MS applications which offer highly sensitive results and enables MS compatible buffers to be used due to the extremely inert silica gel.

Figure 1 : Comparison of Retentivity




Figure 2 : Comparison of Efficiency





Figure 3 : Rapid LC/MS/MS Analysis of Basic Drugs



Analytical Columns

Particle Size: 1.9 µm	Length / I.D. (mm)	2.1	3.0	
	50	5020-88228	5020-88233	
	100	5020-88230	5020-88235	
	150	5020-88231	5020-88236	
HP Series Particle Size: 3 µm Max. Operating Pressure : 50 MPa (500 Bar)	Length / I.D. (mm)	2.1	3.0	4.6
	50	5020-88210	5020-88216	5020-88222
	100	5020-88212	5020-88218	5020-88224
	150	5020-88213	5020-88219	5020-88225
	250	5020-88214	5020-88220	5020-88226

* End-fittings are 1/16" Parker style.

* For maximum operating pressure information, please refer to page 46.

Particle Size: 3 µm	Length / I.D. (mm)	1.0	1.5			
	30	5020-88160	5020-88166			
	50	5020-88161	5020-88167			
	75	5020-88162	5020-88168			
	100	5020-88163	5020-88169			
	150	5020-88164	5020-88170			
	250	5020-88165	5020-88171			
Particle Size: 3 µm	Length / I.D. (mm)	2.1	3.0	4.0	4.6	
	30	5020-88124	5020-88131	5020-88138	5020-88145	
	50	5020-88125	5020-88132	5020-88139	5020-88146	
	75	5020-88126	5020-88133	5020-88140	5020-88147	
	100	5020-88127	5020-88134	5020-88141	5020-88148	
	125	5020-88253	5020-88254	5020-88255	5020-88256	
	150	5020-88128	5020-88135	5020-88142	5020-88149	
	250	5020-88129	5020-88136	5020-88143	5020-88150	
Particle Size: 5 µm	Length / I.D. (mm)	1.0	1.5			
	30	5020-88038	5020-88044			
	50	5020-88039	5020-88045			
	75	5020-88040	5020-88046			
	100	5020-88041	5020-88047			
	150	5020-88042	5020-88048			
	250	5020-88043	5020-88049			
	Particle Size: 5 µm	Length / I.D. (mm)	2.1	3.0	4.0	4.6
		30	5020-88001	5020-88008	5020-88015	5020-88022
		50	5020-88002	5020-88009	5020-88016	5020-88023
		75	5020-88003	5020-88010	5020-88017	5020-88024
		100	5020-88004	5020-88011	5020-88018	5020-88025
		125	5020-88249	5020-88250	5020-88251	5020-88252
		150	5020-88005	5020-88012	5020-88019	5020-88026
		250	5020-88006	5020-88013	5020-88020	5020-88027

* End-fittings are 1/16" Parker style.

* For maximum operating pressure information, please refer to page 46.

Cartridge Guard Column E

I.D. of the Analytical Column Applicable (mm)	Length (mm)	I.D. (mm)	Replacement Cartridge E Guard Column (2 EA.)		Cartridge E Holder / Cartridge Set (2 Cartridge E Guard Columns & 1 Holder)	
			Particle Size		Particle Size	
			3 µm	5 µm	3 µm	5 µm
1.0	10	1.0	5020-88199	5020-88105	5020-88200	5020-88106
1.5, 2.1		1.5	5020-88201	5020-88107	5020-88202	5020-88108
2.1, 3.0		3.0	5020-88197	5020-88103	5020-88198	5020-88104
4.0, 4.6	20	4.0	5020-88195	5020-88101	5020-88196	5020-88102
2.1, 3.0		3.0	5020-88205	5020-88111	5020-88206	5020-88112
4.0, 4.6		4.0	5020-88203	5020-88109	5020-88204	5020-88110
Holder for Cartridge Guard Column E				For 10 mm Length		5020-08500
				For 20 mm Length		5020-08550

* End-fittings are 1/16" Waters-compatible.

* For maximum operating pressure information, please refer to page 46.