

Alltech® Custom Packed Metal and PTFE Columns

Before ordering a custom-packed GC column, consult the list of popular packed columns on pages 244–245 to see if the column you need is there. Each column comes complete with brass nuts and brass ferrules. For assistance or a quote, please contact Technical Support.

Coiling Schematics




Figure 1




Figure 2




Figure 3




Figure 4




Figure 5




Figure 6




Figure 7




Figure 8




Figure 9

Custom Metal and PTFE Columns — Ordering Instructions

Specify Coiling

- 1) Instrument Make and Model: _____
- 2) Configuration Number: (see coiling schematics) _____
- 3) Dimensions:
 - A= _____ in/mm (circle one)
 - B= _____ in/mm (circle one)
 - C= _____ in/mm (circle one)
 - D= _____ in/mm (circle one)
 - E= _____ in/mm (circle one)
 - F= _____ in/mm (circle one)
 - X= _____ in/mm (circle one)
- 4) Column Length: _____ ft/m (circle one)
- 5) Column i.d.: _____ in/mm (circle one)
- 6) Column o.d.: _____ in/mm (circle one)
- 7) Tubing Material: (circle one)
 - Stainless Steel AT-Steel Nickel PTFE
 - TFE-Coated SS Aluminum Copper
- 8) Solid Support: _____ Mesh Size: _____

Specify Packing

Choose either a Pretested Packing (p368–369) or create your own custom (refer to p370–371) by specifying information below:

- 9) Stationary Phase A: _____ % Loading: _____
- Stationary Phase B*: _____ % Loading: _____
- Stationary Phase C*: _____ % Loading: _____
- 10) Mirror Image: _____ Please check if a mirror image is required.
- 11) Injection Port Design: (circle one)
 - on-column not on-column
- 12) Stainless Steel Fittings: _____ Please check if stainless steel fittings are required (PN: C-14010N)
- 13) Preconditioning: _____ Please check if column preconditioning is required (PN: C-14005)
- 14) Other Special Instructions: (please be specific)

*Additional charge applies to dual- and triple-coated packings (Part No. C-14020).

Custom Packed Columns

Description	Part No.
Custom Packed Columns (Specify all information from Table 1)	C-5000
Empty Columns	
Empty Column with Fittings and Ferrules	C-5100
Additional Charges (per Column)	
Preconditioned Column	C-14005
Stainless Steel Fittings	C-14010N
Dual- and Triple-Coated Packings	C-14020

more info

Certain expensive stationary phases and packings may require a surcharge. Be sure to state coiling instructions.

Alltech® GC Packed Columns

Custom Deactiglas® Glass

- Silane treated for inertness
- Precision bore glass for reproducibility
- Brass fittings and vespel/graphite ferrules included
- Available packed or empty




Alltech® Deactiglas® GC Columns are fabricated from precision-bore (2mm ±0.05mm or 4mm ±0.1mm i.d.) glass. All columns are 1/4" o.d. unless otherwise specified. Columns longer than 2 meters may have joints. Several additional services such as stainless steel fittings and preconditioning of the packing are available. Certain expensive stationary phases and packings may require a surcharge.

Due to space limitations we are not able to list all of the different manufacturers columns in our catalog. However, they are still available and can be specified at time of order placement. If you do not see your instrument on the following pages, please call Customer Service.

Special Services for Glass Columns

Description	Part No.
Stainless Steel Fittings (Per Column)	C-14010N
Preconditioned Packing (Per Column)	C-14005

Key to Configurations for Glass Columns







Deactiglas® Column Ordering Instructions:

1. Choose column configuration part no. below.
2. Choose either a pretested packing (p247–248) or create a custom packing (refer to p249–250) by specifying information below:
 Stationary Phase A _____ % Loading _____
 Stationary Phase B _____ % Loading _____
 Stationary Phase C _____ % Loading _____

Alltech® GC Packed Columns

Custom Deactiglas® Glass

Deactiglas® Glass GC Columns (continued)

Chromatograph	Configuration	Length	i.d.	Packed Part No.	Empty Part No.
Agilent 5880, 5890, 5987, 6890* 	Configuration A On-Column Injection A1 = 11.02" (280mm) A2 = 9.05" (230mm) S = 9" (229mm) FID	3ft	2mm	130962	130961
		4ft	2mm	C6082	6082
		4ft	4mm	C6085	6085
		6ft	2mm	C6087	6087
		6ft	4mm	C6089	6089
		2m	2mm	C6336	6336
		2m	4mm	C6343	—
		8ft	2mm	C6347	6347
		1.8m	3mm	C6346	6346
		10ft	2mm	C6091	6091
		10ft	4mm	131552	131551
			Configuration B On-Column Injection A1 = 11.02" (280mm) A2 = 7.09" (180mm) S = 9" (229mm) TCD	6ft	2mm
6ft	4mm			C6391	6391
1.8m	3mm			C6346	6346
	Configuration C Not On-Column Injection A1 and A2 = 9.05" (230mm) S = 9" (229mm) All Detectors Use with Liner or Disposable Insert	3ft	4mm	C6433	6433
		4ft	2mm	C6436	6436
		4ft	4mm	C6440	6440
		6ft	2mm	C6442	6442
		2m	2mm	C6447	6447
		2m	4mm	C6449	6449
PerkinElmer® 115, 300, 900, 910, 990, F30, 2000, 2100, 3920, Sigma Series 	On-Column Injection A1 = 12.62" (321mm) A2 = 7.09" (180mm) S = 8.75" (222mm)	6ft	4mm	C6151	6151

gc columns | packed




more info

For explanation of A1, A2, and S, refer to key on page 240.

Alltech® GC Packed Columns

Custom Deactiglas® Glass

Deactiglas® Glass GC Columns (continued)

Chromatograph	Configuration	Length	i.d.	Packed Part No.	Empty Part No.
<i>PerkinElmer® 8300, 8400, 8700 Series, 9000 Series, Autosystem</i>					
 <p>2534A</p>	Not On-Column Injection A1 = 6.75" (171mm) A2 = 6.75" (171mm) S = 6.5" (165mm)	6ft	2mm	C6479	6479
		6ft	4mm	C6482	6482
		2m	2mm	C6468	6468
<i>Shimadzu® 8A, RIA, GC-8A8IF (No Fittings Included)</i>					
 <p>2539</p>	A1 = 9" (229mm) A2 = 9" (229mm) S = 6" (152mm) 5mm o.d. Tubing				
		2m	3mm	C6493	6493
<i>Shimadzu® 7A, 9A, 12A, 14A, 14B, 15A, 16A (No Fittings Included)</i>					
 <p>2538</p>	A1 = 10.69" (272mm) A2 = 12.81" (327mm) S = 1.57" (40mm) 5mm o.d. Tubing	1m	3mm	C6497	6497
		1.7m	3mm	C6495	6495

gc columns | packed



more info

For explanation of A1, A2, and S, refer to page 240.

Alltech® GC Packed Columns

Custom Deactiglas® Glass

Deactiglas® Glass GC Columns (continued)

Chromatograph	Configuration	Length	i.d.	Packed Part No.	Empty Part No.	
<i>Varian® 3300, 3400, 3600, 3700, 4400, 4600, 6000 Vista Series</i>						
 <p>Injector A to Detector A A1 = 9.5" (241mm) A2 = 8.1" (206mm) S = 5.5" (140mm) FID</p>						
			6ft	4mm	C6192	6192
<i>Varian® Universal Columns 3300, 3400, 3600, 3700, 4400, 4600, 6000 Vista Series</i>						
 <p>A1 = 9.5" (241mm) A2 = 2.25" (57mm) S = N/A FID, TCD, ECD</p>						
			6ft	2mm	C6206	6206

more info

For explanation of A1, A2, and S, refer to page 240.

more info

Universal Glass Column Adapters are available for the Varian® 3000 and Varian® 3700/Vista Series. Please contact Customer Service for part numbers and prices.

gc columns | packed

Hi-EFF™ Pre-Columns


- Extend column life
- Reduce sample clean-up
- Perform reaction chromatography
- Perform subtraction chromatography

Hi-EFF™ Pre-Columns

Description	Length	Qty.	Empty Part No.
1/8" SS	6"	6	5910
1/8" Bottom-Drilled SS Union	—	ea	5920
1/4" Bottom-Drilled SS Union	—	ea	5921

Deactiglas® Column Inlet Liners

- Prolong the life of your glass columns
- Silane-treated for inertness
- Available for 2mm or 4mm i.d. columns
- Inserts directly into the end of your column



Deactiglas® Column Inlet Liners

Description	Qty.	Part No.
For 2mm i.d. Columns	10	6301
For 4mm i.d. Columns	10	6303
Glass Wool Inserter/Remover	2	6292

Alltech® Popular GC Packed Columns

These columns are 6' x 1/8" o.d. 304 premium grade, precleaned stainless steel (unless otherwise indicated). Each column is made with preconditioned, pretested packing and comes complete with fittings, packing quality control chromatogram, column tags, and conditioning instructions.

Popular Packed Column Ordering Instructions:

1. Choose the part no. for your column from pages 244–245.
2. Specify make and model of GC.*
3. If your column is not listed, see page 239 for custom columns.

*If no coiling instructions are provided, we will supply in a nominal 6" coil.

Popular GC Packed Columns

Packing Description	Dimensions*	Part No.
Alumina F-1, 60/80	5' x 1/8" SS Tubing	5664PC
10% Carbowax® 20M on Chromosorb® W-AW, 80/100		12005PC
10% Carbowax® 20M on Chromosorb® W-HP, 80/100		12106PC
0.20% Carbowax® 1500 on Graphpac™-GC, 80/100		8546PC
Chemipack C18, 80/100		2768PC
Chromosorb® 101, 80/100		2406PC
Chromosorb® 107, 80/100	6' x 1/8" TFE Tubing	9783
Chromosorb® 750, 80/100	4' x 1/8" SS Tubing	14480
Gas Chrom™ 220, 80/100		2484PC
Gas Chrom™ 254, 80/100		2486PC
HayeSep® D, 80/100		14487
HayeSep® D, 100/120	20' x 1/8" Ni Tubing	27082PC
HayeSep® D, 100/120	30' x 1/8" SS Tubing	27083PC
HayeSep® D, 100/120	10' x 1/8" SS Tubing	28301PC
HayeSep® DB, 100/120	30' x 1/8" SS Tubing	2836PC
HayeSep® Q, 80/100	8' x 1/8" SS Tubing	28010PC
HayeSep® Q, 80/100		2801PC
HayeSep® Q, 100/120		14489
HayeSep® P, 60/80	8' x 1/8" SS Tubing	2803PC
HayeSep® P, 80/100		2804PC
HayeSep® T, 60/80	3' x 1/8" SS Tubing	14491
HayeSep® T, 80/100		2813PC

*6' x 1/8" stainless steel unless otherwise specified.

Alltech® Popular GC Packed Columns (continued)

Popular GC Packed Columns** (continued)

Packing Description	Dimensions*	Part No.
Molecular Sieve 5A (Washed), 60/80		14494
Molecular Sieve 5A (Washed), 80/100		5605PC
Molecular Sieve 5A (Washed), 80/100	10' x 1/8" SS Tubing	14495
Molecular Sieve 13X (Washed), 80/100	3' x 1/8" SS Tubing	14492
Molecular Sieve 13X (Washed), 80/100		5773PC
Molecular Sieve 13X (Washed), 80/100	10' x 1/8" SS Tubing	14493
3% OV™-101 on Chromosorb® W-HP, 80/100		12019PC
3% OV™-101 on Chromosorb® W-HP, 100/120		12700PC
10% OV™-101 on Chromosorb® W-HP, 80/100		12703PC
3% OV™-17 on Chromosorb® W-HP, 80/100		12719PC
Porapak® N, 80/100		2716PC
Porapak® P, 80/100		14499
Porapak® Q, 50/80		2700PC
Porapak® Q, 80/100	12' x 1/8" SS Tubing	27012PC
Porapak® Q, 80/100		2701PC
Porapak® Q, 100/120		2702PC
Porapak® QS, 80/100		2719PC
Porapak® T, 80/100		2713PC
10% SE-30 on Chromosorb® W-HP, 80/100		12423PC
Silica Gel Grade 12, 60/80	18' x 1/8" SS Tubing	5651PC
Silica Gel Grade 12, 80/100	8' x 1/8" SS Tubing	14530
10% Silar™ 10C on Chromosorb® W-HP, 100/120		12430PC
Tenax® TA, 60/80	4' x 1/8" SS Tubing	14538
Tenax® TA, 60/80		4900PC
Tenax® TA, 80/100		4901PC
Unibeads™ 2S, 60/80		2760PC
Unibeads™ 3S, 60/80		14541
VZN-1, 60/80	23' x 1/8" SS Tubing	84623PC
10% Silicone UCW-98 on Chromosorb® W-HP, 80/100		8493PC


*6' x 1/8" stainless steel unless otherwise specified. **To order, use the part number shown and specify the instrument for which it should be coiled.

Alltech® CTR Concentric Packed Columns

For Gas Analysis

- CTR I is ideal for separating O₂, N₂, CH₄, CO, CO₂ at ambient conditions
- CTR III is ideal for argon analysis in the presence of O₂
- Useful where two separate chromatographic runs are required to analyze one sample
- Individually tested

A CTR column is essentially a column within a column. This permits you to use two different packings for the analysis of your sample. The diagram shows the construction of the CTR column.



Cut-away view of CTR column.


CTR I—Oxygen, Nitrogen, Methane, Carbon Monoxide, and Carbon Dioxide

Until the CTR, no single column operating at room temperature has been designed to permit, in one analysis, all these components in air. The analysis is usually accomplished by separate runs or with expensive valves. The CTR I accomplishes this separation in less than nine minutes.

Calibration Gas for CTR I Column—This special mixture is useful for calibrating Part No. 8700 CTR I columns. It contains approximately 15% CO₂ + 7% CO + 7% O₂ + 4.5% CH₄ with the balance N₂. It is supplied in a 4L push button can with a needle applicator.

CTR I Column

Description	Temp. Limit max.	Part No.
CTR I Column	275°C	8700
Calibration Gas for CTR I	—	9799



CTR III Argon Analysis Column


The CTR III is special designed for the separation of argon in the presence of oxygen. The outer column contains a special adsorbent that removes oxygen allowing the elution of argon and nitrogen. The inner column allows the combined elution of a peak containing oxygen and argon and later nitrogen. This separation is accomplished at ambient temperatures in ten minutes.

more info

For GC Applications, see the GC Applications Section, pages 462–493.

CTR III Column

Description	Temp. Limit max.	Part No.
CTR III Argon Analysis Column	275°C	8725
<i>Additional Fittings</i>		
Tube End Connector, 1/4" to 1/8" Tube	—	14169





Pretested GC Packings

Pretested packings are guaranteed and supplied with a test chromatogram. Preconditioning service is available for 20g or 25g of packing. The packing is heated at 25°C below the upper temperature limit for 3 hours with a flow of oxygen-free nitrogen.

Pretested GC Packings

Packing Description	Temp. Limits	Qty.	Part No.
Alumina F-1, 60/80	300°C	100g	5664
1% AT™-1000 on Carbowax™ 1, 60/80	250°C	15g	1733
5.0% Carbowax® 20M on Chromosorb® W-HP, 80/100	60–225°C	20g	16147
10.0% Carbowax® 20M on Chromosorb® W-AW, 80/100	60–225°C	20g	12005
10.0% Carbowax® 20M on Chromosorb® W-HP, 80/100	60–225°C	20g	12106
Carbosphere®, 60/80	400°C	5g	5680
Carbosphere®, 80/100	400°C	5g	5682
Chemipack C18, 80/100	320°C	15mL	2768
Chromosorb® 101, 80/100	275/325°C	25g	16149
30% DC-200 (500 cstk) on Chromosorb® P-AW, 60/80	20–200°C	25g	12713
Gas Chrom™ 254, 80/100	275°C	25g	2486

Pretested GC Packings (continued)

Pretested GC Packings (continued)

Packing Description	Temp. Limits	Qty.	Part No.
Molecular Sieve 5A (washed), 60/80	350°C	50g	5604
Molecular Sieve 5A (washed), 80/100	350°C	50g	5605
Molecular Sieve 13X (washed), 60/80	350°C	50g	57730
Molecular Sieve 13X (washed), 80/100	350°C	50g	57732
3% OV™-1 on Chromosorb® W-HP, 80/100	100–350°C	20g	12716
3% OV™-17 on Chromosorb® W-HP, 100/120	20–350°C	20g	12017
3% OV™-17 on Chromosorb® W-HP, 80/100	20–350°C	20g	12719
Porapak® PS, 80/100	250°C	15g	16223
Porapak® PS, 100/120	250°C	15g	16224
Porapak® Q, 100/120	250°C	15g	16228
Silica Gel Grade 12, 60/80	350°C	100g	5651
Silica Gel Grade 12, 80/100	350°C	100g	5653
Tenax® TA, 60/80	350°C	10g	04916
Tenax® TA, 80/100	350°C	10g	04917
Tenax® GR, 60/80	350°C	10g	4937
Unibeads™ 1S, 60/80	200°C	25g	2752
Unibeads™ 2S, 60/80	200°C	25g	2760
Unibeads™ 3S, 60/80	200°C	25g	2764
Unibeads™ 3S, 80/100	200°C	25g	2766
VZN-1, 60/80	50°C	10g	8462

Custom Coated Packings

Before ordering a custom coated packing, please check the list of stock pretested packings on pages 247–248 (called "Pretested" on these pages). You may find the packing you need in stock and at a savings.

Custom Coated Packing Ordering Instructions:

1. Choose your solid support from Table 1.
2. Specify mesh range.
3. Choose the stationary phase from Table 3 or 4. It is possible to coat up to three stationary phases on one support).
4. Specify % loading for each stationary phase. Note maximum % loading guidelines for each support in Table 2.

Table 1—Base Price for Custom-Coated Packings

Solid Support	20 Grams Part No.	50 Grams Part No.	100 Grams Part No.
Chromosorb® G, W (HP)	C-3004	C-3036	C-3068
Chromosorb® G, W, or P (NAW, AW, DMCS)	C-3006	C-3038	C-3070
Chromosorb® T	C-3008	C-3040	C-3072
Chromosorb® 101-108	C-2031	C-2032	—
Chromosorb® 750	C-3012	C-3044	C-3076
Gas Chrom™ S	C-3018	C-3050	C-3082
Gas Chrom™ R	C-3020	C-3052	C-3084
Gas Chrom™ 220, 254	C-3022	C-3054	C-3086
Graphpac™-GC and -GB (15 grams), Carbograph™	C-2200	—	—
Porapak®	C-2041	C-2042	—
Tenax® or HayeSep® (15 grams)	C-3028	—	—

Table 2—Maximum Stationary Phase Loading for Solid Supports

Max. Load	Supported Material	Max. Load	Supported Material	Max. Load	Supported Material
5%	Chromosorb® G	12%	Chromosorb® 750	2%	Porapak®
15%	Chromosorb® W	15%	Gas Chrom™ S, A, P, Z, and Q	2%	Super Q
25%	Chromosorb® A	30%	Gas Chrom™ R	12%	T-Port-F
30%	Chromosorb® P	2%	Gas Chrom™ 220 and 254	2%	HayeSep®
12%	Chromosorb® T	0.5%	Glass Beads	2%	Tenax®
2%	Chromosorb® 101–108	8%	Graphpac™, Carbograph™	15%	Supelcoport

Stationary Phases

Table 3—GC Stationary Phases

Description	Solvent	Temp. Limits min.–max.	USP Code	Qty.	Part No.
Apiezon® L	Chloroform	50–300°C	—	25g	5017
bis-(2-Ethoxyethyl) Adipate (BEEA)	Chloroform	150°C	—	25g	5046
Carbowax® 20M	Chloroform	60–225°C	G16	10g	5264
DC-200, 350 cstk (Methyl)	Chloroform	20–250°C	—	100g	5251
OV™-101 (Methyl Fluid)	Chloroform	20–350°C	G1	20g	5309
SE-30 (Methyl Gum)	Chloroform	75–300°C	—	100g	5271
Sebaconitrile	Chloroform	150°C	—	20g	5244
Silar™ 10C (100% Cyanopropyl Silicone) ²	Chloroform	50–250°C	G5	5g	5250
Tetrahydroxyethylene Diamine (THEED)	Methanol	125°C	—	25g	5347

¹U.S. Patent #3,239,997, ²U.S. Patent #4,063,911.

Custom Coated Packings

Additional Stationary Phases for Custom Packings and Columns

gc columns | packed

Table 4—Stationary Phases

AT™-1000	bis-(2-Ethoxyethyl) Sebacate (BEES)	OV™-351
Apiezon® H	Glycerin	Poly-A® 103
Apiezon® J	Halocarbon 25-55 Grease	Poly-A® 135
Apiezon® K	Halocarbon Oil 14-25	Poly-I 110
Apiezon® N	Hexadecanol	Poly-S 176
BBBT	Hi-EFF™ 1AP (Diethyleneglycol Adipate)	Polyox 600M
Bentone 34	Hi-EFF™ 2AP (Ethyleneglycol Adipate)	Polyvinylpyrrolidone
BHXBT	Hi-EFF™ 2BP (Ethyleneglycol Succinate)	Reoplex 400
BMBT	Hi-EFF™ 3BP (Neopentylglycol Succinate)	SE™-52
BMEA	Hi-EFF™ 3CP (Neopentylglycol Sebacate)	SF-96
Carbowax® 540	Hi-EFF™ 8AP (Cyclohexanedimethanol Adipate)	Sorbitol
Carbowax® 600	Igepal® CO-880	SPAN 80
Carbowax® 1000	Igepal® CO-990	Squalane
Carbowax® 1450	Kel-F® Oil No. 3	SUPEROX® 4
Carbowax® 20M-TPA	Kel-F® Oil No. 10	SUPEROX® 20M
CS-10	Kel-F® Wax	TCEPE
N,N-bis-(2-Cyanoethyl) formamide (BCEF)	Krytox® 143AD	TCEP
DC-200 (500cstk)	LAC 3-R-728	β,β-Thiodipropionitrile
DC-410	LAC 4-R-886	Triton® X-100
DC-550	LAC 5-R-737	Triton® X-305
DC-702	Lexan (Polycarbonate Resin)	TWEEN 80
Di-n-butyl Maleate	Neopentylglycol Succinate	UC L-45
Didecyl Phthalate	Octyldecyl Adipate	UCON 50-HB-280-X
Di(2-ethylhexyl) Sebacate	OV™-7	UCON 75-H-90,000
Di(2-ethylhexyl) Phthalate	OV™-17 Vinyl	UCON LB-135
Diglycerol	OV™-22	UCON LB-1715
Diisodecyl Phthalate	OV™-105	UCON LB-550X
Dimethyl Sulfolane	OV™-202	UCON LB-1800-X
Dinonyl Phthalate	OV™-210	UC W-98
ECNSS-M	OV™-215	Versamid® 900
EGSS-X	OV™-225	Versilube® F-50
Ethofat 60-25	OV™-330	Versilube® F-50

Additional Charges for Custom Coated Packings

Description	Part No.
Multiple Stationary Phases (per phase)	C-2000
Expensive Stationary Phases: Add the cost of the stationary phase used to base price	—
Acetone Washing	C-2700

technical assistance

Need help selecting the right packing for your analysis?

Contact Tech Support: Phone: 1.800.255.8324 (North America)

Email: contact.alltech@grace.com

Online: www.discoverysciences.com

related products

Looking for GC packing accessories?

See pages 262–263.



GC Supports

Chromosorb® Diatomite Supports

Chromosorb® A, G, P, and W are available in both untreated and treated forms, as well as in a choice of various mesh ranges.

Chromosorb® Specifications				
Chromosorb® Series	Type	Surface Area	Packed Density	Recommended Liquid Phase Load
A	Flux-calcined	2.7m ² /g	0.48g/cc	25%
G	Flux-calcined	0.5m ² /g	0.58g/cc	5%
P	Calcined	4.0m ² /g	0.47g/cc	30%
W	Flux-calcined	1.0m ² /g	0.24g/cc	15%
750	Flux-calcined	0.75m ² /g	0.40g/cc	12%

Chromosorb® Treatment Specifications	
Abbreviation	Treatments
NAW	Untreated (non-acid washed)
AW	Acid washed (HCL is used)
DMCS	Dimethyldichlorosilane (silanized)
HP	High performance QC'ed (acid washed, silanized, flux-calcined)

Chromosorb® Diatomite Supports

Description	Qty.	Part No.
<i>Chromosorb® A-NAW</i>		
20/30	454g	2431
<i>Chromosorb® G-AW-DMCS</i>		
45/60	225g	23921

Chromosorb® Diatomite Supports (continued)

Description	Qty.	Part No.
<i>Chromosorb® P-NAW</i>		
30/60	454g	2304
45/60	454g	2305
60/80	454g	2306
80/100	454g	2307
100/120	454g	2308
<i>Chromosorb® P-AW</i>		
30/60	454g	2314
45/60	454g	2315
60/80	454g	2316
80/100	454g	2317
<i>Chromosorb® P-AW-DMCS</i>		
80/100	454g	2325
100/120	454g	2326
<i>Chromosorb® W-NAW</i>		
30/60	150g	23331
45/60	150g	23341
60/80	150g	23351
80/100	150g	23361
100/120	150g	23371
<i>Chromosorb® W-AW</i>		
30/60	150g	23451
60/80	150g	23471
100/120	150g	23491
<i>Chromosorb® W-AW-DMCS</i>		
45/60	150g	23581
60/80	150g	23591
80/100	150g	23601
100/120	150g	23611
<i>Chromosorb® W-HP</i>		
80/100	150g	23651
100/120	150g	23661
120/140	150g	24421

Chromosorb® 750

- Designed specifically for biomedical and pesticide analyses

Chromosorb® 750 is the most inert and highly efficient support material in the series. This support material is acid-washed, DMCS treated, and carefully screened.

Chromosorb® 750

Description	Qty.	Part No.
80/100	100g	2400

GC Supports and Packings

Chromosorb® T

Chromosorb® T is a support material made from PTFE and is useful for the analysis of polar compounds such as water, hydrazine, sulfur dioxide, etc.

Chromosorb® T melts at 327°C, starts to decompose at 290°C, and fuses together under prolonged use at 250°C. The density is 0.42g/cc. The surface area is 7–8m²/g and maximum loading capacity is 12%. When packing, chill the material below 19°C.

Chromosorb® T

Description	Qty.	Part No.
30/60	225g	24371

T-Port-F

T-Port-F is a screened tetrafluoroethylene support useful for analysis of very polar compounds. It has a density of about 0.5g/cc with a maximum temperature limit of 150°C.

T-Port-F

Description	Qty.	Part No.
30/60	50g	2447
80/100	50g	2449

Activated Alumina

Activated alumina is useful for the analysis of light hydrocarbons. Unsaturated hydrocarbons are retained longer than saturated ones. Grade F-6 is an indicating alumina which turns from blue to pink upon adsorption of water. The maximum temperature limit is 300°C.

Activated Alumina

Description	Qty.	Part No.
<i>Grade F-1</i>		
60/80	100g	5664
80/100	100g	5666
<i>Stock Packed Column</i>		
5ft x 1/8" SS Packed with Alumina Grade F-1, 60/80	ea	5664PC

Gas Chrom™

Gas Chrom™ Porous Polymers 220 and 254 are low in residual monomers and contaminants, condition faster and have lower column bleed.

Gas Chrom™ Specifications

	Gas Chrom™ 254	Gas Chrom™ 220
Pore Size:	3000nm	10nm
Surface Area:	1000m ² /g	450m ² /g
Equivalents:	Porapak® P	Porapak® Q

Gas Chrom™

Description	Qty.	Part No.
<i>Gas Chrom™ Porous Polymer</i>		
254, 80/100	25g	2486
<i>Stock Packed Columns</i>		
6ft x 1/8" SS Column with Gas Chrom™ Porous Polymer 220, 80/100	ea	2484PC
6ft x 1/8" SS Column with Gas Chrom™ Porous Polymer 254, 80/100	ea	2486PC

Glass Beads

Glass beads can be loaded up to 0.5 wt. % before becoming tacky. Two types are available: untreated (regular) or dimethyldichlorosilane (DMCS)-treated.

Glass Beads

Description	Qty.	Part No.
Regular, 60/80 mesh	125g	5420
Regular, 80/100 mesh	125g	5422
Regular, 100/120 mesh	125g	5424
DMCS Treated, 80/100 mesh	125g	5428
DMCS Treated, 100/120 mesh	125g	5430

GC Packings

Silica Gel

Silica Gel is used for the analysis of fixed gases and light hydrocarbons. It is also useful for dehydration of gases and liquids.

Grade 42 Indicating Silica Gel changes from blue to pink as it becomes saturated with water. It conforms to MIL-D-3716 Type IV Grade H. It will absorb water up to about 40% by weight. The maximum temperature limit is 350°C.

Silica Gel		
Description	100g Part No.	1kg Part No.
<i>Grade 12</i>		
40/60	5650	05532
60/80	5651	05533
80/100	5653	05534
100/120	5655	05535
<i>Grade 15</i>		
35/60	5648	05545
<i>Indicating Silica Gel, Grade 42</i>		
6/16	05559	05560
<i>Stock Packed Column</i>		
18ft x 1/8" SS Packed with Silica Gel, 60/80, ea		5651PC

Unibeads™ S

The Unibeads™ S series consists of spherical porous silica beads. Unibeads™ 2S are comparable to Porasil® A. Unibeads™ 3S are similar to Porasil® B. The smaller the pore size, the greater the retention for hydrocarbons.

Unibeads™ S		
Description	Qty.	Part No.
1S, 60/80, 22A	25g	2752
1S, 80/100, 22A	25g	2758
2S, 60/80, 60A	25g	2760
2S, 80/100, 60A	25g	2762
3S, 60/80, 100A	25g	2764
3S, 80/100, 100A	25g	2766
<i>Stock Packed Column</i>		
6ft x 1/8" SS Packed with Unibeads™ 1S, 60/80	ea	2752PC
6ft x 1/8" SS Packed with Unibeads™ 2S, 60/80	ea	2760PC

Chemipack® C18

Chemipack® C18 is a chemically bonded GC packing. Because it is bonded, it has high thermal stability, zero bleed and stability to water. These properties make it ideal for GC/MS or high sensitivity GC with any detector. The functional group bonded to the silica is octadecyl.

Chemipack®		
Description	Qty.	Part No.
Chemipack® C18, 80/100, 320°C	15mL	2768
<i>Stock Packed Column</i>		
6ft x 1/8" SS Packed with Chemipack® C18, 80/100, 320°C	ea	2768PC

Carbosphere™—Carbon Molecular Sieve

Carbosphere™ exhibits properties similar to molecular sieves. It has a surface area of about 1000m²/g with a pore size of about 13Å. Carbosphere™ separates oxygen, nitrogen, methane, carbon monoxide, and carbon dioxide in the same run. It also separates air, carbon dioxide, and the C1–C2 hydrocarbons. At ambient temperatures, hydrogen will elute before oxygen.

Carbosphere™ resolves trace amounts of methane and acetylene in ethylene and even trace amounts of methane, acetylene, and ethylene in ethane. The carbon molecular sieves also separate oxides of nitrogen and sulfur. Because Carbosphere™ is very nonretentive for water, rapid separations of organics, such as formaldehyde and methanol in water are possible.

It is important to remove every trace of oxygen from the carrier gas used on Carbosphere™ columns. This is essential if you are doing trace analysis or using temperatures above 200°C. An oxygen purifier such as Oxy-Trap™ or Indicating Oxy-Trap™ is ideal for removal of oxygen from carrier gas. Carbosphere™ can be affected by contaminants in the air and in carrier gas. Store under pure nitrogen or argon gas.

Carbosphere™		
Description	Qty.	Part No.
60/80	5g	5680
80/100	5g	5682
<i>Stock Packed Columns</i>		
6ft x 1/8" SS Packed with Carbosphere™, 60/80	ea	5680PC
6ft x 1/8" SS Packed with Carbosphere™, 80/100	ea	5682PC
10ft x 1/8" SS Packed with Carbosphere™, 80/100	ea	56821PC

related products
Need an oxygen trap?
See page 271.



GC Packings

Molecular Sieves

- For drying applications

Molecular sieves are synthetic alkali metal aluminosilicates with various cations. They are used for the separation of fixed gases and drying of liquid or gas streams. They can be reactivated by heating at 250°C for 12 hours; 300°C for 4 hours; or 350°C for 2 hours.

Molecular Sieves		
Description	100g Part No.	1kg Part No.
<i>Molecular Sieve 3A</i>		
1/8" Pellets	05252	—
1/16" Pellets	05254	—
40/60	05306	—
60/80	05307	—
80/100	05308	—
100/120	05309	—
<i>Molecular Sieve 4A</i>		
8/12 Beads	05256	—
1/8" Pellets	05258	05259
1/16" Pellets	05260	05261
40/60	5620	—
60/80	5622	—
80/100	5624	—
100/120	5625	—
<i>Molecular Sieve 5A</i>		
1/8" Pellets	05264	—
1/16" Pellets	5633	—
80/100	5630	—
100/120	5632	—
<i>Molecular Sieve 13X</i>		
8/12 Beads	05268	—
1/8" Pellets	05270	05271
1/16" Pellets	05272	05273
40/60	5634	—
60/80	5636	—
80/100	5638	—

Washed Molecular Sieves

- For chromatography applications

The usual gas chromatographic grades of molecular sieves have a fine dust adhering to each particle which cannot be removed simply by careful sieving. This dust is responsible for poorer separations and higher pressure drops than molecular sieves which are essentially dust-free.

The washed molecular sieves offered here have been carefully cleaned of dust particles by washing with distilled water. The cleaned product is then activated and sealed in glass jars. Molecular sieves are typically reactivated at a temperature of 300°C for four hours.

Washed Molecular Sieves		
Description	Qty.	Part No.
<i>Washed Molecular Sieve 5A</i>		
40/60	50g	5602
60/80	50g	5604
80/100	50g	5605
100/120	50g	5606
<i>Washed Molecular Sieve 13X</i>		
40/60	50g	57728
60/80	50g	57730
80/100	50g	57732
100/120	50g	57734
<i>Stock Packed Columns</i>		
6ft x 1/8" SS Packed with 80/100 Washed 5A	ea	5605PC
6ft x 1/8" SS Packed with 80/100 Washed 13X	ea	5773PC

technical assistance

Contact Tech Support: Phone: 1.800.255.8324 (North America)
 Email: contact.alltech@grace.com
 Online: www.discoverysciences.com

Graphitized Carbon Blacks

Carbograph™ Packings

- High speed analysis
- Extremely hydrophobic—ideal for aqueous injections

Equivalent to the Carbopack™ packings. Carbograph™ packings can be modified with a stationary phase giving them unique selectivities.

Carbograph™ 1

Equivalent to Carbopack™ B. It has a surface area of 100m²/g. Carbograph™ 1-TD has a coarse mesh of 20/40 which is suitable for thermal desorption work. Maximum temperature 500°C.

Carbograph™ 1SC

Similar to Carbopack™ BHT, specifically for the separation of sulfur compounds. It has a mesh of 40/60, and a maximum temperature of 225°C.

Carbograph™ 2

Equivalent to Carbopack™ C. It has a specific surface of 10m²/g for the separation of volatile organic compounds. Carbograph™ 2-TD has a coarse mesh of 20/40 which is suitable for thermal desorption work. Maximum temperature 500°C.

Carbograph™

Description	Qty.	Part No.
<i>Uncoated</i>		
Carbograph™ 1-TD, 20/40, >500°C	15g	1740
Carbograph™ 1, 60/80, >500°C	10g	1722
Carbograph™ 1 80/120, >500°C	10g	1724
Carbograph™ 2-TD, 20/40, >500°C	15g	1745
Carbograph™ 2, 60/80, >500°C	10g	1726
Carbograph™ 2, 80/100, >500°C	10g	1728
Carbograph™ 1-SC, 40/60, 225°C	10g	1734
<i>Coated</i>		
1% AT™-1000 on Carbograph™ 1 60/80, 225°C	15g	1733
0.1% AT™-1000 on Carbograph™ 2 80/100, 225°C	15g	1737
4% Carbowax® 20M on Carbograph™ 1 DA 80/120, 200°C	15g	17271
0.2% Carbowax® 1500 on Carbograph™ 2 60/80, 175°C	15g	1725
0.2% Carbowax® 1500 on Carbograph™ 2 80/100, 175°C	15g	1720
5% Carbowax® 20M on Carbograph™ 1 AW 80/120, 225°C	15g	1729
4% Carbowax® 20M + 0.8% KOH on Carbograph™ 1 60/80, 220°C	15g	1735

Graphpac™ -GC and -GB Packings

- Graphitized carbon
- Inert, non-specific surface
- Unique separation properties
- Optimized phase loading
- Hydrophobic—ideal for aqueous injections

Graphpac™-GC has an average surface area of 10–13m²/g and interacts with samples almost completely by non-specific adsorption, which results in unique separations not achievable with fused silica open tubular (FSOT) capillary columns. Alcohols, hydrocarbons, and acids are separated in series according to carbon number.

The carbons possess partitioning capabilities as uncoated packings, use stationary phases to modify the selectivity of these carbons.

Graphpac™

Description	Qty.	Part No.
<i>Graphpac™-GC</i>		
Uncoated, 60/80	15g	8536
Uncoated, 80/100	15g	8538
Custom—Specify Phase and % Loading	15g	C-2200
<i>Stock Packed Columns</i>		
6ft x 1/8" SS Packed with 0.2% Carbowax® 1500, 175°C, 80/100	ea	8546PC

Porous Polymers

Tenax®-TA (GC) Polymers

- Stable baseline after conditioning
- Short retention times
- Low affinity for water
- Maximum temperature limit is 350°C

Tenax®-TA is a porous polymer that is based on 2,6-Diphenyl-p-phenylene Oxide. It has replaced Tenax®-GC. Tenax®-TA can be used as a packing and as a trapping material. Both the EPA and NIOSH specify the use of Tenax® in their standard methods. Tenax® is particularly useful for the analysis of high boiling compounds such as alcohols, polyethylene glycols, diols, phenols, monoamines and diamines, ethanolamines, aldehydes, ketones, and chlorinated aromatics.

Tenax®-TA (GC) Specifications	
Specific Surface Area:	35m ² /g
Pore Volume:	2.4g/cc
Average Pore Size:	200nm
Density:	0.25g/cc

Tenax®-TA (GC)

Description	Qty.	Part No.
20/35	10g	04914
20/35	100g	049141
35/60	10g	04915
35/60	100g	049151
60/80	10g	04916
60/80	100g	049161
60/80	500g	049162
80/100	10g	04917
80/100	100g	049171

Pretreated Tenax®-TA

Cleaned by Solvent Extraction	—	C-4196
-------------------------------	---	---------------

Stocked Packed Columns

6ft x 1/8" SS Packed with Tenax®-TA Porous Polymer, 60/80	ea	4900PC
6ft x 1/8" SS Packed with Tenax®-TA Porous Polymer, 80/100	ea	4901PC

Tenax®-GR Polymers

- Superior properties as trapping material
- Maximum temperature limit is 350°C

Tenax®-GR contains 23% graphitized carbon as an integral part of the material. This is not an admixture; the graphitized carbon is co-precipitated with the Tenax® polymer. The resulting material gives higher breakthrough volumes for most materials, yet is less water retentive than Tenax®-TA. When using this material for packing GC columns, better peak symmetry is observed. Surface area is 24.1m²/g.

Tenax®-GR

Description	Qty.	Part No.
20/35	10g	4933
20/35	100g	49331
35/60	10g	4935
60/80	10g	4937
60/80	500g	49372
80/100	10g	4939

tech tip

Conditioning

We highly recommend that Tenax®, like all other porous polymers, be preconditioned before being packed into GC columns. This removes most of the shrinkage of the polymer and minimizes the need for further conditioning. Tenax® can be expected to shrink further in use.

GC Packings

HayeSep® Porous Polymer

- Pre-purified porous polymers for GC
- Less bleed and shrinkage than other porous polymers
- Cleaned and conditioned to 200°C (165°C for N & T)

HayeSep® polymers were developed to solve room temperature separation of N₂, O₂, Ar and CO; ppm H₂O in the presence of HCl, Cl₂ etc. The HayeSep® D series better separates the fixed gases to the C₂ hydrocarbons.

HayeSep® Specifications			
Series	Polymer Composition	Surface Area	Polarity
Q	DVB*	582m ² /g	Non-Polar
P	DVB*/Styrene	165m ² /g	Slightly Polar
R	DVB*/N-Vinyl-2-pyrrolidinone	344m ² /g	Moderately Polar
T	EGDM**	250m ² /g	Very Polar
N	DVB*/EGDM**	405m ² /g	Very Polar
A	DVB*/EGDM** (high purity)	526m ² /g	Moderately Polar
D	DVB* (high purity)	795m ² /g	Non-Polar

*Divinylbenzene. **Ethyleneglycoldimethacrylate.

Comparison of HayeSep® “D” Formulations†

Table 1—Retention Time (in minutes)							
	Air	CH ₄	CO ₂	C ₂ H ₂	C ₂ H ₄	C ₂ H ₆	H ₂ O
D	0.9	1.7	3.1	5.8	6.1	8.4	8.6
DOX	1.0	1.8	3.2	5.9	6.3	8.5	8.3
DB	0.9	1.6	3.1	6.1	6.6	8.7	8.1

†Retention data based on 10ft x 1/8" stainless steel at 45°C, 30mL/min.

HayeSep®

Description	Mesh	Temp. Limit max.	Qty.	Part No.
HayeSep® A	80/100	165°C	75cc	2819
	100/120	165°C	75cc	2820
	20/140	165°C	75cc	2821
HayeSep® D	60/80	290°C	75cc	2828
	80/100	290°C	75cc	2829
	100/120	290°C	75cc	2830

HayeSep® (continued)

Description	Mesh	Temp. Limit max.	Qty.	Part No.
HayeSep® DB	60/80	290°C	75cc	2834
	100/120	290°C	75cc	2836
HayeSep® N	60/80	165°C	75cc	2815
	80/100	165°C	75cc	2816
	100/120	165°C	75cc	2817
HayeSep® P	60/80	250°C	75cc	2803
HayeSep® Q	60/80	275°C	75cc	2800
	80/100	275°C	75cc	2801
HayeSep® R	60/80	250°C	75cc	2806
	80/100	250°C	75cc	2807
HayeSep® T	80/100	165°C	75cc	2813
<i>Stock Packed Columns</i>				
10ft x 1/8" SS Packed with HayeSep® D	100/120	290°C	ea	28301PC
20ft x 1/8" NI Packed with HayeSep® D	100/120	290°C	ea	27082PC
30ft x 1/8" SS Packed with HayeSep® D	100/120	290°C	ea	27083PC
30ft x 1/8" SS Packed with HayeSep® DB	100/120	290°C	ea	2836PC
8ft x 1/8" SS Packed with HayeSep® P	60/80	250°C	ea	2803PC
6ft x 1/8" SS Packed with HayeSep® P	80/100	250°C	75cc	2804PC
6ft x 1/8" SS Packed with HayeSep® Q	80/100	275°C	ea	2801PC
8ft x 1/8" SS Packed with HayeSep® Q	80/100	275°C	ea	28010PC
6ft x 1/8" SS Packed with HayeSep® T	80/100	165°C	ea	2813PC

GC Packings

Century Series Chromosorb®

Century Series Chromosorb® porous polymers have a rigid structure and distinct pore size. They are packed into columns in the normal manner and do not require a liquid coating. We recommend preconditioning all porous polymers before packing columns.

Chromosorb® Specifications			
Series	Polymer Composition	Surface Area	Polarity
101	DVB*/Styrene	>50m ² /g	Non-Polar
102	DVB*/Styrene	300–400m ² /g	Slightly Polar
103	Cross-linked Polystyrene**	15–25m ² /g	Non-Polar
106	Cross-linked Polystyrene	700–800m ² /g	Non-Polar
107	Cross-linked Acrylic Ester	400–500m ² /g	Polar

*Divinylbenzene. **The surface is basic.

Chromosorb®

Description	Mesh	Temp. Limit* max.	Qty.	Part No.
<i>Chromosorb® 101</i>	60/80	275/325°C	50g	2405
<i>Chromosorb® 102</i>	60/80	250/300°C	50g	2408
	80/100	250/300°C	50g	2409
<i>Chromosorb® 103</i>	80/100	250/300°C	50g	2412
<i>Chromosorb® 106</i>	60/80	250/275°C	50g	2420
<i>Chromosorb® 107</i>	80/100	250/275°C	50g	2424

*Isothermal/Temp. program.

Porapak®

Porapak® GC packings are cross-linked polymers which can be used directly in GC columns without a stationary phase coating. Acetone washing of Porapaks® improves performance.

Porapak® Specifications			
Series	Polymer Composition	Surface Area	Polarity
Q	DVB*/Ethylvinylbenzene	500–600m ² /g	Slightly Polar
P	DVB*/Styrene	100–200m ² /g	Non-Polar
R	DVB*/Vinyl pyrrolidinone	450–600m ² /g	Mod. Polar
S	DVB*/Vinyl pyridine	300–450m ² /g	Mod. Polar
T	EGDM**	225–350m ² /g	Polar
N	DVB*/Vinyl pyrrolidinone	250–350m ² /g	Very Polar

*Divinylbenzene. **Ethylene glycol dimethacrylate.

Porapak®

Description	Mesh	Temp. Limit max.	Qty.	Part No.
<i>Porapak® Q</i>	50/80	250°C	26g	2700
	80/100	250°C	26g	2701
	100/120	250°C	26g	2702
<i>Porapak® P</i>	50/80	250°C	20g	2703
	80/100	250°C	20g	2704
	100/120	250°C	20g	2705
<i>Porapak® R</i>	50/80	250°C	24g	2706
	80/100	250°C	24g	2707
	100/120	250°C	24g	2708
<i>Porapak® S</i>	50/80	250°C	26g	2709
	80/100	250°C	26g	2710
<i>Porapak® T</i>	50/80	190°C	31g	2712
	80/100	190°C	31g	2713
<i>Porapak® N</i>	50/80	190°C	29g	2715
	80/100	190°C	29g	2716
	100/120	190°C	29g	2717
<i>Porapak® PS</i>	80/100	250°C	20g	2722
	100/120	250°C	20g	2723

tech tip

Temperature limits

Our GC capillary columns are temperature rated. In some cases, we list two maximum operating temperatures, the lower one is for isothermal condition and the higher one for temperature-programmed condition.