


SOLVENT ANALYSIS ON A 0.53 MM ID BP624

SOLVENT ANALYSIS

Column Part No.: 054835
Phase: BP624, 3.0 μm
Column: 50 m x 0.53 mm ID
Initial Temp.: 35 $^{\circ}\text{C}$, 6 min
Rate: 15 $^{\circ}\text{C}/\text{min}$

Final Temp.: 160 $^{\circ}\text{C}$, 2 min
Detector: FID 280 $^{\circ}\text{C}$
Injection Mode: Split (30:1)
Carrier Gas: He, 4 psi




For more information contact our technical customer support team on: techsupport@sge.com

ACRYLIC ACID/ACRYLAMIDE ANALYSIS ON BP21

ACRYLIC ACID/ACRYLAMIDE ANALYSIS

Column Part No.:	054473	Final Temp:	150 °C
Phase:	BP21, 0.5 µm film	Detector:	FID, 280 °C
Column:	12 m x 0.53 mm ID	Injection Mode:	On-Column
Initial Temp:	75 °C, 0.5 min	Carrier Gas:	He, 6 psi
Rate:	10 °C/min		

Notes: When response of acrylic acid does occur, removal of 30 cm from the front of the column will correct this loss. On-column injection is recommended or polymerisation of acrylic acid may occur.




For more information contact our technical customer support team on: techsupport@sge.com

VOLATILE SOLVENT ANALYSIS ON BP624

Column Part No.: 054834
Phase: BP624
Column: 25 m x 0.53 mm ID
Initial Temp.: 40 °C, 2 min
Rate 1: 7 °C/min
Final Temp.: 100 °C


Detector: FID, 280 °C
Injector Mode: Split 30:1
Carrier Gas: He, 2.5 psi
Injection Volume: 0.5 µl
Injection Temp: 200 °C



For more information contact our technical customer support team on: techsupport@sge.com

ETHYLENE GLYCOL AND DIETHYLENE GLYCOL ANALYSIS ON BP21


Column Part No.:	054473	Detector:	FID, 270 °C
Phase:	BP21, 0.5 µm	Injector Mode:	Split 20:1
Column:	12 m x 0.53 mm ID	Carrier Gas:	He, 2 psi
Initial Temp.:	110 °C, 0.5 min	Injection Volume:	0.4 µl
Rate:	15 °C/min	Injection Temp:	240 °C
Final Temp.:	170 °C		



For more information contact our technical customer support team on: techsupport@sge.com

TRIACETIN ANALYSIS ON BP21


Column Part No.:	054473	Detector:	FID, 260 °C
Phase:	BP21, 0.5 µm	Injector Mode:	Split 20:1
Column:	12 m x 0.53 mm ID	Carrier Gas:	He, 2 psi
Initial Temp.:	100 °C, 1 min	Injection Volume:	0.2 µl
Rate:	15 °C/min	Injection Temp:	200 °C
Final Temp.:	170 °C, 2 min		



For more information contact our technical customer support team on: techsupport@sge.com

ETHANOL IMPURITY IN METHANOL ON BP1 ON A 0.53 MM ID COLUMN


Column Part No.:	054095	Detector:	FID, 300 °C
Phase:	BPI, 5.0 µm	Injector Mode:	Split, 20:1
Column:	25 m x 0.53 mm ID	Carrier Gas:	He, 4 psi
Temp:	50 °C, Isothermal	Injection Volume:	1 µl



For more information contact our technical customer support team on: techsupport@sge.com

1,2-BUTADIENE AND 1,3-BUTADIENE IN HEXANE ON THICK FILM BP1


Column Part No.:	054081	Final Temp.:	140 °C, 1 min
Phase:	BP1, 5.0 µm	Detector:	FID, 300 °C
Column:	25 m x 0.32 mm ID	Injector Mode:	Split, 100:1, 240 °C
Initial Temp.:	40 °C, 2 min	Carrier Gas:	He, 6 psi
Rate 1:	10 °C/min	Injection Volume:	1 µL



For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF A STANDARD SOLVENT MIXTURE ON BP20


Column Part No.:	054442	Final Temp.:	90 °C
Phase:	BP20, 1.0 µm	Detector:	FID, 260 °C
Column:	25 m x 0.32 mm ID	Injector Mode:	Split, 100:1, 240 °C
Initial Temp.:	45 °C, 2 min	Carrier Gas:	Hydrogen, 80 kPa
Rate 1:	5 °C/min	Injection Volume:	0.1 µL



For more information contact our technical customer support team on: techsupport@sge.com

SOLVENT IMPURITIES IN DIOXANE (10 PPM) ON BP20


Column Part No.:	054448	Final Temp.:	120 °C
Phase:	BP20, 1.0 µm	Detector:	FID, 280 °C
Column:	25 m x 0.53 mm ID	Injector Mode:	Split, 3:1, 200 °C
Initial Temp.:	40 °C	Carrier Gas:	Hydrogen, 2 psi
Rate:	10 °C/min	Injection Volume:	1.0 µL



For more information contact our technical customer support team on: techsupport@sge.com

SOLVENT IMPURITIES IN DIOXANE (10 PPM) ON BP20


Column Part No.:	054448	Final Temp.:	120 °C
Phase:	BP20, 1.0 µm	Detector:	FID, 280 °C
Column:	25 m x 0.53 mm ID	Injector Mode:	Split, 30:1,
Initial Temp.:	40 °C, 2 min	Carrier Gas:	Hydrogen, 2 psi
Rate:	10 °C/min	Injection Volume:	0.2 µL



For more information contact our technical customer support team on: techsupport@sge.com

SOLVENT MIXTURE ON BP20


Column Part No.:	054436	Final Temp.:	140 °C
Phase:	BP20, 0.5 µm	Detector:	FID, 280 °C
Column:	25 m x 0.32 mm ID	Injector Mode:	Split, 100:1, 240 °C
Initial Temp.:	60 °C, 3.5 min	Carrier Gas:	He, 6 psi
Rate:	15 °C/min	Injection Volume:	1 µL



For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF DIETHYLAMINO ETHOXYETHANOL ON BP20


Column Part No.:	054448	Final Temp.:	180 °C, 1 min
Phase:	BP20, 1.0 µm	Detector:	FID, 280 °C
Column:	25 m x 0.53 mm ID	Injector Mode:	Split, 50:1, 240 °C
Initial Temp.:	150 °C, 1 min	Carrier Gas:	He, 15 kPa
Rate:	10 °C/min	Injection Volume:	0.1 µL



For more information contact our technical customer support team on: techsupport@sge.com

PETROLEUM HYDROCARBONS ON A THICK FILM 0.53 MM ID BP1


Column Part No.:	054095	Detector:	FID, 320 °C
Phase:	BP1, 5.0 µm	Injector Mode:	Split, 100:1, 200 °C
Column:	25 m x 0.53 mm ID	Carrier Gas:	He, 8 mL/min
Initial Temp.:	40 °C isothermal	Injection Volume:	0.1 µL



For more information contact our technical customer support team on: techsupport@sge.com

BIS-ETHYLENE STERAMIDE SAMPLE ON BPX5


Column Part No.:	054118	Final Temp.:	370 °C, 10 min
Phase:	BPX5, 0.25 µm	Detector:	FID, 400 °C
Column:	12 m x 0.32 mm ID	Injector Mode:	On-Column (OCI-5)
Initial Temp.:	60 °C	Carrier Gas:	He, 10 psi
Rate 1:	10 °C/min	Injection Volume:	1 µL



For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF A COMMON SOLVENT MIXTURE ON BP624

Column Part No.:	054832	Carrier Gas Flow:	2.2 mL/min.
Phase:	BP624, 1.8 μ m film	Constant Flow:	On
Alcohol mix:	1000 ppm in Dimethyl Sulfoxide	Average Linear Velocity:	34 cm/sec at 32 $^{\circ}$ C
Column:	30 m x 0.32 mm ID	Injection Mode:	Split
Initial Temp:	32 $^{\circ}$ C, 9 min.	Split Ratio:	100:1
Rate:	30 $^{\circ}$ C/min to 190 $^{\circ}$ C	Injection Volume:	0.2 μ L
Final Temp:	190 $^{\circ}$ C, 0 min.	Injection Temperature:	250 $^{\circ}$ C
Detector Type:	FID	Autosampler:	No
Carrier Gas:	He, 9.6 psi	Liner Type:	4 mm ID Double Taper Liner
		Liner Part Number:	092018



For more information contact our technical customer support team on: techsupport@sge.com

CAPILLARY COLUMN TEST MIX ON SOLGEL-1MS™

Column Part No.:	054795	Injection Mode:	Splitless
Phase:	SolGel-1ms, 0.25 µm film	Purge on Time:	0.5 min.
Sample:	1000 ppm in Pentane	Purge on (Split)	
Column:	30 m x 0.25 mm ID	Vent Flow:	60 mL/min.
Isothermal Temp.:	145 °C, 9 min.	Injection Volume:	1 µL
Detector Type:	Mass Spectrometer	Injection Temp.:	250 °C
Carrier Gas:	He, 34.5 psi	Liner Type:	4 mm ID Double Taper Liner
Carrier Gas Flow:	1.6 mL/min.	Liner Part Number:	092017
Constant Flow:	On	Full Scan / SIM:	Full scan 45-450
Average Linear Velocity:	35 cm/sec at 50 °C		




Table 1.

Kovats Index comparison of SolGel-1 with BP1 and a competitors 100 % Dimethyl Polysiloxane.

No.	Component	Kovats Index		
		SolGel-1	BP1	Competitors 100 % Dimethyl Polysiloxane
1	Pentane			
2	Decane	1000	1000	1000
3	4-Chlorophenol	1165	1167	1164
4	Decylamine	1240	1241	1239
5	Undecanol	1357	1358	1356
6	Acenaphthylene	1433	1434	1434
7	Pentadecane	1500	1500	1500


For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF 22 KETONES ON BPX35

Column Part No.: 054701

Phase: BPX35, 0.25 µm film
 Sample: 300 ppm in dichloromethane
 Column: 30 m x 0.25 mm ID
 Initial Temp: 40 °C, 5 min.
 Rate: 10 °C/min to 170 °C,
 Final Temp: 170 °C, 5 min.
 Detector Type: Mass Spectrometer
 Carrier Gas: He, 25.6 psi
 Carrier Gas Flow: 1.6 mL/min.

Constant Flow: On
 Average
 Linear Velocity: 35 cm/sec at 40 °C
 Injection Mode: Split
 Split Ratio: 80:1
 Injection Volume: 0.5 µL
 Injection Temp.: 250 °C
 Liner Type: 4 mm ID Single Taper Liner
 Liner Part Number: 092017
 Full Scan / SIM: Full scan 45-450




For more information contact our technical customer support team on: techsupport@sge.com

GROB TEST MIX ON SOLGEL-WAX™

Column Part No.: 054796

Phase: SolGel-WAX, 0.25 µm film
 Sample: 200 ppm in dichloromethane
 Column: 30 m x 0.25 mm ID
 Initial Temp.: 40°C, 1 min.
 Rate 1: 6 °C/min to 160 °C,
 Final Temp.: 160 °C, 5 min.
 Detector Type: Mass Spectrometer
 Carrier Gas: He, 25.7 psi
 Carrier Gas Flow: 1.8 mL/min.

Constant Flow: On
 Average
 Linear Velocity: 35 cm/sec at 40 °C
 Injection Mode: Split
 Split Ratio: 50:1
 Injection Volume: 1 µL
 Injection Temp.: 250 °C
 Liner Type: 4 mm ID Single Taper Liner
 Liner Part Number: 092017
 Full Scan / SIM: Full scan 45-450




For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF CIGARETTE LIGHTER FUEL ON A 6 METRE BPX5 COLUMN

Column Part No.: 0541261

Phase: BPX5, 1.0 μ m film
Sample: Butane mixture from lighter
Column: 6 m x 0.32 mm ID
Isothermal Temp.: 30 $^{\circ}$ C, 5 min.
Detector Type: FID
Detector Temp.: 270 $^{\circ}$ C
Carrier Gas: He, 1.94 psi
Carrier Gas Flow: 1.7 mL/min.

Constant Flow: On
Average
Linear Velocity: 35 cm/sec at 30 $^{\circ}$ C
Injection Mode: Split
Split Ratio: 100:1
Injection Volume: 10 μ L
Injection Temp.: 250 $^{\circ}$ C
Liner Type: 4 mm ID Single Taper Liner
Liner Part Number: 092017




For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF GLYCEROL IN BIODIESEL ON SOLGEL-WAX

Column Part No.: 054796

Phase: SolGel-WAX, 0.25 μm film
 Sample: 500 ppm glycerol and Biodiesel
 Column: 30 m x 0.25 mm ID
 Instrument: HP5890
 Initial Temp: 50 $^{\circ}\text{C}$, 2 min.
 Rate 1: 10 $^{\circ}\text{C}/\text{min}$ to 280 $^{\circ}\text{C}$,
 Final Temp: 280 $^{\circ}\text{C}$, 16 min.
 Detector Type: FID
 Detector Temp.: 300 $^{\circ}\text{C}$

Carrier Gas: He, 14 psi
 Carrier Gas Flow: 1.2 mL/min.
 Constant Flow: On
 Average
 Linear Velocity: 29 cm/sec at 50 $^{\circ}\text{C}$
 Injection Mode: Split
 Injection Volume: 1 μL
 Injection Temp.: 250 $^{\circ}\text{C}$
 ATAS Optic 2-200 programmable injector
 Focus Autosampling Robot



For more information contact our technical customer support team on: techsupport@sge.com

Phase:	SolGel-WAX, 0.25µm film	Optic conditions	
Sample:	ethylene glycol in HPLC grade water (standards)	Initial Temp:	60°C
Column:	30m x 0.25 mm ID	Ramp Rate:	16°C/min
Instrument:	Agilent 5890	Final temp:	250°C
Injector	Optic 3 Programmable Injector	Injection Volume:	2 µL
Initial Temp:	100°C, 2 min.	Liner Type:	4 mm ID Fritted
Rate 1:	10°C/min to 230°C,	Transfer column flow:	2.3mL/min
Final Temp:	230°C	Transfer time:	2 min.
Detector Type:	FID	Initial column flow:	2.3 mL/min
Detector temp:	250°C	Column flow:	2.3 mL/min
Carrier Gas:	He,	Split flow:	100 mL/min
		Column Part Number:	054796




Figure 1: 2µL cold split injection of 5ppm ethylene glycol (1) standard




Figure 2: 2µL cold split injection of waste water: (1) ethylene glycol; (2) diethylene glycol

ACKNOWLEDGEMENT

SGE would like to thank ATAS for supplying the chromatograms and data for this application note.

visit www.sge.com or contact our technical Customer Support Team on: techsupport@sge.com

SGE International Pty. Ltd.
email: support@sge.com

SGE, Incorporated (USA)
email: usa@sge.com

SGE (Italia) Srl.
email: sge.italia@tin.it

SGE China Service Centre
email: sge@bj.col.com.cn

SGE India
email: sgeindia@vsnl.com

SGE Europe Ltd. (UK)
email: uk@sge.com

SGE (France) Sarl
email: france@sge.com

SGE Deutschland GmbH
email: germany@sge.com


SGE Japan Inc.
email: japan@sge.com

ANALYSIS OF KETONES ON THICK FILM BPX5

Column Part No.: 054123

Phase: BPX5, 1.0 μm film
 Sample: 300 ppm in dichloromethane
 Column: 60 m x 0.25 mm ID
 Initial Temp.: 40 $^{\circ}\text{C}$, 5 min.
 Rate 1: 1 0 $^{\circ}\text{C}/\text{min}$ to 80 $^{\circ}\text{C}$
 Rate 2: 30 $^{\circ}\text{C}/\text{min}$ to 260 $^{\circ}\text{C}$
 Final Temp: 260 $^{\circ}\text{C}$, 4 min.
 Detector Type: FID
 Detector Temp.: 360 $^{\circ}\text{C}$
 Carrier Gas: He, 27.6 psi

Carrier Gas Flow: 1.9 mL/min.
 Constant Flow: On
 Average
 Linear Velocity: 35 cm/sec at 40 $^{\circ}\text{C}$
 Injection Mode: Split
 Split Ratio: 100:1
 Injection Volume: 0.4 μL
 Injection Tem.: 250 $^{\circ}\text{C}$
 Liner Type: 4 mm ID
 Single Taper Liner
 Liner Part Number: 092017




For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF A COMMON SOLVENT MIXTURE ON A THICK FILM BPX5

Column Part No.: 054123

Phase: BPX5, 1.0 μ m film
Sample: neat
Column: 60 m x 0.25 mm ID
Initial Temp.: 32 $^{\circ}$ C, 5 min.
Rate 1: 20 $^{\circ}$ C/min to 190 $^{\circ}$ C,
Final Temp: 190 $^{\circ}$ C, 2 min.
Detector Type: FID
Detector Temp.: 360 $^{\circ}$ C
Carrier Gas: He, 26.9 psi

Carrier Gas Flow: 1.9 mL/min.
Constant Flow: On
Average
Linear Velocity: 35 cm/sec at 40 $^{\circ}$ C
Injection Mode: Split
Split Ratio: 100:1
Injection Volume: 0.3 μ L
Injection Temp.: 250 $^{\circ}$ C
Liner Type: 4 mm ID Single Taper Liner
Liner Part Number: 092017

**Components**

1. Methanol
2. Ethanol
3. Acetone
4. Isopropyl alcohol
5. Acetonitrile
6. Dichloromethane
7. 2-Butanone
8. Ethyl acetate
9. Benzene

For more information contact our technical customer support team on: techsupport@sge.com