

# Inertsil<sup>®</sup> C8-4


## Physical Properties

- Silica : 3 Series High Purity Silica Gel
- Particle Size : 2 μm, 3 μm, 5 μm
- Surface Area : 450 m<sup>2</sup>/g
- Pore Size : 100 Å (10 nm)
- Pore Volume : 1.05 mL/g
- Bonded Phase : Octyl Groups
- End-capping : Yes
- Carbon Loading : 5 %
- USP Code : L7
- pH Range : 2 ~ 7.5



Many chromatographers prefer a C8 column when an ODS phase shows excessive retention values. Inertsil C8-4 provides the same separation pattern (selectivity) and extreme inertness to any type of compounds just like Inertsil ODS-4, which enables easy method transfer from ODS-4 to C8-4 while other commercially available ODS and C8 columns can show dramatically different selectivity even though they are part of the same brand/series.

**Figure 1 : Comparison of Separation Pattern between ODS and C8 Columns**



### Conditions

Column Size : 5 μm, 250 × 4.6 mm I.D.	1. Uracil	(0.005 mg/mL)
Eluent : A) CH <sub>3</sub> OH	2. Caffeine	(0.04 mg/mL)
B) H <sub>2</sub> O	3. Phenol	(0.08 mg/mL)
A / B = 80/20, v/v	4. <i>n</i> -Butylbenzene	(1.12 mg/mL)
Flow Rate : 1.0 mL/min	5. <i>o</i> -Terphenyl	(0.04 mg/mL)
Col. Temp. : 40 °C	6. <i>n</i> -Amylbenzene	(1.37 mg/mL)
Detection : UV 254 nm	7. Triphenylene	(0.014 mg/mL)
Injection Vol. : 5 μL		

## Analytical Columns

Particle Size: 2 µm	Length \ I.D. (mm)	2.1	3.0
		30	5020-81280
	50	5020-81282	5020-81292
	75	5020-81283	5020-81293
	100	5020-81284	5020-81294
	150	5020-81285	5020-81295

HPSeries Particle Size: 3 µm 50 MPa (500 Bar)	Length \ I.D. (mm)	2.1	3.0	4.6
	30	5020-14071	5020-14074	5020-14077
	50	5020-14072	5020-14075	5020-14078
	75	5020-14073	5020-14076	5020-14079
	100	5020-14051	5020-14054	5020-14057
	150	5020-14052	5020-14055	5020-14058
	250	5020-14053	5020-14056	5020-14059

\* End-fittings are 1/16" Waters-compatible.  
 \* UHPLC compatible end-fittings are also available upon request for UHPLC systems (Ex: UPLC) to avoid dead volume.  
 \* Indicate "UP Type end-fittings" when ordering. (Please note that UP type is not available for a 4.6 mm I.D. column)  
 \* For maximum operating pressure information, please refer to page 46.

Particle Size: 3 µm	Length \ I.D. (mm)	1.0	1.5
	30	5020-81261	5020-81271
	50	5020-81262	5020-81272
	75	5020-81263	5020-81273
	100	5020-81264	5020-81274
	150	5020-81265	5020-81275
	250	5020-81266	5020-81276

Particle Size: 3 µm	Length \ I.D. (mm)	2.1	3.0	4.0	4.6
	30	5020-03971	5020-03978	5020-03985	5020-03992
	50	5020-03972	5020-03979	5020-03986	5020-03993
	75	5020-03973	5020-03980	5020-03987	5020-03994
	100	5020-03974	5020-03981	5020-03988	5020-03995
	125	5020-03977	5020-03984	5020-03991	5020-03998
	150	5020-03975	5020-03982	5020-03989	5020-03996
	250	5020-03976	5020-03983	5020-03990	5020-03997

Particle Size: 5 µm	Length \ I.D. (mm)	1.0	1.5
	30	5020-81221	5020-81231
	50	5020-81222	5020-81232
	75	5020-81223	5020-81233
	100	5020-81224	5020-81234
	150	5020-81225	5020-81235
	250	5020-81226	5020-81236

Particle Size: 5 µm	Length \ I.D. (mm)	2.1	3.0	4.0	4.6
	30	5020-04051	5020-04061	5020-04071	5020-04081
	50	5020-04052	5020-04062	5020-04072	5020-04082
	75	5020-04053	5020-04063	5020-04073	5020-04083
	100	5020-04054	5020-04064	5020-04074	5020-04084
	125	5020-04057	5020-04067	5020-04077	5020-04087
	150	5020-04055	5020-04065	5020-04075	5020-04085
	250	5020-04056	5020-04066	5020-04076	5020-04086

\* End-fittings are 1/16" Waters-compatible.  
 \* For maximum operating pressure information, please refer to page 46.

## Cartridge Guard Column E

I.D. of the Analytical Column Applicable (mm)	Length (mm)	I.D. (mm)	Replacement Cartridge E Guard Column (2 EA.)		Cartridge E Holder / Cartridge Set (2 Cartridge E Guard Columns & 1 Holder)	
			Particle Size		Particle Size	
			3 µm	5 µm	3 µm	5 µm
1.0	10	1.0	5020-19247	5020-19246	5020-19297	5020-19296
1.5, 2.1		1.5	5020-19347	5020-19346	5020-19397	5020-19396
2.1, 3.0		3.0	5020-19147	5020-19146	5020-19197	5020-19196
4.0, 4.6		4.0	5020-19047	5020-19046	5020-19097	5020-19096
2.1, 3.0	20	3.0	5020-19547	5020-19546	5020-19597	5020-19596
4.0, 4.6		4.0	5020-19447	5020-19446	5020-19497	5020-19496
Holder for Cartridge Guard Column E			For 10 mm Length		5020-08500	
			For 20 mm Length		5020-08550	

\* End-fittings are 1/16" Waters-compatible.  
 \* For maximum operating pressure information, please refer to page 46.

# Inertsil<sup>®</sup> C8-3

## Physical Properties

- Silica : 3 Series High Purity Silica Gel
- Particle Size : 2  $\mu\text{m}$ , 3  $\mu\text{m}$ , 5  $\mu\text{m}$ , 10  $\mu\text{m}$
- Surface Area : 450  $\text{m}^2/\text{g}$
- Pore Size : 100  $\text{\AA}$  (10 nm)
- Pore Volume : 1.05 mL/g
- Bonded Phase : Octyl Groups
- End-capping : Yes
- Carbon Loading : 9 %
- USP Code : L7
- pH Range : 2 ~ 7.5




The same base silica gel and bonding technology that is used for Inertsil ODS-3 is also employed for Inertsil C8-3.


The difference between the two phases is just the length of the hydrocarbon ligands. As shown in figure 2, 10  $\mu\text{m}$  particle size columns are also available to meet the requirement of various pharmacopeia methods.

We recommend InertSustain C8 columns for all new method development.

**Figure 1 : Comparison of Retentivity between Inertsil<sup>®</sup> ODS-3 and Inertsil<sup>®</sup> C8-3**



**Figure 2 : Analysis of Losartan Potassium and Hydrochlorothiazide Tablets, Dissolution Test (Based on the Condition of United States Pharmacopeia 36-NF31)**



Analytical Columns

Particle Size: 2 µm	Length \ I.D. (mm)	2.1	3.0	
	30	5020-84930	5020-84935	
	50	5020-84931	5020-84936	
	75	5020-84932	5020-84937	
	100	5020-84933	5020-84938	
HPSeries Particle Size: 3 µm 50 MPa (500 Bar)	Length \ I.D. (mm)	2.1	3.0	4.6
	30	5020-14101	5020-14104	5020-14107
	50	5020-14102	5020-14105	5020-14108
	75	5020-14103	5020-14106	5020-14109
	100	5020-14031	5020-14034	5020-14037
	150	5020-14032	5020-14035	5020-14038
	250	5020-14033	5020-14036	5020-14039

\* End-fittings are 1/16" Waters-compatible.  
 \* UHPLC compatible end-fittings are also available upon request for UHPLC systems (Ex: UPLC) to avoid dead volume.  
 \* Indicate "UP Type end-fittings" when ordering. (Please note that UP type is not available for a 4.6 mm I.D. column)  
 \* For maximum operating pressure information, please refer to page 46.

Particle Size: 3 µm	Length \ I.D. (mm)	1.0	1.5		
	33	5020-84811	5020-84821		
	50	5020-84812	5020-84822		
	75	5020-84813	5020-84823		
	100	5020-84814	5020-84824		
	150	5020-13522	5020-13520		
	250	5020-	5020-		
	Length \ I.D. (mm)	2.1	3.0	4.0	4.6
	33	5020-04811	5020-04821	5020-04831	5020-04841
	50	5020-04812	5020-04822	5020-04832	5020-04842
	75	5020-04813	5020-04823	5020-04833	5020-01910
100	5020-04814	5020-04824	5020-01913	5020-04844	
125	5020-04817	5020-04827	5020-04837	5020-04845	
150	5020-04815	5020-04825	5020-04835	5020-01911	
250	5020-04816	5020-04826	5020-04836	5020-01912	
Particle Size: 5 µm	Length \ I.D. (mm)	1.0	1.5		
	33	5020-84911	5020-84921		
	50	5020-84912	5020-84922		
	75	5020-84913	5020-84923		
	100	5020-84914	5020-84924		
	150	5020-13512	5020-13510		
	250	5020-84916	5020-84926		
	Length \ I.D. (mm)	2.1	3.0	4.0	4.6
	33	5020-04911	5020-04921	5020-04931	5020-04941
	50	5020-04912	5020-04922	5020-04932	5020-04942
	75	5020-04913	5020-04923	5020-04933	5020-04943
100	5020-04914	5020-04924	5020-04934	5020-04944	
125	5020-04917	5020-04927	5020-04935	5020-04945	
150	5020-04915	5020-04925	5020-01902	5020-01900	
250	5020-04916	5020-04926	5020-01903	5020-01901	
Particle Size: 10 µm	Length \ I.D. (mm)	4.6			
	150	5020-01641			
	250	5020-01642			

\* End-fittings are 1/16" Waters-compatible.  
 \* For maximum operating pressure information, please refer to page 46.

Cartridge Guard Column E

I.D. of the Analytical Column Applicable (mm)	Length (mm)	I.D. (mm)	Replacement Cartridge E Guard Column (2 EA.)		Cartridge E Holder / Cartridge Set (2 Cartridge E Guard Columns & 1 Holder)	
			Particle Size		Particle Size	
			3 µm	5 µm	3 µm	5 µm
1.0	10	1.0	5020-19215	5020-19214	5020-19265	5020-19264
1.5, 2.1		1.5	5020-19315	5020-19314	5020-19365	5020-19364
2.1, 3.0		3.0	5020-19115	5020-19114	5020-19165	5020-19164
4.0, 4.6		4.0	5020-19015	5020-19014	5020-19065	5020-19064
2.1, 3.0	20	3.0	5020-19515	5020-19514	5020-19565	5020-19564
4.0, 4.6		4.0	5020-19415	5020-19414	5020-19465	5020-19464
Holder for Cartridge Guard Column E				For 10 mm Length		5020-08500
				For 20 mm Length		5020-08550

\* End-fittings are 1/16" Waters-compatible.  
 \* For maximum operating pressure information, please refer to page 46.

# Inertsil® C8

## Physical Properties

- Silica : 2 Series High Purity Silica Gel
- Particle Size : 5 µm
- Surface Area : 320 m<sup>2</sup>/g
- Pore Size : 150 Å (15 nm)
- Pore Volume : 1.20 mL/g
- Bonded Phase : Octyl Groups
- End-capping : Yes
- Carbon Loading : 10.5 %
- USP Code : L1
- pH Range : 2 ~ 7.5



Inertsil C8 columns have a pore size of 150 Å and it is recommended for rapid analysis of highly hydrophobic compounds. We recommend InertSustain C8 or Inertsil C8-4 columns for all new method development.

## Analytical Columns

	Length \ I.D. (mm)	2.1	3.0	4.0	4.6
Particle Size: 5 µm	150	5020-01221	5020-01222	5020-01223	5020-01224
	250	5020-01225	5020-01226	5020-01227	5020-01228

\* End-fittings are 1/16" Waters-compatible.

\* For maximum operating pressure information, please refer to page 46.

## Cartridge Guard Column E

I.D. of the Analytical Column Applicable (mm)	Length (mm)	I.D. (mm)	Replacement Cartridge E Guard Column (2 EA.)	Cartridge E Holder / Cartridge Set (2 Cartridge E Guard Columns & 1 Holder)
			Particle Size	Particle Size
			5 µm	5 µm
2.1, 3.0	10	3.0	5020-19136	5020-19186
4.0, 4.6		4.0	5020-19036	5020-19086
2.1, 3.0	20	3.0	5020-19536	5020-19586
4.0, 4.6		4.0	5020-19436	5020-19486
Holder for Cartridge Guard Column E		For 10 mm Length		5020-08500
		For 20 mm Length		5020-08550

\* End-fittings are 1/16" Waters-compatible.

\* For maximum operating pressure information, please refer to page 46.

