




PRODUCT LIST

Be selective



Food / Feed Safety



Environment



Cosmetics





Pharmaceutical R&D

AFFINIMIP® SPE PRODUCT LIST FOR MYCOTOXINS ANALYSES

Products	Designation	Description	25c/box	50 c/box
Zearalenone & Fumonisin	AFFINIMIP® SPE FumoZON	3mL for Zearalenone and Fumonisin	FS109-02	FS109-03
Patulin	AFFINIMIP® SPE Patulin	3mL – 100mg for Patulin	FS102-02	FS102-03
		6mL – 200mg for Patulin	FS102-02B-200mg	FS102-03B-200mg
	AFFINIMIP® SPE Patulin & Pectinase kit	Kit of 3mL cartridges for Patulin + 50mL Pectinase enzyme solution	FS102-02K	FS102-03K
	Kit of 6mL - 200mg cartridges for Patulin in dried apple + 50mL Pectinase enzyme solution	FS102-02KB-200mg	FS102-03KB-200mg	
Ochratoxin A	AFFINIMIP® SPE Ochratoxin A	3mL for Ochratoxin A	FS101-02	FS101-03
		6mL for Ochratoxin A	FS101-02B	FS101-03B
DON	AFFINIMIP® SPE Deoxynivalenol	6mL -100mg for Deoxynivalenol in food and babyfood	FS117-02B	FS117-03B
		6mL – 200mg for Deoxynivalenol in feed	FS117-02B-200mg	FS117-03B-200mg
Zearalenone	AFFINIMIP® SPE Zearalenone	3mL for ZON	FS100-02	FS100-03
Pectinase		50 mL Pectinase enzyme solution	REA-001-50mL	

AFFINIMIP® SPE PRODUCT LIST (MISCELLANEOUS)

Products	Designation	Description	25c/box	50 c/box
Bisphenol A and analogues	AFFINIMIP® SPE Bisphenols	3mL for Bisphenols (PP)	FS106-02	FS106-03
		6mL for Bisphenols (PP)	FS106-02B	FS106-03B
		6mL for Bisphenols (Glass)	FS106-02G	FS106-03G
Estrogens	AFFINIMIP® SPE Estrogens	1mL for Estrogens	FS104-02A	FS104-03A
		3mL for Estrogens	FS104-02	FS104-03
		96 well plate for estrogens– 1/pk	FS104-1.96W	
Catechol amines	AFFINIMIP® SPE Catecholamines	3mL for Catecholamines	DG100-02	DG100-03
		1mL for Catecholamines	DG100-02A	DG100-03A
Metanephrines	AFFINIMIP® SPE Metanephrines	3mL for Metanephrines	DG101-02	DG101-03
		1mL for Metanephrines	DG101-02A	DG101-03A
Picloram, Aminopyralid, Clopyralid	AFFINIMIP® SPE Picolinic Herbicides	3mL for Picolinic acid based herbicides	FS115-02	FS115-03
Glyphosate, AMPA	AFFINIMIP® SPE Glyphosate -AMPA	3mL for Glyphosate and AMPA	FS113-02	FS113-03
		6mL for Glyphosate and AMPA	FS113-02B	FS113-03B
NNAL	AFFINIMIP® SPE NNAL	3mL for NNAL	DG103-02	DG103-03
		96 well plate – 1/pk	DG103-1.96W	
Amphetamines	AFFINIMIP® SPE Amphetamines	3mL for Amphetamines derivatives	DG102-02	DG102-03
Chloramphenicol	AFFINIMIP® SPE Chloramphenicol	1mL for Chloramphenicol	FS110-02A	FS110-03A
		3mL for Chloramphenicol	FS110-02	FS110-03
Tamoxifen	AFFINIMIP® SPE Tamoxifen	3mL for Tamoxifen	PH101-02	PH101-03
Tetracyclines	AFFINIMIP® SPE Tetracyclines	1mL for Tetracyclines	FS112-02A	FS112-03A
		3mL for Tetracyclines	FS112-02	FS112-03
Zeranol Residues	AFFINIMIP® SPE Zeranol Residues	3mL for Zeranol Residues	FS105-02	FS105-03
Phenolics	AFFINIMIP® SPE Phenolics	3mL for Phenolic compounds	FS103-02	FS103-03
PAHs	AFFINIMIP® SPE PAHs	3mL for PAHs	FS119-02	FS119-03

Product	Description	Reference	Number of cartridges
Reversible cartridges (2mL) 	SPE cartridges compatible with automates for purification of radiotracers issued from Aromatic Nucleophilic Substitution	RP100-01	10
	SPE cartridges compatible with automates for purification of radiotracers issued from Aromatic Nucleophilic Substitution	RP100-02	25
Reversible cartridges (0.7mL) 	SPE cartridges compatible with automates for purification of radiotracers issued from Aromatic Nucleophilic Substitution	RP100A-01	10
	SPE cartridges compatible with automates for purification of radiotracers issued from Aromatic Nucleophilic Substitution	RP100A-02	25

AttractSPE™ PRODUCT LIST

Format, amount	#/box	AttractSPE™ HLB	AttractSPE™ SCX	AttractSPE™ WCX	AttractSPE™ SAX	AttractSPE™ WAX	AttractSPE™ DVB
1mL, 30mg	100	HLB-100.S.1.30	SCX-100.S.1.30	WCX-100.S.1.30	SAX-100.S.1.30	WAX-100.S.1.30	DVB-100.S.1.30
3mL, 60mg	50	HLB-50.S.3.60	SCX-50.S.3.60	WCX-50.S.3.60	SAX-50.S.3.60	WAX-50.S.3.60	DVB-50.S.3.60
	100	HLB-100.S.3.60	SCX-100.S.3.60	WCX-100.S.3.60	SAX-100.S.3.60	WAX-100.S.3.60	DVB-100.S.3.60
6mL, 200mg	25	HLB-25.S.6.200	SCX-25.S.6.200	WCX-25.S.6.200	SAX-25.S.6.200	WAX-25.S.6.200	DVB-25.S.6.200
	50	HLB-50.S.6.200	SCX-50.S.6.200	WCX-50.S.6.200	SAX-50.S.6.200	WAX-50.S.6.200	DVB-50.S.6.200
	100	HLB-100.S.6.200	SCX-100.S.6.200	WCX-100.S.6.200	SAX-100.S.6.200	WAX-100.S.6.200	DVB-100.S.6.200
6mL, 500mg	25	HLB-25.S.6.500	SCX-25.S.6.500	WCX-25.S.6.500	SAX-25.S.6.500	WAX-25.S.6.500	DVB-25.S.6.500
	50	HLB-50.S.6.500	SCX-50.S.6.500	WCX-50.S.6.500	SAX-50.S.6.500	WAX-50.S.6.500	DVB-50.S.6.500
	100	HLB-100.S.6.500	SCX-100.S.6.500	WCX-100.S.6.500	SAX-100.S.6.500	WAX-100.S.6.500	DVB-100.S.6.500
12mL, 500mg	25	HLB-25.S.12.500	SCX-25.S.12.500	WCX-25.S.12.500	SAX-25.S.12.500	WAX-25.S.12.500	DVB-25.S.12.500
20mL, 1g	25	HLB-25.S.20.1g	SCX-25.S.20.1g	WCX-25.S.20.1g	SAX-25.S.20.1g	WAX-25.S.20.1g	DVB-25.S.20.1g
96 wells Plate, 30mg	1	HLB-1.96W.30	SCX-1.96W.30	WCX-1.96W.30	SAX-1.96W.30	WAX-1.96W.30	DVB-1.96W.30
Reversible 0.7mL, 30mg	25	HLB-25.REV.1.N10	SCX-25.REV.1.N10	WCX-25.REV.1.N10	SAX-25.REV.1.N10	WAX-25.REV.1.N10	DVB-25.REV.1.N10
	50	HLB-50.REV.1.N10	SCX-50.REV.1.N10	WCX-50.REV.1.N10	SAX-50.REV.1.N10	WAX-50.REV.1.N10	DVB-50.REV.1.N10
Reversible 0.7mL, 100mg (200mg for DVB)	25	HLB-25.REV.1.F	SCX-25.REV.1.F	WCX-25.REV.1.F	SAX-25.REV.1.F	WAX-25.REV.1.F	DVB-25.REV.1.F
	50	HLB-50.REV.1.F	SCX-50.REV.1.F	WCX-50.REV.1.F	SAX-50.REV.1.F	WAX-50.REV.1.F	DVB-50.REV.1.F
Reversible 2mL, 225mg	25	HLB-25.REV.2.N10	SCX-25.REV.2.N10	WCX-25.REV.2.N10	SAX-25.REV.2.N10	WAX-25.REV.2.N10	DVB-25.REV.2.N10
	50	HLB-50.REV.2.N10	SCX-50.REV.2.N10	WCX-50.REV.2.N10	SAX-50.REV.2.N10	WAX-50.REV.2.N10	DVB-50.REV.2.N10

AttractSPE™ PRODUCT LIST - removal of interferences

Cartridges format, Sorbent amount	#/box	AttractSPE™ SAX-HCO3	AttractSPE™ PS-H	AttractSPE™ PS-Ag	AttractSPE™ PS-Ba	AttractSPE™ IDA
1mL	100	SAX-HCO3-100.S.1.30		PSAg-100.S.1.30	PSBa-100.S.1.30	IDA-100.S.1.30
3mL, 60mg	25	SAX-HCO3-25.S.3.60	PSH-25.S.3.60	PSAg-25.S.3.60	PSBa-25.S.3.60	IDA-25.S.3.60
	50	SAX-HCO3-50.S.3.60	PSH-50.S.3.60	PSAg-50.S.3.60	PSBa-50.S.3.60	IDA-50.S.3.60
6mL, 200mg	25	SAX-HCO3-25.S.6.200	PSH-25.S.6.200	PSAg-25.S.6.200	PSBa-25.S.6.200	IDA-25.S.6.200
	50	SAX-HCO3-50.S.6.200	PSH-50.S.6.200	PSAg-50.S.6.200	PSBa-50.S.6.200	IDA-50.S.6.200
6mL, 500mg	25	SAX-HCO3-25.S.6.500	PSH-25.S.6.500	PSAg-25.S.6.500	PSBa-25.S.6.500	IDA-25.S.6.500
	50	SAX-HCO3-50.S.6.500	PSH-50.S.6.500	PSAg-50.S.6.500	PSBa-50.S.6.500	IDA-50.S.6.500
96 wells Plate	1	SAX-HCO3-1.96W.30	PSH-1.96W.30			IDA-1.96W.30
Reversible 0.7mL, 30mg	25	SAX-HCO3-25.REV.1.F	PSH-25.S.REV.1.F			IDA-25.REV.1.N10
	50	SAX-HCO3-50.REV.1.F	PSH-50.S.REV.1.F			IDA-50.REV.1.N10
Reversible 0.7mL, 100mg	25	SAX-HCO3-25.REV.2.F	PSH-25.S.REV.2.F	PSAg-25.S.REV.1.F For 400mg	PSBa-25.S.REV.1.F For 400mg	IDA-25.REV.1.F
	50	SAX-HCO3-50.REV.2.F	PSH-50.S.REV.2.F	PSAg-50.S.REV.1.F for 400mg	PSBa-50.S.REV.1.F For 400mg	IDA-50.REV.1.F
Reversible 2mL, 800mg	25	SAX-HCO3-100.S.1.30				IDA-25.REV.2.F
	50	SAX-HCO3-25.S.3.60	PSH-25.S.3.60			IDA-50.REV.2.F

AttractSPE™ PRODUCT LIST - miscellaneous

AttractSPE™ Carbon based SPE - Product list

Product	Vol	Sorbent	25 cartridges/box	50 cartridges/box
AttractSPE™ Carbon	6mL	500mg	Carb-25.S.6.500	Carb-50.S.6.500
AttractSPE™ Carbon/PSA	3mL	250mg/ 250mg	CarbPSA- 25.S.3.250.250	CarbPSA- 50.S.3.250.250
	6mL	500mg/ 500mg	CarbPSA- 25.S.6.500.500	CarbPSA- 50.S.6.500.500
AttractSPE™ Carbon/Amine	6mL	500mg/ 500mg	CarbNH2- 25.S.6.500.500	CarbNH2- 50.S.6.500.500

AttractSPE™ LipRem

Cartridges format, Sorbent amount	#/box	AttractSPE™ LipRem
1mL, 20mg	100	LipRem-100.S.1.20
3mL, 60mg	25	LipRem-25.S.3.50
	50	LipRem-50.S.3.50
6mL, 100mg	25	LipRem-25.S.6.100
	50	LipRem-50.S.6.100
96 wells Plate	1	LipRem-1.96W.20
Reversible 0.7mL, 100mg	25	LipRem-1.REV.1.F
	50	LipRem-1.REV.1.F

SilactSPE™ PRODUCT LIST

		Non polar sorbents			Polar sorbents			
Cartridges format, Sorbent amount	#/box	SilactSPE™ C18	SilactSPE™ C8	SilactSPE™ Phenyl	SilactSPE™ Silica	SilactSPE™ Cyano	SilactSPE™ Florisol	SilactSPE™ Florisol PR
1mL, 50mg	100	C18-100.S.1.50	C8-100.S.1.50	Phe-100.S.1.50	Si-100.S.1.50	CN-100.S.1.50	Flo-100.S.1.50	FloPR-100.S.1.50
1mL, 100mg	100	C18-100.S.1.100	C8-100.S.1.100	Phe-100.S.1.100	Si-100.S.1.100	CN-100.S.1.100	Flo-100.S.1.100	FloPR-100.S.1.100
3mL, 200mg	50	C18-50.S.3.200	C8-50.S.3.200	Phe-50.S.3.200	Si-50.S.3.200	CN-50.S.3.200	Flo-50.S.3.200	FloPR-50.S.3.200
3mL, 500mg	50	C18-50.S.3.500	C8-50.S.3.500	Phe-50.S.3.500	Si-50.S.3.500	CN-50.S.3.500	Flo-50.S.3.500	FloPR-50.S.3.500
6mL, 500mg	50	C18-50.S.6.500	C8-50.S.6.500	Phe-50.S.6.500	Si-50.S.6.500	CN-50.S.6.500	Flo-50.S.6.500	FloPR-50.S.6.500
6mL, 1g	50	C18-50.S.6.1g	C8-50.S.6.1g	Phe-50.S.6.1g	Si-50.S.6.1g	CN-50.S.6.1g	Flo-50.S.6.1g	FloPR-50.S.6.1g
6mL, 2g	50	C18-50.S.6.2g	C8-50.S.6.2g	Phe-50.S.6.2g	Si-50.S.6.2g	CN-50.S.6.2g	Flo-50.S.6.2g	FloPR-50.S.6.2g
12mL, 2g	20	C18-20.S.12.2g	C8-20.S.12.2g	Phe-20.S.12.2g	Si-20.S.12.2g	CN-20.S.12.2g	Flo-20.S.12.2g	FloPR-20.S.12.2g
Reversible 0.7mL, 200mg	25	C18-25.REV.1.200	C8-25.REV.1.200	Phe-25.REV.1.200	Si-25.REV.1.200	CN-25.REV.1.200	Flo-25.REV.1.200	FloPR-25.REV.1.200
Reversible 2mL, 750mg	25	C18-25.REV.2.750	C8-25.REV.2.750	Phe-25.REV.2.750	Si-25.REV.2.750	CN-25.REV.2.750	Flo-25.REV.2.750	FloPR-25.REV.2.750

For other formats, please contact us

SilactSPE™ PRODUCT LIST (continued)

Cartridges format, Sorbent amount	#/box	Polar sorbents			Others sorbents			
		SilactSPE™ Alumina Acidic	SilactSPE™ Alumina Neutral	SilactSPE™ Alumina Basic	SilactSPE™ Amine	SilactSPE™ PSA	SilactSPE™ Carbonate	SilactSPE™ Hydroxy Apatite
1mL, 50mg	100	AluA- 100.S.1.50	AluN- 100.S.1.50	AluB- 100.S.1.50	NH2- 100.S.1.50	PSA- 100.S.1.50	CO3- 100.S.1.50	HAp- 100.S.1.50
1mL, 100mg	100	AluA- 100.S.1.100	AluN- 100.S.1.100	AluB- 100.S.1.100	NH2- 100.S.1.100	PSA- 100.S.1.100	CO3- 100.S.1.100	
3mL, 200mg	50	AluA- 50.S.3.200	AluN- 50.S.3.200	AluB- 50.S.3.200	NH2- 50.S.3.200	PSA- 50.S.3.200	CO3- 50.S.3.200	HAp- 50.S.3.200
3mL, 500mg	50	AluA- 50.S.3.500	AluN- 50.S.3.500	AluB- 50.S.3.500	NH2- 50.S.3.500	PSA- 50.S.3.500	CO3- 50.S.3.500	
6mL, 500mg	50	AluA- 50.S.6.500	AluN- 50.S.6.500	AluB- 50.S.6.500	NH2- 50.S.6.500	PSA- 50.S.6.500	CO3- 50.S.6.500	HAp- 50.S.6.500
6mL, 1g	50	AluA- 50.S.6.1g	AluN- 50.S.6.1g	AluB- 50.S.6.1g	NH2- 50.S.6.1g	PSA- 50.S.6.1g	CO3- 50.S.6.1g	
6mL, 2g	50	AluA- 50.S.6.2g	AluN- 50.S.6.2g	AluB- 50.S.6.2g	NH2- 50.S.6.2g	PSA- 50.S.6.2g	CO3- 50.S.6.2g	
12mL, 2g	20	AluA- 20.S.12.2g	AluN- 20.S.12.2g	AluB- 20.S.12.2g	NH2- 20.S.12.2g	PSA- 20.S.12.2g	CO3- 20.S.12.2g	
Reversible 0.7mL, 200mg	25	AluA- 25.REV.1.200	AluN- 25.REV.1.200	AluB- 25.REV.1.200	NH2- 25.REV.1.200	PSA- 25.REV.1.200	CO3- 25.REV.1.200	HAp- 50.REV.1.F
Reversible 2mL, 750mg	25	AluA- 25.REV.2.750	AluN- 25.REV.2.750	AluB- 25.REV.2.750	NH2- 25.REV.2.750	PSA- 25.REV.2.750	CO3- 25.REV.2.750	

For other formats, please contact us

SilactSPE™ PRODUCT LIST - miscellaneous SPE

SPE for Polycyclic Aromatic Hydrocarbons (PAHs) in soil

Product	Vol	Sorbent	25 cartridges/box	50 cartridges/box
SilactSPE™ CN/SiOH	3mL	500mg/1g	CNSiOH-25.S.3.500.1g	CNSiOH-50.S.3.500.1g
	6mL	500mg/1g	CNSiOH-25.S.6.500.1g	CNSiOH-50.S.6.500.1g
	6mL glass	500mg/1g	CNSiOH-25.G.6.500.1g	CNSiOH-50.G.6.500.1g

SilactSPE™ SLE

Cartridge volume	Sorbent	25 cartridges/box	50 cartridges/box
1mL	250mg	SLE-25.S.1.250	SLE-50.S.1.250
3mL	500mg	SLE-25.S.3.500	SLE-50.S.3.500
6mL	1g	SLE-25.S.6.1g	SLE-50.S.6.1g
15mL	3g	SLE-25.S.15.3g	SLE-50.S.15.3g
30mL	4.5g	SLE-25.S.30.4g	SLE-50.S.30.4g
70mL	14.5g	SLE-25.S.70.14g	SLE-50.S.70.14g

Fritted cartridges

Cartridge volume	SilactSPE™ Double fritted 100 cartridges	SilactSPE™ Single fritted 100 cartridges
1mL	0-100.S.1.2F	0-100.S.1.1F
3mL	0-100.S.3.2F	0-100.S.3.1F
6mL	0-100.S.6.2F	0-100.S.6.1F
15mL	0-100.S.15.2F	0-100.S.15.1F
25mL	0-100.S.25.2F	0-100.S.25.1F
60mL	0-100.S.60.2F	0-100.S.60.1F

Qcleanup™ - PRODUCT LIST

Qcleanup™ EXTRACTION SALTS

QuEChERS methods	Description	Pouches / box	Product reference
Original method	4g MgSO ₄ 1g NaCl	50	EXT.ORL.50
EN 15662	1g Trisodium citrate Dihydrate 0.5g Disodium hydrogencitrate sesquihydrate 1g NaCl and 4g MgSO ₄	50	EXT.EN.50
AOAC 2007.01	1.5g Sodium Acetate and 6g MgSO ₄	50	EXT.AOAC.50

Qcleanup™ DISPERSIVE SPE PRODUCTS

Method	Description	Nber/box	Product reference
For General Fruits & Vegetables			
EN 15662	150mg MgSO ₄ + 25mg PSA	100 tubes of 2mL	dSPE.EN.GFV.100.2
	900mg MgSO ₄ + 150mg PSA	50 tubes of 15mL	dSPE.EN.GFV.50.15
AOAC 2007.01	150mg MgSO ₄ + 50mg PSA	100 tubes of 2mL	dSPE.AOAC.GFV.100.2
	1200mg MgSO ₄ + 400mg PSA	50 tubes of 15mL	dSPE.AOAC.GFV.50.15
For Pigmented Fruits & Vegetables			
EN 15662	150mg MgSO ₄ + 25mg PSA + 2.5mg CB	100 tubes of 2mL	dSPE.EN.PFV.100.2
	900mg MgSO ₄ + 150mg PSA + 15mg CB	50 tubes of 15mL	dSPE.EN.PFV.50.15
AOAC 2007.01	150mg MgSO ₄ + 50mg PSA + 50mg CB	100 tubes of 2mL	dSPE.AOAC.PFV.100.2
	1200mg MgSO ₄ + 400mg PSA + 400mg CB	50 tubes of 15mL	dSPE.AOAC.PFV.50.15
For Highly Pigmented and Fatty Fruits & Vegetables			
EN 15662	150mg MgSO ₄ + 25mg PSA + 7.5mg CB	100 tubes of 2mL	dSPE.EN.HPFV.100.2
	900mg MgSO ₄ + 150mg PSA + 45mg CB	50 tubes of 15mL	dSPE.EN.HPFV.50.15
AOAC 2007.01	150mg MgSO ₄ + 50mg PSA + 50mg CB + 50mg C18	100 tubes of 2mL	dSPE.AOAC.HPFV.100.2
	1200mg MgSO ₄ + 400mg PSA + 400mg CB + 400mg C18	50 tubes of 15mL	dSPE.AOAC.HPFV.50.15
For Fatty and waxed Fruits & Vegetables			
EN 15662	150mg MgSO ₄ + 25mg PSA + 25mg C18	100 tubes of 2mL	dSPE.EN.FWFV.100.2
	900mg MgSO ₄ + 150mg PSA + 150mg C18	50 tubes of 15mL	dSPE.EN.FWFV.50.15
AOAC 2007.01	150mg MgSO ₄ + 50mg PSA + 50mg C18	100 tubes of 2mL	dSPE.AOAC.FWFV.100.2
	1200mg MgSO ₄ + 400mg PSA + 400mg C18	50 tubes of 15mL	dSPE.AOAC.FWFV.50.15

POCIS PRODUCT LIST

Designation	Definition	Composition	Reference
AFFINIMIP® POCIS GLYPHOSATE	POCIS containing AFFINIMIP® GLYPHOSATE - AMPA for the retention of glyphosate and AMPA	1 POCIS	POCIS.GLY.90.55.A.1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.GLY.90.55.kit.10
		Kit of 50 POCIS + empty fritted cartridges	POCIS.GLY.90.55.kit.50
AFFINIMIP® POCIS EDC	POCIS containing AFFINIMIP® Estrogens and AFFINIMIP® Bisphenols for the retention of endocrine disrupters such as natural/synthetic estrogens, Bisphenols...	1 POCIS	POCIS.EDC.90.55.A.1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.EDC.90.55.kit.10
		Kit of 50 POCIS + empty fritted cartridges	POCIS.EDC.90.55.kit.50
Attract POCIS Pesticides	POCIS containing mixture of sorbent for the retention of several pesticides	1 POCIS	POCIS.PEST.90.55.A.1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.PEST.90.55.kit.10
		Kit of 50 POCIS + empty fritted cartridges	POCIS.PEST.90.55.kit.50
Attract POCIS HLB	POCIS containing Attract HLB for the retention of pharmaceutical drug residues	1 POCIS	POCIS.HLB.90.55.A.1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.HLB.90.55.kit.10
		Kit of 50 POCIS + empty fritted cartridges	POCIS.HLB.90.55.kit.50
CANISTER – 3 POCIS	Canister for 3 POCIS . Requires a holder	1 canister	CAN-3P.A.1
HOLDER – 3 POCIS	Holder for 3 POCIS	1 holder	HOLD-3P.A.1



POCIS



CANISTER – 3 POCIS



HOLDER – 3 POCIS

On-line SPE columns and SPE ACCESSORIES

On-line SPE columns – Product list

Product	Product reference	Nber column	I.D. (mm)	Lenght (mm)
On-line AttractSPE™ HLB columns	OnlineSPE-HLB-1.2.20	1	2.1	20
	OnlineSPE-HLB-1.5.20	1	4.6	20
On-line AFFINIMIP® PHENOLICS columns	OnlineSPE-PHE-1.2.20	1	2.1	20
	OnlineSPE-PHE-1.5.20	1	4.6	20
On-line AFFINIMIP® ESTROGENS columns	OnlineSPE-EST-1.2.20	1	2.1	20
	OnlineSPE-EST-1.5.20	1	4.6	20

SPE ACCESSORIES – Product list

SPE Accessories	Designation	Definition	Reference
Manifold	SPE Vaccum Manifold	12-port model	ACC-MAN1
SPE Adapter & Reservoir kit	SPE Adapter & Reservoir kit	Kit of 12 reservoirs 60ml and adapters for use with 1,3 & 6 mL cartridges	ACC-AR1
Mini-Vap	Mini Evaporator/Concentrator	6 port Mini-Vap Evaporator/Concentrator for use with 1 to 250mL containers	ACC-VAP1
Mini PUMP	Mini vacuum pump	Laboport diaphragm vacuum mini pump, 5.5L/min	ACC-PUMP
Vacuum pump trap	SPE Vacuum pump trap kit	1L trap kit	ACC-TRAP