

Septa for GC Injectors

- Center point guides the needle for easy penetration
- Reduce needle bending
- Precision molding assures accurate fit
- Available in selected sizes of our premium septa
- Packaged in glass vials for high purity



Premium Septa Advanced Green 3			Premium Septa BTO		
<p><i>The Advanced Green 3 Septum was specifically created to combine significantly longer injection life, low bleed and low injection port adhesion. The result is a green septum you can use for all your critical analysis. Septa are 3mm thick, except Shimadzu-style plug.</i></p>			<p><i>Septa BTO® is <u>B</u>leed and <u>T</u>emperature <u>O</u>ptimized for today's most demanding GC and GC/MS applications. Septa BTO is formulated to extend low-bleed and outstanding mechanical properties of premium GC septa to the highest-temperature applications (400°C with most modern GCs).</i></p>		
Part No.	Description	Qty.	Part No.	Description	Qty.
SA5	5mm CenterGuide	pk.50	SB5	5mm CenterGuide	pk.50
SA9	9mm CenterGuide	pk.50	SB64	6.4mm (1/4")	pk.50
SA95	9.5mm (3/8")	pk.50	SB9	9mm CenterGuide	pk.50
SA11	11mm (7/16") CenterGuide	pk.50	SB95	9.5mm (3/8")	pk.50
SA11-100	11mm (7/16") CenterGuide	pk.100	SB10	10mm	pk.50
SA115	11.5mm CenterGuide	pk.50	SB11	11mm (7/16") CenterGuide	pk.50
SA127	12.7mm (1/2") CenterGuide	pk.50	SB11-100	11mm (7/16") CenterGuide	pk.100
SA17	17mm	pk.50	SB115	11.5mm CenterGuide	pk.50
SAS	for Shimadzu	pk.50	SB127	12.7mm (1/2") CenterGuide	pk.50
			SB17	17mm CenterGuide	pk.50
			SBS	for Shimadzu	pk.50
<p>Premium Septa Marathon</p> <p><i>The Marathon septum typically achieves 400 injections without failure, when used with a rounded-tip (Agilent-style) needle and autosampler or needle guide. The CenterGuide versions of Marathon are pre-pierced to channel the needle through the same spot with each injection, which helps eliminate the common type of coring that occurs on the first injection. The Marathon septum is also ideal for standard manual-injection GC and GC/MS. Septa are 3mm thick, except Shimadzu-style plug septa.</i></p>			<p>Premium Septa Silicone Pure</p> <p><i>The Silicone Pure Septum was specifically created to combine significantly longer injection life, low bleed and low injection port adhesion. Suitable for inlet temperatures of up to 340°C.</i></p>		
Part No.	Description	Qty.	Part No.	Description	Qty.
SM5-50	5mm CenterGuide	pk.50	SS95	9.5mm	pk.25
SM9-50	9mm CenterGuide	pk.50	SS10	10mm	pk.25
SM11-50	11mm (7/16") CenterGuide	pk.50	SS11	11mm	pk.25
SM11-100	11mm (7/16") CenterGuide	pk.100	SS125	12.5mm	pk.25
SM115-50	11.5mm CenterGuide	pk.50	SS17	17mm	pk.25
SM127-50	12.7mm (1/2") CenterGuide	pk.50	SSS	for Shimadzu	pk.25
SM17-50	17mm	pk.50			
SMS-50	for Shimadzu	pk.50			

Septa Pick

The Septa Pick is useful for removing small pieces of septa lodged in the threads and corners of the septum cap or injection port. It can also be used to remove small bits of graphite ferrule from the cap and fittings of column connectors.

Part No.	Description	Qty.
SP	Septa Pick	ea.



SETonic Microsyringes
see page 79-82.