

Filtration and microfiltration products for applications in beer elaboration and laboratories

The following filter papers and microfiltration products are ideal to elaborate, clarify and prepare beer before further analysis.



PROCESS	TECHNIQUE	TYPE OF FILTER	REF. CHM
Solid - Liquid Separation	Malt and beer Analysis (EBC)	Qualitative Low Ash filter paper	F1001- F1002- F1003- F1004- F1005-
	Filtration (Funnel/E.Buchner)	Qualitative General Purpose Filter Paper	F1093-
Determination of suspended solids		Glass Microfiber Filter	GF3-
Aqueous samples clarification	Malta Color. Spectrophotometry.	Cellulose Acetate Syringe Filter, 0,45um	SCA045025Q
Nitrogen Cuantification	Nitrogen Determination	Weighing boats	
Micro organisms Retention		Cellulose Nitrate membrane. White with black grid. Sterile individually packed. 0,45um, 47mm.	MNW045047H-SG
Yeast count	Microbiological Analysis	Cellulose Nitrate membrane. Grey with white grid. Sterile individually packed. 0,45 and 0,65um, 47mm.	MNB045047H-SW MNB065047H-SW
Micro organisms Incubation		Petri dishes with absorbent Pad. (Sterile and non-sterile version)	SSP050T-S

Filtration and microfiltration products for applications in beer elaboration and laboratories